



Edition: 11.24
Valid until: 11.25

BasiX²

Invoice to:
Name
Street
Town: Postal code:
E-Mail:
Tel: Fax:

Order date
Purchase order:
Number required:

- ☐ Order
☐ Calculation
☐ Built 4 me



Deliver to:
Name
Street
Town
E-Mail:
Tel:

Marked for:
Sunrise order no:

Basic Price:
Reimbursement code:
Max. user weight: 125 kg
USER WEIGHT kg
Mandatory information to ensure correct configuration
Art. no.: 074100-000
Dimensional tolerances +/- 1 cm

☐ = Standard ☐ = Cost option ☐ = No cost option

Sunrise Medical Customer service

Address: Kahlbachring 2-4, 69254 Malsch, Germany
Tel. +497253 980-0
Fax. +497253 980-222

E-mail: export@sunrisemedical.de
Website: www.sunrisemedical.eu

BASIC PRICE 845 €

Standard measures (all measurements are in mm and without cushion):

Seat height: 475 mm (450 mm with 12"); adj. from 425 to 505 mm
Seat depth: 435 mm as standard; adjustable to 410 to 460 mm
Back height: 435 mm as standard; adjustable to 410 and 460 mm
Colour: brilliant silver , upholstery black, fork black
Crash tested

Standard features:

Seat upholstery padded Nylon, black
Back upholstery padded Nylon, tension adjustable, black
Backrest with fixed push handles
Standard armrest, flip back, detachable with short & long arm pad
Hanger 70° with composite angle adj footplate, heel loop
Castor 8" x 1 1/4" solid
Fork, aluminium, black
Rear wheel 24" x 1", aluminium handrim, with quick release axle, solid tyre or 12" (UK only)
Push to lock knee lever brake
Step tube integrated in frame
Colour: brilliant silver
Crash tested

BasiX² Std Back, Fixed Pushhandles ,8" Caster, 24" Rear Wheel

074100-001	☐	SW 380 mm	-	
074100-002	☐	SW 405 mm	-	
074100-003	☐	SW 430 mm	-	
074100-004	☐	SW 455 mm	-	
074100-005	☐	SW 480 mm	-	
074100-006	☐	SW 520 mm	-	

BasiX² Std Back, Fixed Pushhandles, 8" Caster, 24" Rear Wheel, Drumbrake

074101-001	☐	SW 380 mm	293 €	
074101-002	☐	SW 405 mm	293 €	
074101-003	☐	SW 430 mm	293 €	
074101-004	☐	SW 455 mm	293 €	
074101-005	☐	SW 480 mm	293 €	
074101-006	☐	SW 520 mm	293 €	

BasiX² Std Back, Integrated Height Adj. Pushhandles, 8" Caster, 24" Rear Wheel, Cushion, Crutch Holder

074108-002	☐	SW 405 mm	238 €	
074108-003	☐	SW 430 mm	238 €	
074108-004	☐	SW 455 mm	238 €	
074108-005	☐	SW 480 mm	238 €	
074108-006	☐	SW 520 mm	238 €	

BasiX ² Std Back, Fixed Pushhandles, Height adj. Armrest, 8" Caster, 24" Rear Wheel, 3" Travelwheels			
074102-001	☐	SW 380 mm	250 €
074102-002	☐	SW 405 mm	250 €
074102-003	☐	SW 430 mm	250 €
074102-004	☐	SW 455 mm	250 €
074102-005	☐	SW 480 mm	250 €
074102-006	☐	SW 520 mm	250 €
SEAT Cushion 50 mm			
064335-001	☐	SW 380 mm	30 €
064335-002	☐	SW 405 mm	30 €
064335-003	☐	SW 430 mm	30 €
064335-004	☐	SW 455 mm	30 €
064335-005	☐	SW 480 mm	30 €
064335-006	☐	SW 520 mm	30 €
BACKREST Options			
064707-001	☐	Swing away stabilizer bar SW 380 mm (not with folding back)	132 €
064707-002	☐	Swing away stabilizer bar SW 405 mm (not with folding back)	132 €
064707-003	☐	Swing away stabilizer bar SW 430 mm (not with folding back)	132 €
064707-004	☐	Swing away stabilizer bar SW 455 mm (not with folding back)	132 €
064707-005	☐	Swing away stabilizer bar SW 480 mm (not with folding back)	132 €
064707-006	☐	Swing away stabilizer bar SW 520 mm (not with folding back)	132 €
064909-001	☐	Headrest (with swing away stabilizer bar 380 mm) (not with folding back)	433 €
064909-002	☐	Headrest (with swing away stabilizer bar 405 mm) (not with folding back)	433 €
064909-003	☐	Headrest (with swing away stabilizer bar 430 mm) (not with folding back)	433 €
064909-004	☐	Headrest (with swing away stabilizer bar 455 mm) (not with folding back)	433 €
064909-005	☐	Headrest (with swing away stabilizer bar 480 mm) (not with folding back)	433 €
064909-006	☐	Headrest (with swing away stabilizer bar 520 mm) (not with folding back)	433 €
064920-001-ET	☐	Extension back SW 405 mm (not with int. height adj. push handles and SW 38.52)	130 €
064920-001-ET	☐	Extension back SW 430 mm (not with int. height adj. push handles and SW 38.52)	130 €
064920-001-ET	☐	Extension back SW 455 mm (not with int. height adj. push handles and SW 38.52)	130 €
064920-001-ET	☐	Extension back SW 480 mm (not with int. height adj. push handles and SW 38.52)	130 €
ARMREST			
064475-100	☐	Height adjustable armrest with short & long pad (sliding) (not with recliner)	202 €
LEGREST & FOOTPLATES			
064897-011	☐	Elevating Leg Rest (Left) SW 380-405 mm (with curved calf pad)	148 €
064897-021	☐	Elevating Leg Rest (Right) SW 380-405 mm (with curved calf pad)	148 €
064897-101	☐	Elevating Leg Rest (Pair) SW 380-405 mm (with curved calf pad)	296 €
064897-012	☐	Elevating Leg Rest (Left) SW 430-455 mm (with curved calf pad)	148 €
064897-022	☐	Elevating Leg Rest (Right) SW 430-455 mm (with curved calf pad)	148 €
064897-102	☐	Elevating Leg Rest (Pair) SW 430-455 mm (with curved calf pad)	296 €
064897-013	☐	Elevating Leg Rest (Left) SW 480-520 mm (with curved calf pad)	148 €
064897-023	☐	Elevating Leg Rest (Right) SW 480-520 mm (with curved calf pad)	148 €
064897-103	☐	Elevating Leg Rest (Pair) SW 480-520 mm mm (with curved calf pad)	296 €
064280-002	☐	Amputee support pad (Left)	216 €
064280-001	☐	Amputee support pad (Right)	216 €
067725	☐	Cast support (left and right usable)	130 €
CASTOR			
064936	☐	8" x 50 solid (always order 2 units)	74 €
WHEEL LOCKS			
000861253	☐	Lever extension (pair)	23 €
ACCESSORIES			
064780-001	☐	reimbursement kit SW380mm (antitip (pair)/cushion/belt/crutch holder)	184 €
064780-002	☐	reimbursement kit SW405mm (antitip (pair)/cushion/belt/crutch holder)	184 €
064780-003	☐	reimbursement kit SW430mm (antitip (pair)/cushion/belt/crutch holder)	184 €
064780-004	☐	reimbursement kit SW455mm (antitip (pair)/cushion/belt/crutch holder)	184 €
064780-005	☐	reimbursement kit SW480mm (antitip (pair)/cushion/belt/crutch holder)	184 €
064780-006	☐	reimbursement kit SW520mm (antitip (pair)/cushion/belt/crutch holder)	184 €
067708	☐	Crutch holder (UniX)	32 €
062553-033	☐	Transparent Spoke guard (24")	192 €
064559-100	☐	3" travel wheels (pair)	245 €
760848050.01	☐	Slide on tray (SW 380-410)	257 €
760848050.02	☐	Slide on tray (SW430)	298 €
760848050.03	☐	Slide on tray (SW 455-480)	298 €
064723-ET	☐	Anti tips (Pair)	298 €
064723	☐	Anti tips (Unit)	62 €
064931	☐	Step tube (one unit)	20 €
069935	☐	Positioning belt	38 €
Breezy BasiX² - Basic Price			
Option price			
Packaging and transport			
Exclusive VAT			
VAT		€	
Total price		€	
		For more information visit our webpage www.sunrisemedical.eu	
		© Sunrise Medical 2025	