

BREEZY STYLE

LIGHTNESS,
VERSATILITY
AND STYLE



LIGHTNESS, VERSATILITY AND STYLE

The aluminium chair with the
best reliability and style!



 **STYLE**



Configuration with options



THE STEEL WHEELCHAIR IN WHICH YOU CAN TRUST

RIGIDITY AND RESISTANCE

More than 25 years of experience in design and manufacture of manual wheelchair endorse us. Designed with an **aluminium tube with 25mm**, our Breezy Style is the chair in which you can trust.

In its manufacture we have **only use the highest quality components** to guarantee a strong model, responding day to day requirements.

Details as the **robust armrest attachments to the frame**, the efficient design of the **brakes with an ergonomic handle** or the **great resistance of the anchoring system of the footrest platform** are the proof of this.

And there are still more! The **aluminium plate** and the **cross brace**, without screw and anchors, equip the chair with a great resistance and rigidity to be able to enjoy it as the first day, for a long time. We are so sure that we offer **5 year warranty on frame and cross brace!**



Aluminium plate



Robust armrest



Brakes with an ergonomic system



Great resistance on the footrest platform

CONFIGURABLE OPTIONS

Breezy Style allows 3 different height from the seat to the floor (51,5, 49 o 46,5 cm). In addition allows adjustments in:

- Backrest height between 40 and 45 cm, relocating the upholstery.
- Padded armrest adjustable in depth

THE PERFECT CHOICE

Breezy Style is a **versatile** chair, and it can respond to any user needs. Its numerous specially-made options are from being able to choose between 6 seat width options, 3 types of backrest (Standard, Half-folding or Reclinable), big or small rear wheel and different frame colours, why have we settle for less?



3 BACKREST



Standard backrest



Reclinable backrest



Half-folding backrest

2 REAR WHEEL SIZES

All backrest are available with big rear wheel (24") to make easier the displacement for the user and also with small rear wheel (12").



Navy Blue



Apple Green



Brilliant Red



Selenium Grey

COLOUR RANGE

4 frame colours in each model: navy blue, apple green, brilliant red, selenium grey.

COMFORT AND FUNCTIONALITY

EXTENSIVE EQUIPMENT

- Vented breathable backrest upholstery
- Side padded protector
- 70° removable and folding footrest
- Padded and justable armrest
- Ergonomic brake handle
- Rear wheel with spokes and aluminium rings

LESS WEIGHT

Breezy Style is the solution for a standard chair as light as possible. Thanks to its aluminium chassis, Breezy Style, has a **low transport weight** and **folds easily** to carry it anywhere.



Breezy Style with half-folding backrests, folded

DIFFERENT OPTION AVAILABLE

A large option range to be adapted to the user needs:

- Anatomic seats and backrest
- Adjustable armrest in height and depth
- Adjustable aluminium footplates
- Flip-up hangers
- Footboard
- Headrest
- Drum brakes
- IV pole
- Crutch holder
- Brake lever extension
- Double amputee kit
- Anti-tip tubes
- Transit wheels
- Tray
- Belt
- Back pack ...



STANDARD BACKREST

STANDARD CONFIGURATION

- Aluminium wheelchair with standard backrest
- Nylon seat upholstery
- Vented breathable backrest upholstery
- Detachable armrest and can be folded up with with the padded panels
- Padded side protector
- Standard footrest with detachable 70 degree hanger that swing outwards and inwards
- Cross wise spokes and Handrim Alu anodized silver
- Colours: Navy blue, green apple, brilliant red and selenium grey



Standard Configuration



**BREEZY STYLE WITH STANDARD
BACKREST AND 12" REAL WHEELS**



RECLINING BACKREST

STANDARD CONFIGURATION

- Aluminium wheelchair with reclining backrest till 60°
- Nylon upholstery
- T armrest, detachable in depth and height
- Standard footrest with detachable 70 degree hanger that swing outwards and inwards
- Cross wise spokes and Handrim Alu anodized silver
- Colours: Navy blue, green apple, brilliant red and selenium grey



Standard configuration with upholstery
recline headrest (optional)



BREEZY STYLE OPTIONAL CONFIGURATION

- 12" solid rear wheels
- Anatomic seat
- Anatomic backrest
- Skay headrest
- Elevating legrests
- Anti-tip tubes



HALF-FOLDING BACKREST

STANDARD CONFIGURATION

- Aluminium wheelchair with half-folding backrest for easier transportation and storage
- Nylon upholstery
- Detachable padded armrest, adjustable in depth and height
- Side padded protector
- Standard footrest with detachable 70 degree hanger that swing outwards and inwards
- Cross wise spokes and Handrim Alu anodized silver
- Colours: Navy blue, green apple, brilliant red and selenium grey



Standard configuration
(Optional Quick Release)



Breezy Style with half-folding
backrest, folded



**BREEZY STYLE WITH
HALF-FOLDED AND 12"
SOLID REAR WHEEL**



OPTIONS AND ACCESSORIES



LEGRESTS AND FOOTPLATES



Footboard



Aluminium angle adjustable footplates



Elevating legrest

ARMRESTS



Height adjustable armrest

BACKRESTS



Height adjustable push handles



Back pack

HEADRESTS



Skay headrest



Anatomic headrest



Fabric headrest for reclining backrest



ANATOMICAL OPTIONS AND BELTS



Anatomical backrest

Anatomical seat

Separating wedge

Belt

Anatomical seat with toilet bucket

Breezy Style with reclining backrest and anatomical seat and backrest

WHEELS



8" x 2" Pneumatic front castors



8" x 2" Pneumatic solid castors



6" Solid front wheels



12" Rear wheels (pneumatic or solid)



Different options in rear wheels:
22" or 24", with spokes, sport, pneumatic, solid.



Quick release axle on rear wheels

BRAKES



Brake lever extension



24" drum brakes

ACCESSORIES



One hand operation



Crutch holder



Anti-tip tubes



Tray



Doble amputee kit



Air pump



IV pole

Other options: Short footrest, transit wheels, oxygen bottle holder, spokes protectos, anti-folding system...



TECHNICAL SPECIFICATION

Seat width	38 / 40,5 / 43 / 45,5 / 48 / 52 cm	Total length reclining	116,5 cm (85,5 cm w/o footrest)
Seat depth	41,5 cm (45 w/reclining backrest)	Total length reclining 12" rear wheel	104 cm (75,5 cm w/o footrest)
Seat height	51,5 cm (49 ó 46,5 cm w/6" front wheels)	Total height	88 cm (132,5 cm w/reclining backrest and headrest)
Backrest height	45 cm (50 cm w/reclining backrest)	Seat angle	3°
Total width	Seat width + 18 cm (24" wheels) Seat width + 14,5 cm (12" wheels)	Backrest angle	10° (0° a 60° w/reclining backrest)
Total length w/24" rear wheels	103 cm (74,5 cm w/o footrest)	Colours	Green apple, brilliant red, navy blue and selenium grey
Total length w/12" rear wheels	98 cm (69,5 cm w/o footrest)	Approved user weight	125 kg
		Max. weight	From 13,9 kg (complete)*
		Folded width	26 cm

*Breezy Style Standard configuration, 38 cm size
Tolerance: +/- 1 cm

For further information on the full specification, options and accessories, please refer to the order form.
All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



Sunrise Medical S.L

Polígono Bakiola, 41
48498 Arrankudiaga
Vizcaya, España
Tel.: 902 14 24 34
Fax: 94 648 15 75



Ver 1. May 2018



www.SunriseMedical.eu