

Cirrus G5



Rev:11/2024
Valid:11/2025

ORDER FORM

Errors and imissions excepted.

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address	
C.-ID:			C.-ID:			
Name:			Name:			
Street:			Street:			
Town/PC:			Town/PC:			
E-Mail:			E-Mail:			
Phone:			Phone:			
Order Information			User Weight			
Order	Calculation	Built-4-Me®	User Weight Capacity: 135 kg			
Purchase Order No.:			User Weight: kg			
Art. No.: 35001000		Reimbursement Code: xxx				
Order Date:	Day	Month	Year	Base Price € X.XXX,-		

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat Width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether).
- b) On the wheelchair: From outside seat frame tube to outside seat frame tube.

Seat Depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee (approx. 30 - 60 mm). Care must taken with adj. upholstery.
- b) On the wheelchair: From back tube to front seat upholstery.

Front Seat Height:

- a) For the user: When sitting with 90° knee angle: from back of knee to bottom of shoe sole.
 - b) On the wheelchair: From front of seat tube to the ground. Min. 30 mm longer than user measurement.
- Remember to include the cushion height.

Rear Seat Height:

On the Wheelchair: From back of seat tube to the ground.

Note: Max. seat height difference rear to front is 110 mm

Back Height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.

Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably

Back angle:

Important: The back height depends on the degree of disability. 0° is the angle 90° to the ground. Negative angles are <90° and positive angles are >90°.

Centre of Gravity:

On the wheelchair: From the back tube to the middle of the axle plate.

Technical Data

Max. Overall Height:	1100 mm
Max. Overall Length:	1150 mm
Overall Width:	SW + 240 mm
Weight:	48 kg

Product is crash tested according to ISO 7176-19 and therefore approved for transportation in a vehicle in combination with crash test option. For transport without a user, fix the chair safely in the vehicle (e.g. car or plane).

Seat Width

	SW39	SW42	SW45	SW48	SW51	SW54
mm	390	420	450	480	510	540

Seat Depth

	SD44	SD47	SD50			
mm	440 (390-490)	470 (420-520)	500 (450-550)			

Front Seat Height (1)

	FSH375	FSH40	FSH425	FSH45	FSH475	FSH50	FSH525
mm	375	400	425	450	475	500	525
5" / 16" (4)			—	—	—	—	—
6" / 16" (4)	—	—				—	—
8" / 16" (4)	—	—	—			—	—
5" / 20"			—	—	—	—	—
6" / 20"	—	—					—
8" / 20"	—	—	—				
5" / 22"	—		—	—	—	—	—
6" / 22"	—	—					—
8" / 22"	—	—	—				
5" / 24"	—	—	—	—	—	—	—
6" / 24"	—	—					—
8" / 24"	—	—	—				

Seat (1)

€

CIR020086	Seat, tilting (tilt range 0° - 33°)	—
CIR020088	Locking mechanism tilt	XXX,-
CIR020089	Seat plate, fixed, +/- 5 cm seat depth	—
CIR020091	Seat plate, without pelvic hole (2)	—
CIR020090	Seat, sliding, +/- 5 cm seat depth (3)	XXX,-
CIR020092	Seat plate, with pelvic hole	—

(1) Seat height will be 25 mm higher with CIR020090 Seat, sliding, +/- 5 cm seat depth

(2) Not in combination with CIR020090 Seat, sliding, +/- 5 cm seat depth

(3) Not in combination with CIR020089 Seat plate, fixed, +/- 5 cm seat depth

(4) Only in combination with CIR070008 Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes

20" Rear wheel positioning		5", 6", 8" Front wheel configurations										25									
Fixed seat	Sliding seat	20" 5"		20" 5"		20" 5"		20" 6"		20" 6"		20" 6"		20" 8"		20" 8"		20" 8"		20" 8"	
375	400	a b c		a b c		a b c		a b c		a b c		a b c		a b c		a b c		a b c		a b c	
400	425	1 2 3		1 2 3		1 2 3															
425	450					3 3 3															
450	475					4 4 4		4 4 4				4 4 4									
475	500							5 5 5		5 5 5		5 5 5		5 5 5		5 5 5					
500	525									5 5 5						5 5 5		5 5 5		5 5 5	
525	550																	5 5 5		5 5 5	

22" Rear wheel positioning			7", 6", 5" Front wheel configurations										26	
														
	Fixed seat	Sliding seat												
			22" 7"	22" 6"	22" 5"	22" 7"	22" 6"	22" 5"	22" 7"	22" 6"	22" 5"	22" 7"	22" 6"	22" 5"
			a b c	a b c	a b c	a b c	a b c	a b c	a b c	a b c	a b c	a b c	a b c	a b c
	401	495	1 1 1											
475	475		2 2 2	2 2 2										
450	475				3 3 3	3 3 3								
475	500				4 4 4	4 4 4								
500	525						5 5 5	5 5 5						
525	550								5 5 5	5 5 5				

[illegible]

16" Rear wheel positioning			5", 6", 8" Front wheel configuration												28											
			a b c			a b c			a b c			a b c			a b c			a b c			a b c					
			375			400			375			400			375			400			375			400		
			400			425			400			425			400			425			400			425		
			425			450			425			450			425			450			425			450		
			450			475			450			475			450			475			450			475		
475			500			475			500			475			500			475			500					

Seat Cushions		€
CIR020093	Seat cushion, Trevira	XXX,-
CIR020094	Seat cushion, Nylon, PU coating	–
CIR020110	Omit seat cushion	–
SEDEO PRO Seat Cushions		
CIR020073	SEDEO PRO, seat cushion, profiled, height 70 mm (6) (8)	XXX,-
CIR020075	SEDEO PRO, seat cushion, support, height 90 mm (6)	XXX,-
CIR020078	SEDEO PRO, seat cushion, stabilizer, height 90 mm (6)	XXX,-
CIR020076	SEDEO PRO, seat cushion, support memory, height 90 mm (7)	XXX,-
SEDEO PRO Seat Cushion Cover		
CIR020518	Seat cushion cover, incontinence protection	XXX,-
CIR020080	Seat cushion cover, comfort	XXX,-
JAY Cushions		€
JAY020006	JAY Easy Visco F (1)	XXX,-
JAY020009	JAY Ion (1) (2)	XXX,-
JAY020024	JAY Union (2) (3)	XXX,-
JAY020015	JAY Balance, Fluid (4)	XXX,-
JAY020020	JAY J2 Deep Contour (5)	XXX,-

(1) Medium skin protection
(2) UK only
(3) Medium-high skin protection
(4) High skin protection

(5) Very high skin protection, only for SW 390/450/510 mm
(6) Not with SD44
(7) SD50 only
(8) Only in combination with Seat, sliding CIR020090

Seat cushions & Seat plate



// CIR020093



// CIR020073



// CIR020075



// CIR020076



// CIR020078



// JAY020006



// JAY020020




// JAY020009



// JAY020024



// CIR20092



// CIR20091

Seat Accessories			€	
CIR020070	Incontinence cover (1)		XXX,-	
CIR090202	Abduction pommel, small, 70 x 60 x 100 mm		XXX,-	
CIR090302	Abduction pommel, large, 90 x 90 x 100 mm		XXX,-	
Positioning Belt with lock				
CIR090135	Positioning belt, carseat buckle		XXX,-	
CIR090156	Upholstery positioning belt, small		XXX,-	
CIR090136	Upholstery positioning belt, large		XXX,-	
Positioning Belts Kinetic Balance (Sweden only)				
CIR090137	Positioning belt Kinetic Balance, XS		XXX,-	
CIR090144	Positioning belt Kinetic Balance, S		XXX,-	
CIR090150	Positioning belt Kinetic Balance, M		XXX,-	
CIR090151	Positioning belt Kinetic Balance, L		XXX,-	
CIR090152	Positioning belt Kinetic Balance, XL		XXX,-	
CIR090138	Positioning belt Kinetic Balance, 4 point fixation, XS		XXX,-	
CIR090153	Positioning belt Kinetic Balance, 4 point fixation, S		XXX,-	
CIR090154	Positioning belt Kinetic Balance, 4 point fixation, M		XXX,-	
CIR090155	Positioning belt Kinetic Balance, 4 point fixation, L		XXX,-	
JAY Pelvic Belts				
JAY090071	2 Point Anterior Pelvic Belt 40 mm Push Button - Single Pull 190 mm Pad		XXX,-	
JAY090072	2 Point Anterior Pelvic Belt 40 mm Push Button - Single Pull 240 mm Pad		XXX,-	
JAY090077	JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 140 x 40 mm		XXX,-	
JAY090078	JAY Medial 4-Point Pelvic Belt, Push Button Buckle, Pad 165 x 40 mm		XXX,-	
JAY090102	JAY V-Contour Pelvic Belt, Side Release Buckle, Pad, Medium		XXX,-	
JAY090103	JAY V-Contour Pelvic Belt, Side Release Buckle, Pad, Large		XXX,-	
Backrest				€
Backrests				
CIR030123	Backrest recline (0° - 35°), gas strut		—	
	BH55	Back height 550 mm	—	
	BH65	Back height 650 mm	XXX,-	
CIR030134	Backrest, reclining, 0° - 35°, gas strut, for JAY, SEDEO and solid backs		—	
CIR030124	Locking mechanism recline		XXX,-	
Push Handles				
CIR030126	Push bar		—	
CIR030070	Push bar, angle-adjustable (2)		XXX,-	
CIR030204	Push handles, height-adjustable		XXX,-	
Back Cushions				
CIR030127	Back cushion, Trevira		XXX,-	
CIR030128	Back cushion, Nylon, PU coating		—	
CIR030129	Back cushion, Unistabil (3)		XXX,-	
CIR030209	Omit back cushion		—	
SEDEO LITE Back Cushions (4)				
CIR030130	SEDEO LITE, back cushion, profiled, height 500 mm, comfort cover		XXX,-	
CIR030131	SEDEO LITE, back cushion, profiled, adjustable straps, height 570 mm, comfort cover		XXX,-	
CIR030132	SEDEO LITE, back cushion, deep contour, passive, adjustable straps, height 570 mm, comfort cover		XXX,-	
CIR030133	SEDEO LITE, back cushion, deep contour, active, adjustable straps, height 570 mm, comfort cover		XXX,-	

(1) Not with SEDEO & JAY seat cushions

(2) Not in combination with CIR030134 Backrest, reclining, 0° - 35°, gas strut, for JAY and solid backs

(3) Only available height 550 mm, only with CIR030134 Backrest, reclining, 0° - 35°, gas strut, for JAY, SEDEO LITE and solid backs; Sweden only

(4) Only with CIR030134 Backrest, reclining, 0° - 35°, gas strut, for JAY, SEDEO LITE and solid backs

Seat accessories



// CIR090137/144/150/151/152



// JAY090071



// CIR090302/202



//JAY090077



// JAY090102

Backrest



// CIR30127



// CIR030126



// CIR030070



// CIR030204

Headrest			€
CIR030064	Headrest Standard, small, 139 x 210 x 40 mm		XXX,-
	CIR030137	Receiver, vertical and horizontal angle adjustable	XXX,-
	CIR030138	Cover, padded	XXX,-
CIR030066	Headrest, lateral guidance, small, 140 x 230 x 70		XXX,-
	CIR030137	Receiver, vertical and horizontal angle adjustable	XXX,-
	CIR030138	Cover, padded	XXX,-
CIR030135	Head-/neckrest, comfort, tension adjustable, Trevira		XXX,-
	CIR030137	Receiver, vertical and horizontal angle adjustable	XXX,-
CIR030136	Head-/neckrest, comfort, tension adjustable, PU coated		XXX,-
	CIR030137	Receiver, vertical and horizontal angle adjustable	XXX,-

Whitmyer Headrests

WHT030406	Whitmyer, Headrest Pad, Adjust-A-Plush, Large, (LxH) 460 x 127 mm (1)		XXX,-
WHT030412	Whitmyer, Headrest Pad, Contoured Cradle, Medium (LxH) 250 x 158 mm (1)		XXX,-
WHT030413	Whitmyer, Headrest Pad, Contoured Cradle, Large (LxH) 330 x 190 mm (1)		XXX,-
WHT030414	Whitmyer, Headrest Mount, Cirrus G5, vertical and angle adjustable		XXX,-

JAY J3 Backs (For complete product range see JAY order forms)

JAY J3 Back, Shallow Contour (SC)	Short (S)		Medium (M)		Tall (T)		
Upper Thoracic (UT)	JAY030012	420 mm	JAY030019	460 mm	JAY030025	500 mm	XXX,-
Shoulder Height (SH)	JAY030015	530 mm	JAY030022	570 mm	JAY030028	610 mm	XXX,-
JAY J3 Back, Mid Deep Contour (MDC)	Short (S)		Medium (M)		Tall (T)		
Upper Thoracic (UT)			JAY030042	460 mm			XXX,-
Shoulder Height (SH)	JAY030043	530 mm					XXX,-
JAY J3 Back, Deep Contour (DC)	Short (S)		Medium (M)		Tall (T)		
Upper Thoracic (UT)	JAY030014	420 mm	JAY030021	460 mm	JAY030027	500 mm	XXX,-
Shoulder Height (SH)	JAY030017	530 mm	JAY030024	570 mm	JAY030030	610 mm	XXX,-

JAY J3 Back Headrests

JAY030004	JAY Flip Back Headrest		XXX,-
JAY030006	JAY Std. Headrest		XXX,-
WHT030415	Whitmyer Headrest, COBRA XTRA, flip back (ONLY WITH JAY J3 BACK)		XXX,-
WHT030416	Whitmyer Headrest, AXYS, Standard (ONLY WITH JAY J3 BACK)		XXX,-

JAY J3 Lateral Thoracic Pads

JAY090024	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 127 mm, Right		XXX,-
JAY090015	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 127 mm, Left		XXX,-
JAY090025	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 152 mm, Right		XXX,-
JAY090016	JAY J3 Lateral Thoracic Pad, (LxH) 152 mm x 152 mm, Left		XXX,-

JAY J3 Lateral Thoracic Brackets

JAY090028	Lateral Thoracic, Swing Away Bracket, Right		XXX,-
JAY090019	Lateral Thoracic, Swing Away Bracket, Left		XXX,-

Back Accessories

CIR030139	Positioning wedge kit		XXX,-
CIR030143	Lateral, thoracic support, height- and depth-adjustable, swing away, right		XXX,-
CIR030142	Lateral, thoracic support, height- and depth-adjustable, swing away, left		XXX,-
CIR030144	Cover, lateral thoracic support, padded (2)		XXX,-
JAY090087	Classic Style Anterior Trunk Support - Dynamic, Medium		XXX,-
JAY090088	Classic Style Anterior Trunk Support - Dynamic, Large		XXX,-
JAY090104	Contour Style Anterior Trunk Support - Dynamic, Medium		XXX,-
JAY090105	Contour Style Anterior Trunk Support - Dynamic, Large		XXX,-

(1) Including WHT030414 Whitmyer, Headrest Mount, Cirrus G5, vertical and angle adjustable

(2) Not in combination with SEDEO LITE backs

Headrests



// CIR030064



// CIR030066



// CIR30135



// WHT030412



// WHT030406



// WHT030414

JAY J3 Backs



// Shallow Contour (SC)



// Deep Contour (DC)



// Mid Deep Contour (MDC)



Back accessories



// CIR030143



// JAY090087



// JAY090104

Electric back angle and seat tilt adjustment €

CIR020071	Seat tilt adjustment, powered	XXX,-
CIR030145	Back angle adjustment, powered	XXX,-
CIR020072	Lock for powered adjustment	XXX,-

Armrests €

CIR040042	Single-post sideguard, armpad (400 mm), tool height-adjustable, removable	–
	CIR040031 Armpad hemiplegia, integral foam, right	XXX,-
	CIR040032 Armpad hemiplegia, integral foam, left	XXX,-
	CIR020095 Seat width extension (+ 30 mm)	XXX,-
CIR040043	Single-post sideguard, armpad (440 mm), height-adjustable, reclining, flip back, removable (1)	–
	CIR040031 Arm pad hemiplegia, integral foam, right	XXX,-
	CIR040032 Arm pad hemiplegia, integral foam, left	XXX,-
CIR040000	Without sideguard	XXX,-

Armrest accessories

CIR040044	Cover arm pad, padded, PU (pair)	XXX,-
CIR040045	Cover, side panel, half-padded, PU (pair) (2)	XXX,-
CIR040046	Cover, side panel, complete padded, PU (pair) (2)	XXX,-

Tray tables

CIR090131	Tray table, transparent, slide on, width-adjustable 390 - 540 mm	XXX,-
	CIR090125 Partial padding for tray table (3)	XXX,-
CIR090139	Tray table, Half-Side, slide on, right (4)	XXX,-
CIR090140	Tray table, Half-Side, slide on, left (5)	XXX,-

(1) Not with CIR030134 Backrest, reclining, 0° - 35°, gas strut, for JAY, SEDEO LITE & solid backs

(2) Only with Single-post sideguard, armpad (440 mm), height-adj., reclining, flip back, removable

(3) Not with Arm pad hemiplegia, integral foam

(4) Not with Arm pad hemiplegia, integral foam, right, Optioncode CIR040031

(5) Not with Arm pad hemiplegia, integral foam, left, Optioncode CIR040032

Footrest & Footplates

€

Hanger aluminium, angle adjustable (70°/80°/90°)

CIR050117	Hanger, aluminium, tool angle-adjustable 70°/80°/90°, swing-away, right	–
CIR050118	Hanger, aluminium, tool angle-adjustable 70°/80°/90°, swing-away, left	–

Hanger aluminium, flip-up (0° - 90°)

CIR050010	Hanger, aluminium, flip-up (0° - 90°), swing away, right	–
CIR050009	Hanger, aluminium, flip-up (0° - 90°), swing away, left	–

Amputee support

CIR050013	Amputee support, padded, right	XXX,-
CIR050014	Amputee support, padded, left	XXX,-

Cast support

CIR050094	Cast support, right	XXX,-
CIR050093	Cast support, left	XXX,-

Omit hanger and footplate

CIR050000	Omit hanger and footplate	XXX,-
-----------	---------------------------	-------

Footplates

CIR050119	Footboard, divided, composite, angle-adjustable, flip-up (sideways), black	–
CIR050120	Footboard, divided, composite, depth- and angle-adjustable, flip-up (sideways), black	XXX,-
CIR050121	Footboard, platform, aluminium, angle-adjustable, black (1)	XXX,-

Footbox (1)

CIR050122	Footbox, padded, PU coated	XXX,-
-----------	----------------------------	-------

Footrest accessories

CIR050144	Clinic kit for lower leg length adjustment, without tools	XXX,-
CIR050145	Padding for footboard, PU coated (1)	XXX,-
CIR050123	Calfpad, padded, not split (1)	XXX,-
CIR050124	Knee pad, PU coated (1)	XXX,-
CIR050134	Universal padding, PU coated (2)	XXX,-
CIR050133	Footplate extension, aluminium, black (3)	XXX,-

(1) Only for CIR050009-10 Hanger, aluminium, flip-up (0° - 90°), swing away left/right

(2) E.g. for footplate, amputation support, leg support

(3) Only CIR050119 Footboard, divided, composite, angle-adjustable, flip-up (sideways), black & CIR050120 Footboard, divided, composite, depth-and angle-adjustable, flip-up (sideways), black

Armrests & Accessories



// CIR090131



// CIR090140



// CIR040044



// CIR040045



// CIR040046



// CIR040042



// CIR040031

Footrests & Footplates



// CIR050117



// CIR050010



// CIR050013



// CIR050094



// CIR050119



// CIR050122



// CIR050145



// CIR050124



// CIR050134



// CIR050123

Center Of Gravity (COG)		€
COG040	-40 mm (Safe)	-
COG00	0 mm (Standard)	-
COG40	40 mm (Active)	-

Frame accessories		€
Tip assist (step tube)		
CIR090003	Tip assist, pair	-
Anti-tip tubes swing away		
CIR090006	Anti-tip tubes, swing away, pair	-
CIR090060	Crutch holder, right	XXX,-
CIR090039	Crutch holder, left	XXX,-

Castor Wheels		€
CIR080301	Castor, solid (hard), rim: composite, black	-
CIR080305	Castor, pneumatic, rim: composite, black (1)	-

Rear Wheels		€
CIR070073	Rear wheel, 36 spokes (silver), cross-wise spoked (2)	-
CIR070008	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes (7)	XXX,-

Tyres		€
CIR070076	Tyre, 1", pneumatic, tread (3)	-
CIR070077	Tyre, 1 3/8", pneumatic, tread (4)	-
CIR070074	Tyre, 1", puncture proof, tread (3)	XXX,-
CIR070075	Tyre, 1 3/8", puncture proof, tread	XXX,-

Handrims		€
CIR070202	Handrims, aluminium, anodised, silver (4)	-
CIR070204	Handrims, stainless steel, silver (4)	XXX,-
CIR070207	Handrims, Supergrip® (4)	XXX,-
CIR070208	Handrims, Maxgripp® (4) (5)	XXX,-
CIR070251	Handrim, cover, silicon (4) (6)	XXX,-

(1) Only for 6" & 8" castors
(2) Only for 20", 22" and 24" rear wheels
(3) Not with drumbrake, 16" wheel
(4) Not with 16" wheel

(5) Not with 20" wheel
(6) Not with CIR070207 Handrims, Supergrip® and CIR070208 Handrims, Maxgripp®
(7) 16" rear wheel with composite rim

Frame accessories



// CIR090003



// CIR090006



// CIR090060

Rear wheels & Castors



// CIR070008



// CIR070073



// CIR080301



// CIR070208



// CIR070207

Spoke Guards (1) €

CIR070831	Spoke guards, transparent	XXX,-
-----------	---------------------------	-------

Wheel locks (Choose only one option) €

CIR060009	Brake, pull to lock, brake lever long	-
-----------	---------------------------------------	---

One arm wheel lock

CIR060021	Brake, pull to lock, one-arm, right actuation	XXX,-
CIR060022	Brake, pull to lock, one-arm, left actuation	XXX,-

Accessories €

CIR090000	Tool kit	XXX,-
CIR090149	Air pump, powered	XXX,-
CIR090025	Air pump	XXX,-
CIR090030	Back pack	XXX,-
CIR090083	IV support with multiposition clamp	XXX,-
CIR090107	Oxygen bottle holder	XXX,-
CIR090148	Universal headrest fixation	XXX,-

Colours €

CIR100134	Matt-Anthracite	-
-----------	-----------------	---

Wheel locks



// CIR060009



// CIR060021

Accessories



// CIR090148



// CIR090083



// CIR090107

Comments for Adaptations



Cirrus G5	€
-----------	---

Base price	
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total price	

Date	Signature
------	-----------

Day Month Year

**Sunrise Medical HCM**

Phone: +31 (0) 492 593 888

Fax: +31 (0) 492 546 497

E-Mail: customerservice@sunrisemedical.nl

www.SunriseMedical.eu