

EASY 200 / 300

The Benchmark in active

mobility for many years...



Its high level of flexibility and individual adjustment account for the success of the Easy 300 and today, more than ever before, these features take centre stage.

Easy 300 Swing away version

A total concept, which will convince you again and again: three frame lengths, the very best in adjustments and compact construction – the Easy 200 is the classic in folding wheelchairs



Fixed Front version

Highly Configurable



Excellent range of improved sideguards, backrests, wheels, footrests, **swing-away legrests (Easy 300)** and designed to meet all users needs. The new selection of upholsteries and options gives the Easy 200 / 300 a fantastic look.

Easy handling

Easy to fold with retaining strap and **quick release legrest (Easy 300).**Good driving characteristics with centre of gravity adjustment and rear wheel camber option.



New selection of back upholsteries with a wide range of different versions and an improved adjustability



New sideguard adapter for improved locking and sideguard flip up function



New des-sideguard with height adjustbale armrest



Flip up and swing away footboard offers larger space for transfers (optional)



Elevating legrest with angle adjustable footrests for comfortable seating



Folding push-handles with a wide range of push handle options available

Easy 200 / 300




















Technical data sheet

28 individual frame colours 3 coloured anodized packages (silver, black and orange)		Maximum user weight Weight of heaviest part	140 kg 7.60 kg

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.





PB_EASY200_EN_EU_REV1 PB_EASY300_EN_EU_REV1

Sunrise Medical GmbH & Co. KG

Kahlbachring 2-4 69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0 Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de

www.SunriseMedical.eu