

Order form

EMINEO K

Valid from	01-01-2021	Version	1.0
Dealer		Date	
Reference No.			
Delivery address		Quantity	
		PRICE LIST ID	0249

- Seat depth 37 cm (possible adjustments 34 till 40 cm)
- Seat height : 39,7 cm
- Frame colour : grey
- Brake : Pull to lock, foldable
- Armrest : height adjustabel with side guards
- Velcro back : angle adjustable +7 (97°)
Adjustable from 0° to +20° (90° - 110°)
Back height 55cm, stepless adjsutable 50 - 60 cm
- Push bar
- Anti tip RHS









The standard configuration below is indicated with 'basic'. If an alternative option is required please mark clearly by ticking the associated box.

THE CODES WHICH ARE MENTIONED BEHIND THE OPTIONS, ARE NOT ARTICLE NUMBERS. ALWAYS USE THESE CODES IN COMBINATION WITH THE PRICE LIST-ID NUMBER MENTIONED ABOVE.



CHASSIS

		Option	Excl.
	<input type="checkbox"/> 1 Emineo KID SW 33	Y01053	Basic
	<input type="checkbox"/> 1 Emineo KID SW 36	Y01054	0

WHEELS

		Option	Excl.
	<input type="checkbox"/> 1 One tool Wheelblock and Castorhousing	Y01570	Basic
	<input type="checkbox"/> 2 Basic Wheelblock and Castorhousing	Y01569	0
	<input type="checkbox"/> 3 Front wheel solid	Y00829	Basic
	<input type="checkbox"/> 3 Front wheels 4" (only i.c.w. 20" rear wheels)	Y00882	0
	<input type="checkbox"/> 3 Front wheel 5"	Y00827	Basic
	<input type="checkbox"/> 4 Rear wheel 20"	Y00824	0
	<input type="checkbox"/> 4 Rear wheel 22"	Y00825	Basic
	<input type="checkbox"/> 4 Rear wheels solid 37 mm	Y00149	65
	<input type="checkbox"/> 4 Rear wheels high pressure 25 mm	Y00880	Basic
	<input type="checkbox"/> 4 Rear wheels solid 25 mm	Y00881	65
	<input type="checkbox"/> 4 Pushrim aluminium	Y00685	Basic
	<input type="checkbox"/> 5 Pushrim grip-friction	Y00686	111
	<input type="checkbox"/> 6 Silicon cover for pushrim, set	Y00021	97
	<input type="checkbox"/> 7 Rear wheels solid with drumbrake	Y00830	136
	<input type="checkbox"/> 7 Rear wheels pneumatic with drumbrakes	Y00150	108
	<input type="checkbox"/> 8 12,5" solid rear wheels	Y00135	99



BACKRESTS

		Option	Excl.
	<input type="checkbox"/> 1 Velcro back, angle adj. +7 (97°) Adjustable from 0° to +20° (90°-110°) Back height 55 cm, stepless adjustable 50-60cm	Y00625	Basic
	<input type="checkbox"/> Backrest adjustment with gas spring	Y00387	98
	<input type="checkbox"/> 2 Positioning Cushion, Small, 15 mm	Y00681	36
	<input type="checkbox"/> 2 Positioning Cushion, Large, 30 mm	Y00682	36



The product picture shown may deviate from the standard configuration.

BACKRESTS , CONTINUATION

			Option	Excl.
		<input type="checkbox"/> 3 Push bar, height 55 cm, stepless adjustable 50-60 cm	Y01115	Basic
		<input type="checkbox"/> 4 Height adjustable push handles, set	Y00124	67

SEAT

	Option	Excl.
<input type="checkbox"/> Seatplate standard	Y00709	Basic
<input type="checkbox"/> Seat tilt adj. handle on armrest, right	Y00887	0






BRAKES

	Option	Excl.
	<input type="checkbox"/> 1 Pull to lock, foldable	Y00602 Basic





















ARMRESTS

	Option	Excl.
	<input type="checkbox"/> 1 Armrest height adjustable 23 - 33 cm, right	Y00075 Basic
	<input type="checkbox"/> 1 Armrest height adjustable 23-33 cm, left	Y00062 Basic

LEGRESTS

	Option	Excl.
	<input type="checkbox"/> 1 Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, right	Y00741 43
	<input type="checkbox"/> 1 Legrest hanger, adjustable 70°, 80° or 90°, set up 80°, left	Y00740 43
	<input type="checkbox"/> 2 Angle adjustable legrest hanger, right	Y00737 68
	<input type="checkbox"/> 2 Angle adjustable legrest hanger, left	Y00736 68
	<input type="checkbox"/> 3 Calf support SW 36-39, right	Y00640 32
	<input type="checkbox"/> 3 Calf support SW 36-39, left	Y00636 32
	<input type="checkbox"/> 4 Footplate angle adjustable without heelstrap, right	Y00890 32
	<input type="checkbox"/> 4 Footplate angle adjustable without heelstrap, left	Y00889 32
	<input type="checkbox"/> 5 Footplate angle adjustable with heelstrap, right	Y00892 39
	<input type="checkbox"/> 5 Footplate angle adjustable with heelstrap, left	Y00891 39

ACCESSORIES

	Option	Excl.
	<input type="checkbox"/> 1 Calf strap	Y00635 11
	<input type="checkbox"/> 2 Anti-tip stabiliser with Step-tip for 20" and 22" rear wheels, right	Y00868 Basic
	<input type="checkbox"/> 2 Anti-tip stabiliser with Step-tip for 20" and 22" rear wheels, left	Y00867 48
	<input type="checkbox"/> 3 Anti-tip stabiliser swing forward with leg, right	Y00614 58
	<input type="checkbox"/> 3 Anti-tip stabiliser swing forward with leg, left	Y00612 58
	<input type="checkbox"/> 4 Crutchholder Basic	Y00651 21
	<input type="checkbox"/> Front pouch/satchel	Y00019 21
	<input type="checkbox"/> 5 Rucksack	Y00006 32
	<input type="checkbox"/> 6 Abduction pommel, small	Y00540 80
	<input type="checkbox"/> 7 Side support, right	Y01434 61
	<input type="checkbox"/> 7 Side support, left	Y01435 61
	<input type="checkbox"/> Cover Side support, small	Y00650 42
	<input type="checkbox"/> 8 Positioning belt with plastic buckle	Y01330 12
	<input type="checkbox"/> 9 Positioning belt with steel buckle	Y01331 14
	<input type="checkbox"/> 10 Padding for positioning belt, small	Y01332 19
	<input type="checkbox"/> 11 Padding for positioning belt, large	Y01333 21
	<input type="checkbox"/> 12 Amputation axle plates	Y00082 63
	<input type="checkbox"/> 13 Spoke protectors 20"	Y00032 87
	<input type="checkbox"/> 13 Spoke protectors 22"	Y00033 87
	<input type="checkbox"/> 14 Headrest Basic	Y00665 68
	<input type="checkbox"/> 15 Headrest Basic, Advanced pad	Y00666 141
	<input type="checkbox"/> 16 Headrest Advanced	Y00667 178
	<input type="checkbox"/> 17 Headrest Advanced, Basic Pad	Y00668 141
	<input type="checkbox"/> 18 Cover headrest Basic pad	Y00648 38
	<input type="checkbox"/> 18 Cover headrest Advanced pad	Y00649 58
	<input type="checkbox"/> 19 Electrical inflator	Y00609 81
	<input type="checkbox"/> 20 Seat cushion	Y00476 47
	<input type="checkbox"/> Seat cushion depth 34	Y01055 0
	<input type="checkbox"/> Seat cushion depth 40	Y01056 0

SEAT HEIGHT MATRIX

(**seat cushion excluded)

Rear wheel	Seat height back **	Front wheel	Fork	Seat height front basic **	Comments
20"	38 (adjustable 36-44)	4"	92 mm	38	
20"	38 (adjustable 36-44)	5"	92 mm	38 (38-42)	
22"	40 (adjustable 38-45)	5"	92 mm	40 (40-44)	

COMMENTS

Sunrise Medical HCM can tailor solutions to meet your requirements. Please contact your Sunrise Medical HCM contact person for all additional wishes that are not described above.