

ATTITUDE Hybrid



Rev:11/2024
Valid:11/2025

ORDER FORM

Errors and imissions excepted.

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address
C.-ID:			C.-ID:		
Name:			Name:		
Street:			Street:		
Town/PC:			Town / PC:		
E-Mail:			E-Mail:		
Phone:			Phone:		
Order Information			User Weight		
Order	Calculation	Built-4-Me®	Max. user weight: 100 kg		
Purchase Order No.:			User Weight: kg		
Art. No.: 47500000 / 47500001 IC			Reimbursement Code: xxx		
Order Date:	Day	Month	Year		

Base Price € X.XXX,-

Attitude - Model Selection

€

ATT010802	ATTITUDE Hybrid 10 speed cassette and derailleur; trigger shifter; 17 Ah Lithium-ion battery; motor support up to max 25 kph; front and rear lamp; docking system; stand with wheels, handle; v-crank, disk brake; parking brake; speedometer integrated in control display.	XXX,-
ATT072102	Connect Display - coloured Display, includes Smartphone connectivity via bluetooth	XXX,-
FMG-MJ02-3187-02	B4Me option: Tetraplegia operation for Neodrives power unit Build 4 me - different lead-times and discount structures	—

Technical Data Overview

Max. user weight	100 kg	Width of V-crank	470 mm
Front frame width	260 – 460 mm ⁽¹⁾	Length of the crank	170 mm
Total weight	21,5 kg	Brake	disc brake (180 mm)
Wheel size	20"	Parking brake	yes
Shifter	thumb shifter	Frame colors	24
Gears	10		

Wheelchair Model - which will be used in combination with the ATTITUDE

€

Sunrise Wheelchairs		Other Wheelchairs	
Easy Life T (fixed front frame only)	—	Küschall K-serie; KSL Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's Fits Kueschall round-tube frames until build date 03/ 2020	XXX,-
Easy Life RT (fixed front frame only)	—	Küschall Advance, Champion Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's Fits Kueschall round-tube frames until build date 03/ 2020	XXX,-
Easy Life (fixed front frame only)	—	Otto Bock Avantgarde CS/CLT Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's Only possible for frames without inset	XXX,-
Easy Life R (fixed front frame only)	—	Otto Bock Motus CS Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's Only possible for frames without inset	XXX,-
Easy 200	—	Meyra Avanti Pro Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-
Easy Max (fixed front frame only)	—	Meyra Hurricane Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-
Neon (fixed front frame only)	—	Meyra Nano Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-
Xenon/Xenon ² (fixed front frame only)	—	Meyra Smart F Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-
Argon ²	—	PRO ACTIV Traveler Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-
Helium Helium with handbike axle is only approved with frame reinforcement kit for add-on bikes	—	PRO ACTIV Speedy Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-
Nitrum Nitrum with handbike axle is only approved with frame reinforcement kit for add-on bikes	—	Progeo Joker Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-
RGK Tiga Please note: V-Frames are only suitable up to 5°; please check with your authorized dealer if there are any further NoGo's	XXX,-	TiLite TR, TRA	XXX,-
QS5 X (fixed front frame only)	-	TiLite Aero T, Z, ZR, ZRA	XXX,-
		Panthera S2, S3, U2, U3, Build 4 me - different lead-times and discount structures	XXX,-
		Panthera U2 Light, U3 Light Build 4 me - different lead-times and discount structures	XXX,-

Attitude manual adjustment and set-up


A1T080022




// ATTITUDE Manual with Standard fork: crank position range (adjustment via spacers in the upper docking clamp)

Attitude manual fork: Angle and depth adjustability

A1T080035



// ATTITUDE Manual / Hybrid - crank position range with option A1T080035 "Fork, angle and depth-adjustability of crank-position"



// **Depth adjustability:** 100mm (+/- 50mm compared to standard fork) continuously adjustable

Angle adjustment: from 0° (parallel to the ground) to 28° continuously adjustable

B4ME-Options

FVG-MA03-3201-01



// To move the pedal position 2 cm closer or up to 6 cm further

Order process information

ATT090213	ATTITUDE only - without a related wheelchair order	
ATT090214	ATTITUDE with wheelchair order - on a separate wheelchair order form	
FMG- MA03-3201- 01	B4Me option: Adaptation link -2 to + 6 cm Build 4 me - different lead-times and discount structures	

Wheelchair front frame width - measured from outside front frame tube to outside front frame tube

€

260 mm	380 mm
280 mm	400 mm
300 mm	420 mm
320 mm	440 mm
340 mm	460mm
360 mm	

Fork

€

ATT080034	Fork, Standard	–
ATT080035	Fork, angle and depth-adjustability of crank-position	XXX,-

Tyres

€

ATT070104	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	–
ATT070110	Tyre, pneumatic, tread, Mountainbike	–

Gear and shifters

€

ATT070045	10 speed cassette with derailleur for Hybrid, Shifter, trigger	–
FMG- MJ02-3187-01	B4Me option: Tetraplegia operation for universal shifter Build4me - different lead-times and discount structures; only in combination with ATT070044	–

Crank and handle - special length, widths and shapes via B4Me

€

ATT080022	V-crank, standard (Length: 170 mm, Width: 470 mm)	–
ATT080023	Handle, vertical, 10° angle, round, diameter 30 mm	–
	ATT080027 Quad, hand support	XXX,-
	ATT080028 Quad, arm support	XXX,-
ATT080024	Handle, vertical, 10° angle, round, diameter 40 mm	–
ATT080025	Handle, vertical, 25° angle, oval, diameter 37 x 21 mm	XXX,-
ATT080026	Handle, vertical, 25° angle, oval, diameter 47 x 32 mm	XXX,-

Disk brake - disc-diameter 180 mm, parking brake: V-brake

€

ATT060050	Brake, manual actuation not with gear hub, 8 speed, coaster brake	–
ATT060053	Tyre, pneumatic, tread, Mountainbike	XXX,-
	ATT060102 Brake lever mounted to right handle	–
	ATT060101 Brake lever mounted to left handle	–
ATT060052	2 Driving brakes plus parking brake parking brake as rim brake and not as disk brake, not with gear hub, 8 speed, coaster brake	XXX,-
ATT060051	Brake, backpedal actuation not with ATT080035 Fork, angle and depth-adjustability of crank-position, not with gear hub, 8 speed, coaster brake	XXX,-

Fork

A1 T080034



// Fork for Attitude Manual, angle and depth-adjustability of crank-position

Tyres

A1 T070104



// Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution

A1 T070110



// Tyre, pneumatic, tread, Mountain-bike

Shifter and handel

A1T070045



// Shifter, trigger

FM G-MJ02-3187-01



// Tetra operation

A1T080023

// Tyre, pneuma Handle, vertical,
10° angle, round, diameter 30 mm
tic, tread, Big Apple

A1T080027



// Quad, hand support

ATT080028



// Quad, arm support

A1T080025



// Handle, vertical, 25° angle, oval

Disk brake

ATT060050



// Performance brake lever

Battery - charger always included

€

ATT160020	Battery, Lithium-ion, 17 Ah	–
ATT090088	Additional battery, Lithium-ion, 17 Ah	–
ATT090091	Additional charger, Hybrid	

Accessories

€

ATT090074	Front- & rear lamp	–
ATT090080	Speedometer	–
ATT090095	Mudguard	XXX,-
ATT090067	Rear view mirror	XXX,-
ATT090081	Bottle holder	XXX,-
ATT090082	Low rider carrier for saddle bags	XXX,-
ATT090024	High pressure pump (0 - 11 bar)	XXX,-
ATT090084	Additional weight, removable, 5 kg, to improve traction	XXX,-
ATT090085	Bell	XXX,-
ATT090041	Multitool	XXX,-
ATT090096	Protection foil	XXX,-

Chain guard Colours You can choose separate colours for the frame, the cranks and the chain guard.

	Frame	Chain guard	Crank		Frame	Chain guard	Crank
02 Brilliant Silver	ATT100102	ATT100802	ATT100902	29 Black	ATT100129	ATT100829	ATT100829
03 Silver Matt	ATT100103	ATT100803	ATT100903	31 Matt Black	ATT100131	ATT100831	ATT100931
06 Sun Yellow	ATT100106	ATT100806	ATT100906	32 Signal White	ATT100132	ATT100832	ATT100932
14 Blue	ATT100114	ATT100814	ATT100914	33 Pink	ATT100133	ATT100833	ATT100933
15 Brilliant Marine	ATT100115	ATT100815	ATT100915	34 Dark Chocolate Metallic	ATT100134	ATT100834	ATT100934
16 Matt Marine	ATT100116	ATT100816	ATT100916	35 Skyblue	ATT100135	ATT100835	ATT100935
19 Brilliant Purple	ATT100119	ATT100819	ATT100919	36 Cyan-Metallic	ATT100136	ATT100836	ATT100936
21 Red Ruby	ATT100121	ATT100821	ATT100921	37 Marsred (Orange)	ATT100137	ATT100837	ATT100937
22 Brilliant Red	ATT100122	ATT100822	ATT100922	40 Kiwi	ATT100140	ATT100840	ATT100840
23 Matt Red	ATT100123	ATT100823	ATT100923	41 Flamingo-Glitter	ATT100141	ATT100841	ATT100941
24 Tutti Frutti Orange	ATT100124	ATT100824	ATT100924	42 Matt Anthracite	ATT100142	ATT100842	ATT100942
27 Graphite-metallic	ATT100127	ATT100827	ATT100927	46 Mamba Green	ATT100146	ATT100846	ATT100946

Sticker colours

€

ATT100550	Black	–
ATT100551	Silver	–

Accessories

A1T060051



// Brake, backpedal actuation

A1T090095



// Mudguard, bolt on, black

A1T090041



// Multitool

A1T090067



// Rear view mirror

A1T090084



// Additional weight, removable, 5 kg, to increase traction

A1T090085



// Bell

A1T090081



// Bottle holder

A1T090024



// High pressure pump (0 - 11 bar)

A1T090082



// Low rider carrier for saddle bags

A1T090081



// Connect Display
 connects with your smartphone
 (download of neodrive app necessary)
 Records trip data such as average speed, maximum speed, trip distance, and total distance.

Comments

Please find some examples for already existing B4me options below:

FMG-MJ02-3187-02 - Tetra operation (ATTITUDE Hybrid)

FMG-MA03-3201-01 Adaptation link -2 to + 6 cm (ATTITUDE Manual, Hybrid & Power)

**ATTITUDE Manual**

€

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date**Signature**

Day Month Year

**Sunrise Medical GmbH**

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: export@sunrisemedical.de

www.SunriseMedical.eu