






After undocking only two small pins (170 g) remain on the wheelchair.

# READY TO GO IN SECONDS

The Attitude Junior's clever docking mechanism attaches to your wheelchair in a matter of seconds...and it's just as quick to remove. Ideal for explorers who don't have time to lose when they want to change between cycling or wheelchair use.

The docking mechanism on the Attitude Junior bike is so simple, it just locks securely to the front of the wheelchair – without any need for parts or cluttered wiring around the knees. So there's nothing coming between them and the comfort that they have come to expect from their wheelchair.



- 1 Lock in the docking pins.
- 2 Swing the upper docking connector downwards.
- 3 Lock the connection by pulling the red lever.





## WHERE ADVENTURES BEGIN







#### // LET IT ROLL

The 8 speed hub-gearing is perfect for all terrains. The powerfull backpedal brake allows for smoother speed control.



## // HORIZONTAL HANDLES

The horizontal push handle position is ideal for children and supports the paediatric arm ergonomy.



### // ERGONIMIC HANDLE

The ergonomic handle helps to manouvre the bike easily with just one hand while undocking or storing it away.



#### // PARKING BRAKE

Easy operating additional parking brake with lock-in function as standard on all Attitudes.















Kahlbachring 2-4 69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0 Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de

www.SunriseMedical.eu

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at

www.SunriseMedical.eu.



