

ATTITUDE Junior



Rev:11/2024
Valid:11/2025

ORDER FORM

Errors and imissions excepted.

Billing Address Or Customer-ID			Shipping Address	Identical To Billing Address
C.-ID:			C.-ID:	
Name:			Name:	
Street:			Street:	
Town/PC:			Town/PC:	
E-Mail:			E-Mail:	
Phone:			Phone:	
Order Information			User Weight	
Order	Calculation	Built-4-Me®	Max. user weight: 100 kg	
Purchase Order			User Weight: kg	
Art. No.: 47500001		Reimbursement Code: xxx		
Order Date:	Day	Month	Year	

Base Price € X.XXX,-

Attitude - Model Selection

€

ATT010804	ATTITUDE Junior 8 speed hub gearing, back-pedal brake and parking brake, grip-shift, rear and front light, angle- and depth-adjustability, 16" front wheel, docking system for fixed-front frames, length adjustable crank (2 positions)	
-----------	--	--

Technical Data Overview

Max. user weight	100 kg	Shifter	Twist grip shifter
Front frame width	260 – 460 mm	Gears	8
Weight of the attachment adapters on the wheelchair	85 g each	Brake	Coaster brake
Total weight	from 21,5 kg	Parking brake	Rim brake
Wheel size	16"	Frame colors	24

Wheelchair Model - which will be used in combination with the ATTITUDE Junior

€

Sunrise Medical Pediatric Wheelchairs		Sunrise Medical Wheelchairs	
Easy Life T fixed front frame only	—	Easy Life fixed front frame only	—
Easy Life RT fixed front frame only	—	Easy Life R fixed front frame only	—
Youngster 3 please see next page for NoGos in combination with the Youngster	—	Xenon/Xenon ² fixed front frame only	—
		Argon ²	—
		Helium	—
		RGK Tiga please note: V-Frames are only suitable up to 5°, Please check with your authorized dealer/ manufacturer if there are any further NoGos	XXX,-

Other Wheelchairs

Sorg Vector fixed front frame only, available from June 2017, please note: V-Frames are only suitable up to 5°	XXX,-
Sorg Jump fixed front frame only, available from June 2017, please note: V-Frames are only suitable up to 5°	XXX,-
PRO ACTIV Speedy / Buddy please note: V-Frames are only suitable up to 5°, Please check with your authorized dealer/ manufacturer if there are any further NoGos	XXX,-

Wheelchair front frame width - measured from outside front frame tube to outside front frame tube

260 mm	380 mm
280 mm	400 mm
300 mm	420 mm
320 mm	440 mm
340 mm	460mm
360 mm	

Fitting-kits for Youngster 3 Necessary to convert wheel lock attachment on the Youngster- see pictures on the next page

€

ATT060150	Retro-fit kit wheel-lock attachment for Youngster (Knee-lever-to Standardbrake),incl. Standardbrake left & right	XXX,-
ATT060151	Retro-fit kit wheel-lock mounting brackets for Youngster with Standardbrake	XXX,-
ATT060152	Retro-fit kit wheel-lock mounting brackets for Youngster with Safari brake	XXX,-


Order process information

€

ATT090213	ATTITUDE only - without a related wheelchair order	—
ATT090214	ATTITUDE with wheelchair order - on a separate wheelchair order form	—

Youngster 3: Wheel-lock conversion

To dock the Attitude junior to a Youngster3, the wheel lock attachment needs to be converted- the upper clamp of the Attitude collides with the standard wheel-lock mounting bracket of the Youngster3. If the Youngster3 is equipped with knee-lever brakes, those need to be swapped for Standard-brakes.



ATT060150 **Retro-fit kit wheel-lock attachment for Youngster (Knee-lever-to Standardbrake),incl. Standardbrake left & right**
(Necessary if your Youngster3 is equipped with knee-lever brakes. This kit includes two Standardbrakes (left & right) and mounting brackets)

ATT060151 **Retro-fit kit wheel-lock mounting brackets for Youngster with Standardbrake**
(Necessary if your Youngster3 is already equipped with Standardbrakes. This kit includes mounting brackets)

ATT060152 **Retro-fit kit wheel-lock mounting brackets for Youngster with Safaribrake**
(Necessary if your Youngster3 is already equipped with Safaribrakes. This kit includes mounting brackets)

Youngster 3: NoGos

Short frame (seat-depth ≤ 300 mm)

wheel size	COG Position			Wheel lock		
	1 (passiv)		2 (active)	knee-lever	Standard	Safari
24"	X		X	X	X	X
22"	✓		X	✓ (1)	✓ (2)	✓ (3)
20"	✓		✓	✓ (1)	✓ (2)	✓ (3)

Long frame (seat-depth > 300 mm)

wheel size	COG Position			Wheel lock		
	1 (passiv)	2 (middle)	3 (active)	knee-lever	Standard	Safari
24"	✓	✓	X	✓ (1)	✓ (2)	✓ (3)
22"	✓	✓	✓	✓ (1)	✓ (2)	✓ (3)
20"	✓	✓	✓	✓ (1)	✓ (2)	✓ (3)

(1) knee-lever brake needs to be converted to Standard-brake with retro-fit kit ATT060150

(2) wheel-lock mounting brackets need to be dismantled and swapped for brackets included in retro-fit kit ATT060151

(3) wheel-lock mounting brackets need to be dismantled and swapped for brackets included in retro-fit kit ATT060152


ATTITUDE JUNIOR fork: Angle and depth-adjustability





Angle adjustment:

from 0° (parallel to the ground) to 28° continuously adjustable

Depth adjustability:

100mm (+/- 50mm compared to standard fork) continuously adjustable



Fork 			€
ATT080036	Fork for Attitude Junior, angle and depth-adjustability of crank-position		–
Tyres 			€
ATT070110	Tyre, pneumatic, tread, Mountainbike		–
ATT070111	Tyre, pneumatic, tread, Big Apple		–
Gear and shifters			€
ATT070043	Gear hub, 8 speed, coaster brake with grip-shift		–
Crank and handle - special length, widths and shapes via B4Me 			€
ATT080037	Junior crank: length adjustable (2 positions), short: length 110/130mm		–
ATT080038	Junior crank: length adjustable (2 positions), long: length 130/150mm		–
ATT080039	Junior Handle, horizontal, round, diameter 22mm Grip-shift mounted to the frame		–
ATT080023	Handle, vertical, 10° angle, round, diameter 30 mm Grip-shift mounted to the right handle		–
	ATT080028 Quad, hand support		XXX,-
	ATT080027 Quad, arm support		XXX,-
ATT080025	Handle, vertical, 25° angle, oval, diameter 37 x 21 mm Grip-shift mounted to the right handle		XXX,-
Brake			€
	Brake, backpedal actuation		–
	Additional parking brake with lever mounted to the frame		–
Accessories 			€
ATT090074	Front- & rear lamp		–
ATT090080	Speedometer		–
ATT090095	Mudguard		XXX,-
ATT090067	Rear view mirror		XXX,-
ATT090081	Bottle holder		XXX,-
ATT090024	High pressure pump 0 - 11 bar		XXX,-
ATT090085	Bell		XXX,-
ATT090041	Multitool		XXX,-
ATT090096	Protection foil		XXX,-

Fork

ATT080036



// Fork for Attitude Junior, angle and depth-adjustability of crank-position

Tyres

ATT070110



// Tyre, pneumatic, tread, Mountain-bike

ATT070111



// Tyre, pneumatic, tread, Big Apple

Crank and handle

ATT070043



// Gear hub, 8 speed, coaster brake with grip-shift

ATT080036



// Junior crank short: length 110/130mm (2 positions)
// Junior Handle, horizontal, round, diameter 22mm

ATT070110



// Handle, vertical, 10° angle, round, diameter 30 mm

ATT070111



// Quad, hand support

ATT070110



// Quad, arm support

ATT070111



// Handle, vertical, 25° angle, oval, diameter 37 x 21 mm

Accessories

ATT090095



// Mudguard

ATT090067



// Rear view mirror

ATT090081



// Bottle holder

ATT090024



// High pressure pump (0 - 11 bar)

ATT090085



// Bell

ATT090041



// Multitool

Chain guard Colours You can choose separate colours for the frame, the cranks and the chain guard.

	Frame	Chain guard	Crank		Frame	Chain guard	Crank
02 Brilliant Silver	ATT100102	ATT100802	ATT100902	29 Black	ATT100129	ATT100829	ATT100829
03 Silver Matt	ATT100103	ATT100803	ATT100903	31 Matt Black	ATT100131	ATT100831	ATT100931
06 Sun Yellow	ATT100106	ATT100806	ATT100906	32 Signal White	ATT100132	ATT100832	ATT100932
14 Blue	ATT100114	ATT100814	ATT100914	33 Pink	ATT100133	ATT100833	ATT100933
15 Brilliant Marine	ATT100115	ATT100815	ATT100915	34 Dark Chocolate Metallic	ATT100134	ATT100834	ATT100934
16 Matt Marine	ATT100116	ATT100816	ATT100916	35 Skyblue	ATT100135	ATT100835	ATT100935
19 Brilliant Purple	ATT100119	ATT100819	ATT100919	36 Cyan-Metallic	ATT100136	ATT100836	ATT100936
21 Red Ruby	ATT100121	ATT100821	ATT100921	37 Marsred (Orange)	ATT100137	ATT100837	ATT100937
22 Brilliant Red	ATT100122	ATT100822	ATT100922	40 Kiwi	ATT100140	ATT100840	ATT100840
23 Matt Red	ATT100123	ATT100823	ATT100923	41 Flamingo-Glitter	ATT100141	ATT100841	ATT100941
24 Tutti Frutti Orange	ATT100124	ATT100824	ATT100924	42 Matt Anthracite	ATT100142	ATT100842	ATT100942
27 Graphite-metallic	ATT100127	ATT100827	ATT100927	46 Mamba Green	ATT100146	ATT100846	ATT100946

Sticker colours

€

ATT100550	Black	—
ATT100551	Silver	—

Comments for special build - please apply for a quote; different lead time and discount structure.

ATTITUDE Manual

€

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date
Signature

Day Month Year


Sunrise Medical GmbH

Kahlbachring 2-4
D-69254 Malsch/HD
Phone: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: export@sunrisemedical.de
www.SunriseMedical.eu