Self propelled Easy to operate for the independent user.

Safe and simple transfer The flat seat and the brake position facilitate lateral transfer.

6

0


Etac Clean 24"

This variant of Etac Clean is suitable for those who can operate the chair independently and want to maximise their mobility.

Easy to manoeuvre, for the independent user

The Etac Clean 24-inch wheels are positioned to make it easy to manoeuvre the chair without compromising user safety or comfort. It is easy to transfer to and from the chair thanks to the large, unrestricted seat in front of the wheels. The brake handles are also positioned below the seat so that they do not get in the way during transfer, while remaining easy to reach.

Etac Clean 24" has the same basic benefits as the other members of the Clean family – but allows the user to operate the chair independently.





Etac Clean 24" 55 cm (21¾) with 24" rear wheels, white 55 cm (21¾) with 24" rear wheels, lagoon green

Width 69 cm

Width between arm supports 43.5 cm

Seat height 55 cm

Seat width

48 cm Seat depth

44 cm

Weight 20 kg

Maximum user weight 130 kg

Material

80229276

80229216

Rear wheel adapter: aluminum Rear Wheels: polyamide, solid tires Chassis: Powder-painted steel tubing Seat/back support: Polypropylene Arm supports: Polypropylene Foot support: Polypropylene Castor brackets: Polyamide Castors, lockable, ball-bearing: Polyamide/TPE





 $\begin{array}{l} \textbf{Maintenance} \\ \textbf{Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70 % disinfectant solution. Can be decontaminated at max 85°C / 185°F is so required. \end{array}$

Accessories





Soft back support 80209225 Adjustable with velcro straps. Washing at max 60°C (140°F).

Side Support 80209509 Adjustable in width. Material: Polyester fabric, Velcro, polyethylene. Washing at max. 60°C (140°F).



Crossbar 80209246 Attached to the arm supports. Material: Aluminium, cellular polyethylene, polypropylene.



Safety strap 80209034 Attachments on the rear side of the back support. Material: Polyester.



Positioning belt 80209440 Adjustable in length. Attachment with velcro straps. Material: Polyester.



Chest harness strap 80209442 Adjustable in length. Attach to back support. Always use together with positioning belt.





Comfort Kit

Material: Polyester.

80209237 Includes both Comfort Seat Soft and Comfort Cover. Cleaning max. 60° C (140 °F) or with alcohol wipe.

Comfort Seat Soft

Opening width 15 cm (6"): 80209301 Opening width 18 cm (7"): 80209227

A soft contoured seat in two versions, small opening: W 15 cm (6") or narrow opening: W 18 cm (7"). Material: polyurethane. Cleaning max. 85° C (185°F) or with alcohol wipe.

















Comfort Cover

80209228 A soft back support with integrated arm support pads. Material: polyethylene. Cleaning max. 60° C (140 °F) or with alcohol wipe.

Comfort Arm Pads 80209226

Soft and warm, leaving the handles free to ensure a firm grip when standing up. Material: polyethylene. Cleaning max. 60° C (140 °F) or with alcohol wipe.

Cushion, protect

80209245 Air-inflated ROHO cushion with excellent pressure distribution properties. Adjustable air pressure. Recess width: 21 cm (8¼"). Material: Neoprene. Cleaning max: 60°C (140°F).

Soft pad

Grey: 80209260 Heat-reflecting with a non-slip surface. Recess width: 17 cm (6¾"). Material: polyethylene. Cleaning max max. 60° C (140 °F) or with alcohol wipe.

Soft back pad

Grey: 80209261 Heat-reflecting with a non-slip surface. Material: polyethylene. Cleaning max. 60°C (140°F) or with alcohol wipe.

Seat cover

80209247 When a complete seat is required. Material: Polyurethane. Cleaning max. 85°C (185°F) or with alcohol wipe.

Soap dish (does not fit Clean 24") Grev: 80209266

Attaches to the side of the seat. Built-in shower handle holder. Material: Polypropylene. (Does not fit Clean 24")

Pan Grey: 80209254 Material: polypropylene.

Pan incl. lid with handle Grey: 80209255 Grey lid only, with handle: 80209256 Material: Polypropylene.

Accessories



Pan holder

Grey: 80209257 Fits both pan, bedpan and regular hospital pans with a thin rim.



Bucket with handle and lid 80302026 The bucket holder is not included. Material: Polypropylene.





Attachment for low mounting of pan/bucket Grey: 80209259 Fitted to the pan/bucket holder. Lowers the pan/bucket holder with

6 cm (2¹/₂"). Material: Polypropene.

Disposable pans, 50 pc. 80209280 Strong and firm. 50 pcs.



Bucket holder Grey: 80209258 Fits Bucket with handle and lid and other buckets with a thick rim. Material: Polypropene.



Disposable pan holder kit Grey: 80209229 Fits Disposable pan.