



JAY®

Balance

Rev 2021

Invoice to:

Name:

Street:

Town:

Postal code:

E-Mail:

Tel:

Fax:

Deliver to:

Name:

Street:

Town:

Postal code:

E-Mail:

Tel:

Fax:

Order date:

Purchase order:

Number required:



Marked for:

Sunrise order no:

- ☐ Order
☐ Calculation

Jay® Balance Standard Features

Fluid pad or Air (single or dual valve) insert to optimize immersion and accommodate positioning needs. Positioning accessories available.

Innovative dual-cover system: water-resistant inner cover, easy to clean, combined with 3 types of machine-washable outer covers: stretch, microclimatic and incontinent.

2 Year guarantee; 6 months on cover.

Max. user weight:

150 kg (cushion width 34 to 50 cm)

225 kg (cushion width 56 to 60 cm)

USER WEIGHT kg

Mandatory information to ensure correct configuration

Reimbursement code:Application:

Client at high risk of skin breakdown with or without high postural needs.

Symmetrical to significant asymmetrical posture requiring laterally, forward/rearward stability and/or positioning.

May be unable to reposition or perform an independent pressure relief.

To order a Jay Balance cushion, follow the steps below to build your part number and then choose your accessories

fixed	Width	Depth	Insert	Outer Cover	Pos. Kit	fixed
JBALC	40	42	FL	M	P	EU
Jay Balance Cushion	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	


Example: JBALC4042FLMPEU

CUSHION WIDTH AND DEPTH - STEP 1 + STEP 2

Standard Cushion (width: 34 - 50 cm; max user weight 150 kg)

Bariatric Cushion (width: 56 - 60 cm; max user weight 225 kg)

		Width										
Depth	cm	34	36	38	40	42	44	46	48	50	56	60
	34	○										
	36	○	○	○								
	38	○	○	○	○	○	○					
	40	○	○	○	○	○	○					
	42	○	○	○	○	○	○	○	○	○		
	44		○	○	○	○	○	○	○	○		
	46		○	○	○	○	○	○	○	○	○	○
	48					○	○	○	○	○	○	○
	50						○	○	○	○	○	○
56										○	○	
60										○	○	
PLA Size		A			B					C		

Jay Balance Measuring Guide

Width X cm - Measurement made from outside edge to outside edge of cushion (including cover)

Depth Y cm - Measurement made from front of seat notch to front topside (or bottom of contracture cut) of cushion (including cover)

WIDTH x DEPTH

STEP 1 STEP 2



For more information
visit our webpage

www.sunrisemedical.com

© Sunrise Medical 2021

JAY BALANCE CUSHION

JAY Balance Lead Time - for the below listet cushions the lead time is 5 days; for all others the lead time is 15 days

Seat Width	Seat Depth	Balance Part No.
34	34	JBALC3434FLMPEU
	36	JBALC3640FLMEU
	40	JBALC3440FLMEU
36	36	JBALC3636FLMEU
	40	JBALC3640ASMEU JBALC3640FLMEU JBALC3640FLMPEU
	42	JBALC3642FLMEU JBALC3642FLMPEU
	44	JBALC3644FLMPEU
	46	JBALC3646FLMEU JBALC3646FLMPEU
	46	JBALC3646FLMPEU
40	40	JBALC4040ADMEU JBALC4040ASMEU JBALC4040FLIEU
	42	JBALC4042ADMEU JBALC4042ASMEU JBALC4042FLMPEU
	44	JBALC4044ASMEU JBALC4044FLMEU JBALC4044FLMPEU
	46	JBALC4046ADMEU JBALC4046FLIEU JBALC4046FLMPEU
	40	JBALC4240FLIEU JBALC4240FLMEU JBALC4240FLMPEU
	42	JBALC4242ADMEU JBALC4242ADMPEU JBALC4242ASMPPEU JBALC4242FLIEU JBALC4242FLIPEU JBALC4242FLMEU JBALC4242FLMPEU JBALC4242FLSPEU
	44	JBALC4244ADMEU JBALC4244FLIEU JBALC4244FLIPEU JBALC4244FLMEU JBALC4244FLMPEU
	46	JBALC4246ADMEU JBALC4246FLIEU JBALC4246FLIPEU JBALC4246FLMEU JBALC4246FLMPEU JBALC4246FLSPEU
42	48	JBALC4248FLIEU JBALC4248FLIPEU JBALC4248FLMEU JBALC4248FLMPEU JBALC4248FLSEU
	40	JBALC4440FLMEU JBALC4440FLIEU JBALC4440FLMPEU
	42	JBALC4442FLIPEU JBALC4442FLIEU JBALC4442FLMEU JBALC4442FLMPEU
	44	JBALC4444ADMEU JBALC4444ADMPEU JBALC4444ASMPPEU JBALC4444ASMEU JBALC4444FLIPEU JBALC4444FLIEU JBALC4444FLMEU JBALC4444FLMPEU JBALC4444FLSPEU
	46	JBALC4446ASSEU JBALC4446ADMEU

Seat Width	Seat Depth	Balance Part No.
44	46	JBALC4446ADMPEU JBALC4446ASMPPEU JBALC4446FLSPEU JBALC4446FLIPEU JBALC4446FLIEU JBALC4446FLMEU JBALC4446FLMPEU JBALC4446FLSEU
	48	JBALC4448FLMEU JBALC4448FLMPEU
	50	JBALC4450FLMEU JBALC4450FLIPEU JBALC4450ASMPPEU JBALC4450FLMPEU
	42	JBALC4642FLIEU JBALC4642FLMEU JBALC4642FLMPEU
	44	JBALC4644FLMEU JBALC4644FLMPEU
	46	JBALC4646ADMEU JBALC4646ASIEU JBALC4646FLSPEU JBALC4646ASMPPEU JBALC4646ADMPEU JBALC4646ASMEU JBALC4646FLIEU JBALC4646FLIPEU JBALC4646FLMEU JBALC4646FLMPEU
46	48	JBALC4648FLIEU JBALC4648FLMEU JBALC4648FLMPEU JBALC4648FLSEU
	50	JBALC4650FLIEU JBALC4650FLMEU JBALC4650FLIPEU JBALC4650FLMPEU JBALC4650FLSEU
	42	JBALC4842FLMPEU
	44	JBALC4844FLMPEU JBALC4844FLMEU
	46	JBALC4846FLIEU JBALC4846FLMEU JBALC4846FLMPEU
	48	JBALC4848FLMEU JBALC4848FLMPEU
48	50	JBALC4850FLIEU JBALC4850FLMEU JBALC4850FLMPEU
	42	JBALC5042FLMEU JBALC5042FLMPEU
	44	JBALC5044FLMEU JBALC5044FLMPEU
	46	JBALC5046FLIEU JBALC5046FLMEU JBALC5046FLMPEU
	48	JBALC5048FLMEU JBALC5048FLMPEU
	50	JBALC5050ASMPPEU JBALC5050FLIEU JBALC5050FLMEU JBALC5050FLMPEU
50	46	JBALC5646FLMEU
	50	JBALC5650FLMEU JBALC5650ASMEU JBALC5650FLMPEU
	56	JBALC5656FLMPEU JBALC5656FLMPEU
	56	JBALC5656FLMPEU




For more information
visit our webpage
www.sunrisemedical.com
© Sunrise Medical 2021

JAY BALANCE CUSHION

PELVIC LOADING AREA (PLA) INSERT* - STEP 3

	Select	Part #	Charge
FLUID - Fluid Pad Option	☐	FL	-
AIR Insert Option - Single Valve**	☐	AS	add XX €
AIR Insert Pad Option - Dual Valve**	☐	AD	add XX €



INSERT

STEP 3
(FL, AS or AD)

* The insert pad sizes supplied are based on PLA size and automatically determined when you select the cushion (see matrix)

**Air Bladder Options shipped with instructions, air pump, patch kit and nozzle fittings

OUTER COVER OPTION - STEP 4

	Select	Part #	Charge
Microclimatic	☐	M	-
Stretch	☐	S	-
Incontinent	☐	I	-

OUTER COVER

STEP 4
(M, S or I)

All cushions are delivered with an water resistant inner cover.

POSITIONING KIT *) - STEP 5

If you wish your JAY Balance to have positioning components, please write "P" in the Step 5 box. Otherwise, leave the box empty to get a standard cushion (without positioning components). When ordering a positioning version, you will get the following positioning kit:

	Part #	Charge
Positioning Kit (containing each 1 of the components mentioned below)	P	€
Lateral Thigh 2.5 cm Left		
Lateral Thigh 2.5 cm Right		
Medial Thigh 2.5 cm (1") (each)		
Lateral Pelvic 2.5 cm (2") Left		
Lateral Pelvic 2.5 cm (2") Right		
Lateral Pelvis LH		
Lateral Pelvis RH		
Lateral Thigh LH - 2,5 cm		
Medial Thigh - 2,5 cm		
Lateral Thigh RH - 2,5 cm		

POS. KIT

STEP 5
(P or blank)

*) all positioning components are delivered with a stretchable water resistant inner cover to accommodate these components onto the Balance cushion.

JAY Balance Smart Part Number

	Width	Depth	Insert	Outer Cover	Pos. Kit	
JBALC	40	42	FL	M	P	EU
fixed	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	fixed

Notes

The part # above (JBALC4042FLMPEU) describes a Jay Balance cushion with 40 width and 42 depth, fluid insert, microclimatic outer cover and positioning components.

Additional explanations:

When selecting the PLA insert and the positioning option for your cushion, the water resistant inner cover will be automatically selected (here in this example it will be an inner cover which accommodates fluid insert and positioning components).

If any modifications of your cushion is needed please contact your Sunrise Medical authorized Customer Service.

JAY BALANCE COSTS

	Unit price	Total price
Cushion Price		
Option Price		
Total Cost		

Sunrise Medical Customer service
Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg
Tel. +497253 980-400
Fax. +497253 980-222

For more information
visit our webpage
www.sunrisemedical.com
© Sunrise Medical 2021

Spare Parts			Charge
Cushion base (ww = width; dd = depth)			
☐	BJBALCwwddEU	Base assembly	€
Spare Outer Cover (ww = width; dd = depth)			
☐	CJBALCwwddSEU	Stretch outer cover	€
☐	CJBALCwwddMEU	Microclimatic outer cover	€
☐	CJBALCwwddIEU	Incontinent outer cover	€
Spare Inner Cover (ww = width; dd = depth)			
☐	CIFLJBALCwwddEU	Inner cover for fluid cushion	€
☐	CIFLPJBALCwwddEU	Inner cover for fluid and positioning cushion	€
☐	CIAJBALCwwddEU	Inner cover for air cushion	€
☐	CIAPJBALCwwddEU	Inner cover for air and positioning cushion	€
Solid Seat (ww = width; dd = depth)			
☐	JSIwwdd	Solid Seat insert (wood plate)	€
Fluid & Air inserts *			Charge
PLA A			
☐	FLJBALCAEU	Fluid insert for PLA size A	€
☐	J3CSAAEU	Air insert - single valve - for PLA size A	€
☐	J3CSAADEU	Air insert - dual valve - for PLA size A	€
PLA B			
☐	FLJBALCBEU	Fluid insert for PLA size B	€
☐	J3CSBASEU	Air insert - single valve - for PLA size B	€
☐	J3CSBADEU	Air insert - dual valve - for PLA size B	€
PLA C			
☐	FLJBALCCEU	Fluid insert for PLA size C	€
☐	J3CSCASEU	Air insert - single valve - for PLA size C	€
☐	J3CSCADEU	Air insert - dual valve - for PLA size C	€
* Note: A Balance cushion with Air insert has a different inner cover as with fluid insert. In case you want to change the insert for an existing Balance cushion, you need to order a suitable inner cover.			
JAY Balance Positioning kit *			Charge
	for depth 34-40 cm	for depth 42-60 cm	
Lateral Thigh Left 2.5cm Left	☐ JBALCLTLH-S	☐ JBALCLTLH	€
Lateral Thigh Left 2.5cm Right	☐ JBALCLTRH-S	☐ JBALCLTRH	€
Medial Thigh 2.5cm (each)	☐ JBALCMT-S	☐ JBALCMT	€
Lateral Pelvic 2.5 cm Left	☐ JBALCLPLH-S	☐ JBALCLPLH	€
Lateral Pelvic 2.5 cm Right	☐ JBALCLPRH-S	☐ JBALCLPRH	€
Positioning kit	☐ JBALCPKIT-S	☐ JBALCPKIT	€
(Positioning kit is containing 1 of each of the components above)			
* Note: A Balance Standard cushion has a different inner cover as the Positioning version. In case you want to change from Standard to Positioning (and versa vice) you need to order a suitable inner cover			
JAY Balance Obliquity Pads			Charge
Order 1 or 2 obliquity pads per side, depending on the level of obliquity that needs to be corrected. The pad has hook/loop on it and can be stacked (maximum 2).			
Obliquity pads can be stacked the following way: for PLA size A: order the obliquity pad size A - one side (RH or LH) - and additionally order the size B (RH or LH). For PLA size B, order the obliquity pad size B - one side (RH or LH) - and additionally order the size C (RH or LH). For PLA size C, order the obliquity pad size C - one side (RH or LH) - and additionally order the size D (RH or LH). Please notice that size D can only be ordered in combination with size C otherwise it will not fit in the PLA.			
☐	JBALCOP-ALH	Obliquity pad (left) for PLA A	€
☐	JBALCOP-BLH	Obliquity pad (left) for PLA B	€
☐	JBALCOP_CLH	Obliquity pad (left) for PLA C	€
☐	JBALCOP-DLH	Additional obliquity pad (left) that can be stacked over obliquity pad for PLA C	€
☐	JBALCOP-ARH	Obliquity pad (right) for PLA A	€
☐	JBALCOP-BRH	Obliquity pad (right) for PLA B	€
☐	JBALCOP-CRH	Obliquity pad (right) for PLA C	€
☐	JBALCOP-DRH	Additional obliquity pad (right) that can be stacked over obliquity pad for PLA C	€
Note: Cushions or inner covers purchased before February 2017 do not contain appropriate velcro attachment for the obliquity pads.			
			For more information visit our webpage www.sunrisemedical.com © Sunrise Medical 2021