

JAY Zip Back

Clinical seating, designed for zippie kids!



The JAY Zip backrest was designed through the use of anthropometric data so they are specially sized to meet a pediatric client's unique needs. It is an ultra lightweight, antimicrobial backrest that features X-static™ and Aquaguard™ cover technologies as well as adjustable, one step release hardware to accommodate a wide range of mounting locations and user presentations.

Versatile One Step Release Hardware

- Easy one step release
- Fits on wide range of tubing: 1,9 cm with inserts; 2,2 - 2,3 cm with inserts or 2,5 cm without insert
- One tool adjustments



Convenient & Antimicrobial

Outer Cover

- Machine washable and quick to dry
- Silver impregnated X-static® stretch fabric naturally inhibits bacterial growth
- Extra layer of microclimatic spacer fabric for comfort
- Safe non-toxic material

Inner Cover

- Easy to wipe off
- Aquaguard™ zipper resists moisture
- Anti-wicking seam thread protects foam base



One step release

Backrest width:	20 / 25 / 30 / 35 cm	Back angle:	+20° to -20°
Backrest height:	15 / 20 / 25 / 30 cm	Width adjustability:	+2,5 cm
Transport weight:	1,1 kg	Depth adjustability:	7,5 cm
Max. user weight:	75 kg	Accessories: Headrest, Lateral support, spine align positioning components, chest strap, vanity flap and spare cover	