

**Life**Edition: 11.24
Valid until: 11.25

Invoice to:
Name:
Street:
Town: Postal code:
E-Mail:
Tel: Fax:

Deliver to:
Name:
Street:
Town: Postal code:
E-Mail:
Tel: Fax:

Order date:
Purchase order:
Number required:



Marked for:
Sunrise order no:

☐ Order ☐ Calculation ☐ Built 4 me

• = Standard
○ = No cost option

☐ = Cost option

Art. no.: 07200101

Max. user weight:

140 kg

USER WEIGHT

kg

Mandatory information to ensure correct configuration

Base Price

2.599 €

CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH SELECTION - measured from outside frame tube to outside frame tube (mm)

360	380	400	420	440	460	480	500	520	540
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

SEAT DEPTH SELECTION - measured in mm

360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Short Frame			Medium Frame			Long Frame	

FRONT SEAT HEIGHT SELECTION: in dependence of front castors and frame versions

		Front seat height in mm													
		Fork Type (mm)													
Front castor	4"	Hemi Version	<input type="checkbox"/>	<input type="checkbox"/>											
		Aluminium (143 mm)			<input type="radio"/>	<input type="radio"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
		Aluminium (158 mm)						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		
	5"	Hemi Version	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
		Aluminium (143 mm)			<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>				
		Aluminium (158 mm)					<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		
	6"	Hemi Version			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>								
		Aluminium (143 mm)						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
		Aluminium (158 mm)							<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	7"	Hemi version only with SA frame									<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

REAR SEAT HEIGHT SELECTION: in dependence of rear wheel size

		Rear seat height in mm										
		Rear wheels	400	410	420	430	440	450	460	470	480	490
	22"		<input type="radio"/>	<input type="radio"/>		<input type="radio"/>	<input type="radio"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
	24"				<input type="radio"/>	<input type="radio"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
	25"					<input type="radio"/>	<input type="radio"/>	<input type="radio"/>		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

For more information
www.SunriseMedical.eu
© Sunrise Medical 2025

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair
(approx. 10 - 20 mm altogether)
- b) On the wheelchair: From outside frame tube to outside frame tube.

Seat depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between
front edge of the seat upholstery and back of knee.
(approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

- a) For the user: When sitting with 90° kneeangle:
From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground.
(min. 30 mm longer than user measurement)
Remember to include the cushion height.

Rear seat height:

- b) On the wheelchair: From back of seat tube to the ground.
Note: Max. seat height difference rear to front is 130 mm

Back height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.
Note: The user's ability to sit upright has a major influence on the
back height and it can therefore differ considerably

Back angle:

Important: The back height depends on the degree of disability.
0° is the angle 90° to the ground. Negative angles are < 90°
and positive angles are > 90°.

Quickie Life

● = Standard ○ = No cost option □ = Cost option

€

SEAT UPHOLSTERY AND CUSHIONS

Seat Slings

LFE020014	● Seat sling, nylon, 1 utility bag	Std
LFE020000	○ Without seat sling	-

Seat cushion, cover with zipper, black

LFE020101	□ Seat cushion, foam (medium-soft), cover with zipper, black, 30 mm	116 €
LFE020102	□ Seat cushion, foam (medium-soft), cover with zipper, black, 50 mm	116 €

Jay Cushion - for complete Jay product range see Jay order forms

JAY020002	□ JAY Basic	163 €
JAY020003	□ JAY Soft Combi P	238 €
JAY020001	□ JAY Easy Visco C	349 €
JAY020007	□ JAY Easy Fluid	602 €
JAY020004	□ JAY LITE	551 €

Seat accessories

LFE090018	□ Lap belt with buckle, metal	87 €
-----------	-------------------------------	------

BACKREST

Backrest height range	250-300	300-350	350-400	400-450	450 - 500
Back height set at (mm)	○ 250	○ 300	○ 350	○ 400	○ 450
	○ 275	○ 325	○ 375	○ 425	○ 475
	○ 300	○ 350	○ 400	○ 450	○ 500

LFE030010	● Backrest, fixed (90°)	Std
LFE030101	□ Backrest, half-folding (to the back) (from back height 450 mm-475 mm)	340 €
LFE030031	□ Backrest, fold-down (to the front) (4) (from back height 300 mm-475 mm)	362 €
LFE030012	□ Backrest, angle-adjustable (not with BH 500 mm)	286 €
	○ 78° ○ 82° ○ 86°	-
	○ 90° ○ 94° ○ 98° ○ 102°	-

Push handles

LFE030201	○ Push handles long	-
LFE030203	□ Push handles, fold-down (1)(5) (only from back height 300 mm; from 350 mm with angle adj back, not with JAY Shape)	195 €
LFE030204	□ Push handles, height-adjustable (1) (2) (only from back height 300 mm)	310 €
LFE030205	□ Push handles, height-adjustable for Jay backs (1), (2) (only from back height 300 mm)	310 €
LFE030200	○ Without push handles (1), (2), (5)	-
FMG3217	FMG Pushhandles long (+30mm) for attendee (not with height adj pushhandles)	B4me

Backrest upholstery (black)

LFE030316	○ EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black (EXO Evo version) Colour, binding: ● Black ○ Orange ○ Silver ○ Red ○ Blue	-
LFE030317	□ EXO PRO Backrest sling, padded, adjustable, breathable, spacer fabric, 1 accessory bag, black (EXO Evo version) Colour, binding: ● Black ○ Orange ○ Silver ○ Red ○ Blue	179 €
LFE030000	○ Without backrest upholstery	-

Jay Shape

JAY080001	□ JAY SHAPE Back (SW 360-500 mm, from BH 375 mm, from RSH 440 mm)	301 €
JAY080011	□ JAY SHAPE Lateral supports 10x12 (Pair)	275 €
JAY080010	□ JAY SHAPE Lateral supports 15x15 (Pair)	275 €

Backrest accessories

LFE030402	□ Stabilising bar, automatic folding (1)	172 €
LFE030410	□ Headrest, tall height, depth and angle-adjustable, 230 x 110 mm (1) (2) (3)	412 €
LFE030411	□ Headrest, medium height, depth&angle-adjustable, 200x110 mm (1) (2) (3)	411 €

(1): Not in combination with backrest, half-folding (to the back), R20
 (2): Not in combination backrest backrest, fold-down (to the front), R20
 (3): Only with push handles height adjustable (LFE030204, LFE030205)
 (4): Backheight 350-425 mm with pushhandles fold down; 300-375 mm with push handles long
 (5): Not in combination with R20

For more information
 visit our webpage
www.SunriseMedical.eu
 © Sunrise Medical 2025

SEAT UPHOLSTERY AND CUSHIONS



Seat cushion, cover with zipper, black




Lap belt with buckle, metal



Jay Basic
(picture shown without upholstery)



Jay Soft Combi P
(picture shown without upholstery)



Jay Easy Visco
(picture shown without upholstery)



Jay Easy Fluid
(picture shown without upholstery)



Jay Lite

BACKREST



Angle adjustable back
-12° to 12°



Backrest, Half Folding



Backrest, fold down (to the front)



Push handles long



Fold down push



Height adjustable
push handles



EXO
backrest sling (EXO



EXO PRO
backrest sling (EXO Evo version)



Autofolding
Stabilizer Bar

Quickie Life								
● = Standard ○ = No cost option □ = Cost option €								
JAY BACKS plus ACCESSORIES (max . user weight 136 kg)								
Jay 3 Back - Shallow Contour								
Mid Thoracic:		□ JAY030010	Medium 340mm	□ JAY030033	Tall 380mm	889 €		
Upper Thoracic:		□ JAY030012	Short 420mm	□ JAY030019	Medium 460mm	986 €		
Shoulder Height:		□ JAY030015	Short 530mm	□ JAY030022	Medium 570mm	□ JAY030028	Tall 610mm	1.056 €
Jay 3 Back - Mid Contour								
Mid Thoracic:		□ JAY030011	Medium 340mm	□ JAY030018	Tall 380mm	972 €		
Upper Thoracic:		□ JAY030013	Short 420mm	□ JAY030020	Medium 460mm	□ JAY030026	Tall 500mm	1.085 €
Shoulder Height:		□ JAY030016	Short 530mm	□ JAY030023	Medium 570mm	□ JAY030029	Tall 610mm	1.141 €
JAY J3 Back - Mid Deep Contour (150 mm contour depth)								
Mid Thoracic:		□ JAY030044	Medium - 340 mm	□ JAY030041	Tall 380mm	1.115 €		
Upper Thoracic:		□ JAY030045	Short - 420 mm	□ JAY030042	Medium 460mm	□ JAY030046	Tall - 500 mm	1.115 €
Shoulder Height:		□ JAY030043	Short 530mm	□ JAY030047	Medium - 570 mm	□ JAY030048	Tall - 610 mm	1.141 €
JAY J3 Back - Mid Deep Contour (Extra Deep, 180 mm contour depth)								
Mid Thoracic:		□ JAY030049	Tall - 380 mm	1.115 €				
Upper Thoracic:		□ JAY030050	Medium - 570 mm	1.115 €				
Jay 3 Back - Deep Contour								
Upper Thoracic:		□ JAY030014	Short 420mm	□ JAY030021	Medium 460mm	□ JAY030027	Tall 500mm	1.170 €
Shoulder Height:		□ JAY030017	Short 530mm	□ JAY030024	Medium 570mm	□ JAY030030	Tall 610mm	1.198 €
Jay Easy Back - Posterior Deep Contour								
WHT030415		□ Whitmyer Headrest, COBRA XTRA, flip back	(only with JAY J3 backs)	564 €				
WHT030416		□ Whitmyer Headrest, AXYS, Standard	(only with JAY J3 backs)	338 €				
ARMREST AND SIDEGUARDS								
Armrest								
□ Single post height adjustable armrest		(not in combination with sideguards; only with 0° and 1° camber)						
LFE040050		○ Single-post sideguard, armpad (250 mm), height-adjustable, removable		265 €				
LFE040051		○ Single-post sideguard, armpad (360 mm), height-adjustable, removable		265 €				
LFE040052		□ Single-post sideguard, armpad (330 mm), depth adjustable, tool height-adjustable, removable		265 €				
Desk armrest with long and short pad (1)								
LFE040005		□ Desk sideguard, flip-up, removable		319 €				
LFE040144		□ Armpad for hemiplegic XL (length: 540 mm, width: 135 mm), left		258 €				
LFE040145		□ Armpad for hemiplegic XL (length: 540 mm, width: 135 mm), right		258 €				
LFE040146		□ Armpad for hemiplegic (length: 460 mm, width: 105 mm), left		258 €				
LFE040147		□ Armpad for hemiplegic (length: 460 mm, width: 105 mm), right		258 €				
LFE040001		□ Desk sideguard, armpad (250 mm), fixed, flip-up, removable		319 €				
LFE040144		□ Armpad for hemiplegic XL (length: 540 mm, width: 135 mm), left		258 €				
LFE040145		□ Armpad for hemiplegic XL (length: 540 mm, width: 135 mm), right		258 €				
LFE040146		□ Armpad for hemiplegic (length: 460 mm, width: 105 mm), left		258 €				
LFE040147		□ Armpad for hemiplegic (length: 460 mm, width: 105 mm), right		258 €				
LFE040003		□ Desk sideguard, armpad (310 mm), fixed, flip-up, removable		319 €				
LFE040144		□ Armpad for hemiplegic XL (length: 540 mm, width: 135 mm), left		258 €				
LFE040145		□ Armpad for hemiplegic XL (length: 540 mm, width: 135 mm), right		258 €				
LFE040146		□ Armpad for hemiplegic (length: 460 mm, width: 105 mm), left		258 €				
LFE040147		□ Armpad for hemiplegic (length: 460 mm, width: 105 mm), right		258 €				
LFE040002		□ Desk sideguard, armpad (250 mm), height-adjustable, flip-up, removable		431 €				
LFE040004		□ Desk sideguard, armpad (310 mm), height-adjustable, flip-up, removable		430 €				
LFE040020		□ Release assistance for tetraplegics		-				
Sideguards (not in combination with armrest)								
LFE040107		□ Sideguard, composite, removable, 210 mm high (310 mm long), black		-				
LFE040102		□ Sideguard, aluminium, cold weather protection (in frame colour)	(effective seat width is sw +20 mm)	238 €				
LFE040103		□ Sideguard, aluminium, fender, cold weather protection (in frame colour)	(effective seat width is sw +20 mm)	290 €				
LFE040104		□ Sideguard, carbon fibre, fender, cold weather protection (black)	(effective seat width is sw +20 mm)	301 €				
TRAYS								
Slide on tray								
LFE090011		□ Tray, push-fit, SB 380-480 mm	(SW 380 mm to 480 mm)	306 €				
Swing away tray								
LFE090012		□ Tray, swing-away (to the side), SW 340 - 600 mm, left-mount	(only with desk armrest height adjustable)	389 €				
LFE090013		□ Tray, swing-away (to the side), SW 340 - 600 mm, right-mount	(only with desk armrest height adjustable)	389 €				
FRAME								
LFE010010		○ Frame FF 75 ° (integrated footrest)						
LFE010011		□ Frame FF 85 ° (integrated footrest)						
LFE010016		□ Frame FF 85 ° (integrated footrest), inset (20 mm on each side), active		161 €				
LFE010015		□ Frame SA (removable footrest)						
LFE090102		□ With crash test option		106 €				
(1): Not in combination with Half folding back < BH 400mm. Only in combination with 0° & 1° wheel camber; Fold down backrest only from BH 400 mm								
For more information www.SunriseMedical.eu © Sunrise Medical 2025								

JAY J3 Backs

JAY J3 Back
Shallow Contour



JAY J3 Back
Mid Contour



JAY J3 Back
Mid Deep Contour



JAY J3 Back
Deep Contour



ARMREST AND SIDEGUARDS



Single post height
adjustable armrest



Single-post sideguard,
depth adjustable, tool
height-adjustable



Sideguard composite
detachable



Sideguard carbon with
fender



Sideguard alu with fender



Desk armrest,
height adjustable, short armpad(25cm)



Desk armrest,
fixed height, short armpad(25cm)



Desk armrest
without armpad

Frame



Life SA



Life FF 85°



Life FF 75°



Frame FF 85°
(integrated footrest),
inset (20 mm on each
side), active

Quickie Life			€
● = Standard ○ = No cost option □ = Cost option			
CENTER OF GRAVITY (COG) measured from front of back post; choose only one option			
COG23	○ Centre of gravity 23 mm		-
COG43	● Centre of gravity 43mm		Std
COG63	○ Centre of gravity 63mm		-
COG83	○ Centre of gravity 83mm		-
FRAME ACCESSORIES			
LFE090010	□ Transit wheels	(Not with crutch holder, tip assist and anti tip)	286 €
Tip assist (step tube) & Crutch holder (choose either tip assist or anti-tip tubes, not both)			
LFE090001	□ Tip assist, left		72 €
LFE090002	□ Tip assist, right		72 €
LFE090003	□ Tip assist, pair		143 €
LFE090009	□ Crutch holder, loop	(only in combination with tip assist)	41 €
Anti-tip tubes, swing-away (choose either tip assist or anti-tip tubes, not both)			
LFE090004	□ Anti-tip tubes, swing-away, left		166 €
LFE090005	□ Anti-tip tubes, swing-away, right		166 €
LFE090006	□ Anti-tip tubes, swing-away, pair		333 €
LFE090036	□ Anti-tip tubes, aluminium, lightweight, swing-away (forwards), pair	(not with tip assist, max. user weight: 100 kg)	333 €
PREPARATION FOR ATTITUDE ADD-ON BIKE - ONLY WITH FF-FRAME			
LFE090211	□ Preparation for ATTITUDE, w/o ATTITUDE order	(docking pins assembled)	
LFE090212	○ Preparation for ATTITUDE, with ATTITUDE order	(docking pins assembled; ATTITUDE order separately)	-
EMPULSE R20			
EMO000203	□ Empulse R20 added	(Please use separate OF to order the R20.)	-
EMO000204	□ Empulse R20 mounted	(Please use separate OF to order the R20.)	-
EMPULSE F55 BUNDLING			
EMO000101	○ With F55 order, not mounted to the wheelchair	(Please use separate OF to order the F55.)	
EMO000102	○ With F55 order, mounted to the wheelchair	(Please use separate OF to order the F55.)	
HANGERS & FOOTPLATES			
Aluminum hanger			-
LFE050002	○ Hanger angle 70°		-
LFE050001	○ Hanger angle 80°		-
LFE050041	○ Hanger angle 90°	(not with 6", 7" castor)	-
LFE050006	□ Hanger, aluminium, swing-away, LH	LFE050007 □ Hanger, aluminium, swing-away, RH	87 €
pick both	□ Hanger, aluminium, swing-away, Pair		174 €
Composite hanger			
LFE050002	○ Hanger angle 70°		-
LFE050001	○ Hanger angle 80°		-
LFE050003	□ Hanger, composit, swing-away, LH	LFE050004 □ Hanger, composit, swing-away, RH	-
pick both	□ Hanger, composit, swing-away, Pair		-
Elevating Legrest (ELR)			
LFE050009	□ Hanger, aluminium, flip-up (0° - 90°), swing-away, LH	(1)	154 €
LFE050010	□ Hanger, aluminium, flip-up (0° - 90°), swing-away, RH	(1)	154 €
pick both	Hanger, aluminium, flip-up (0° - 90°), swing-away, Pair	(1) (only from SW 38)	314 €
Amputee Support Pad			
LFE050014	□ Amputee support, left		216 €
LFE050013	□ Amputee support, right		216 €
LFE050012	□ Amputee support, pair		401 €
Footplate			
LFE050020	● Footboard, divided, composit, angle-adjustable, flip-up (sideways), heel loop, black (3)	(pair)	Std
LFE050050	□ Footboard, divided, aluminum, angle-adjustable, flip-up (sideways), heel loop, black	(from LLL 360 mm)	267 €
LFE050102	□ Footboard, platform, aluminum, angle-adjustable, flip-up (sideways), heel loop, black	(2)	197 €
LFE050126	□ Swing-away mechanism for platform footboard	(from LLL 380 mm)	197 €
LFE050062	□ Locking system, platform footplate, aluminium		68 €
Platform footrest, angle adjustable (only up to seat width 500 mm)			
LFE050131	□ Footboard, platform, lightweight, composite, angle & depth-adj. flip up (sideways), calf strap, black		172 €
LFE050132	□ Footboard, platform, lightweight, carbon fibre, angle & depth-adj. flip up (sideways), calf strap, black		333 €
LFE050127	□ Adjustable side protection, to create a side-stop on the footboard	(only in combination with Platform)	72 €
Autofolding (only up to seat width 500 mm)			
LFE050028	□ Footboard, platform, lightweight, composite, auto-folding, angle & depth-adj., calf strap, black (2)		415 €
LFE050027	□ Footboard, platform, lightweight, carbon fibre, auto-folding, angle & depth-adj., calf strap, black (2)		415 €
(1): Not in combination with platform & autofolding footplate			Fore more information www.SunriseMedical.eu © Sunrise Medical 2025
(2): Not in combination with ELR and SA frame			
(3): From SW 360 mm, with inset from SW 400mm			

FRAME ACCESSORIES



Transit wheels



Tip assist



Crutch holder



Anti-tip swing away

FRAME ACCESSORIES




Aluminum hanger, 90°



Composite hanger 70°/80°



Elevating legrest



Amputee Support Pad



Platform footplate lightweight angle & depth adjustable composite flip up



Footboard platform, lightweight, carbon fibre, auto folding



Divided footplate composite



Divided footplate aluminum



Platform, angle and depth adjustable; flip up, aluminium with locking system

F55 BUNDLING - NOGO MATRIX FOR Swing Away with 24" REAR WHEELS

COG	Seat Depth						
	380	400	420	440	460	480	500
40	✓	✓	✓	✓	✓	✓	✓
60	X	✓	✓	✓	✓	✓	✓
70	X	X	✓	✓	✓	✓	✓
80	X	X	✓	✓	✓	✓	✓
90	X	X	✓	✓	✓	✓	✓
100	X	X	X	✓	✓	✓	✓
110	X	X	X	✓	✓	✓	✓
120	X	X	X	X	✓	✓	✓
130	X	X	X	X	✓	✓	✓

F55 BUNDLING - NOGO MATRIX FOR Fixed Front with 24" REAR WHEELS

COG	Seat Depth						
	380	400	420	440	460	480	500
40	X	X	✓	✓	✓	✓	✓
60	X	X	X	✓	✓	✓	✓
70	X	X	X	✓	✓	✓	✓
80	X	X	X	X	✓	✓	✓
90	X	X	X	X	✓	✓	✓
100	X	X	X	X	X	✓	✓
110	X	X	X	X	X	✓	✓
120	X	X	X	X	X	X	✓
130	X	X	X	X	X	X	✓

F35 BUNDLING - NOGO MATRIX FOR Swing Away with 24" REAR WHEELS

COG	Seat Depth						
	380	400	420	440	460	480	500
40	✓	✓	✓	✓	✓	✓	✓
60	X	✓	✓	✓	✓	✓	✓
80	X	X	✓	✓	✓	✓	✓
90	X	X	✓	✓	✓	✓	✓
100	X	X	X	✓	✓	✓	✓
110	X	X	X	✓	✓	✓	✓
120	X	X	X	X	✓	✓	✓
130	X	X	X	X	✓	✓	✓

F35 BUNDLING - NOGO MATRIX FOR Fixed Front with 24" REAR WHEELS

COG	Seat Depth						
	380	400	420	440	460	480	500
40	X	X	X	X	X	X	X
60	X	X	X	X	X	X	X
70	X	X	X	X	X	X	X
80	X	X	X	X	X	X	X
90	X	X	X	X	X	X	X
100	X	X	X	X	X	X	X
110	X	X	X	X	X	X	X
120	X	X	X	X	X	X	X
130	X	X	X	X	X	X	X

Quickie Life			
<input checked="" type="radio"/> = Standard <input type="radio"/> = No cost option <input type="checkbox"/> = Cost option			€
LOWER LEG LENGTH - LLL in mm (Attention: LLL has to be 25 mm lower than seat height front)			
<input type="radio"/> # <input type="radio"/> 350 <input type="radio"/> 360 <input type="radio"/> 370 <input type="radio"/> 380 <input type="radio"/> 390 <input type="radio"/> # <input type="radio"/> 410 <input type="radio"/> 420 <input type="radio"/> 430 <input type="radio"/> 440 <input type="radio"/> 450 <input type="radio"/> # <input type="radio"/> 470 <input type="radio"/> 480 <input type="radio"/> 490 <input type="radio"/> 500			
FORK (compare with your choice on the first page)			
LFE080029	<input checked="" type="radio"/> Castor fork, Aluminium, black		Std
	LFE080042 <input type="checkbox"/> Colour Cast. Fork Silver		53 €
Castor fork alu powder coated			
	<input type="checkbox"/> 16 Matt Marine		53 €
	<input type="checkbox"/> 21 Red Ruby		53 €
	<input type="checkbox"/> 32 Signal White		53 €
	<input type="checkbox"/> 37 Marsred (Orange)		53 €
LFE080101	<input type="checkbox"/> Expansion kit, designed for hemiplegic users (Not in combination with 7" castors; SA frame only)		318 €
CASTOR (compare with your choice on the first page)			
LFE080301	<input type="radio"/> Castor, solid (hard), rim: composite, black		-
LFE080302	<input type="radio"/> Castor, solid (soft), rim: composite, black (only for 5" & 6" castors; max SH 520 mm)		-
LFE080305	<input type="radio"/> Castor, pneumatic, rim: composite, black (only for 6" & 7" castors)		-
LFE080303	<input type="checkbox"/> Castor, solid (hard), light-up, rim: composite, black (Not available for 6"/7")		76 €
Castor with aluminium rim			
LFE080304	<input type="checkbox"/> Castor, solid (hard), rim: aluminium, silver (Only for 4" and 5")		103 €
LFE080306	<input type="checkbox"/> Castor, solid (soft), rim: aluminium, silver (Only for 4" and 5")		155 €
	Colour, rim: <input type="checkbox"/> Black <input type="checkbox"/> Orange <input type="checkbox"/> Red <input type="checkbox"/> Blue		53 €
REAR WHEELS (compare with your choice on the first page; choose only one wheel type)			
LFE070000	<input type="checkbox"/> Without rear wheel		-
LFE070001	<input type="radio"/> Rear wheel, Universal, 36 spokes (silver), cross-wise spokes (only 22"/24")		-
LFE070002	<input type="radio"/> Rear wheel, Design, 36 spokes (silver), straight spokes (only 24")		-
LFE070003	<input type="checkbox"/> Rear wheel, Lightweight, 24 spokes (silver), straight spokes (max. user weight: 125 kg)		286 €
LFE070004	<input type="checkbox"/> Rear wheel, Proton, 24 spokes (black), straight spokes (max. user weight: 100 kg)		395 €
LFE070005	<input type="checkbox"/> Rear wheel, Spinergy, 18 spokes (black), straight spokes, hub: silver (only 24"/25")		712 €
LFE070008	<input type="checkbox"/> Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes (only 22"/24")		398 €
LFE070009	<input type="checkbox"/> Rear wheel, Mountainbike, 36 spokes (silver), straight spokes (only 24")		362 €
LFE070010	<input type="checkbox"/> Mountainbike wheel, replacement kit incl. tyre (extra coarse profile) (only 24")		548 €
One-hand operation (22" & 24"; only with stainless steel handrim, 0° camber)			
LFE070401	<input type="checkbox"/> Rear wheel, one-hand operation, 36 spokes, cross-spoked, left actuation		1.170 €
LFE070402	<input type="checkbox"/> Rear wheel, one-hand operation, 36 spokes, cross-spoked, right actuation		1.170 €
TYRES (choose only one tyre type)			
LFE070101	<input type="radio"/> Tyre, pneumatic, fine profile, Schwalbe Right Run (only 22"/24")		-
LFE070102	<input type="radio"/> Tyre, puncture proof, tread (25" only with Proton, lightweight wheel)		-
LFE070107	<input type="radio"/> Tyre, pneumatic, slick, Schwalbe One (only 24"/25")		55 €
LFE070104	<input type="checkbox"/> Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution (only 24"/25")		56 €
HANDRIM			
LFE070204	<input type="radio"/> Handrims, stainless steel, silver		-
LFE070202	<input type="radio"/> Handrims, aluminium, anodised, silver		-
LFE070201	<input type="radio"/> Handrims, aluminium, powder-coated, clear		-
LFE070203	<input type="checkbox"/> Handrims, aluminium, powder-coated, coloured		72 €
LFE070207	<input type="checkbox"/> Handrims, Supergrip ® (22" only with universal and drum brake wheels)		368 €
LFE070208	<input type="checkbox"/> Handrims, Maxgrepp ® (only 22"/24")		486 €
LFE070212	<input type="checkbox"/> Handrims, Ellipse 3R, black (ergonomic handrims with rubber-strip, 24" and 25" only)		386 €
LFE070200	<input type="radio"/> Without handrims		-
Handrim Mounting			
LFE070307	<input type="checkbox"/> Handrim mount, super narrow (7 mm) (Only with 24" stainless steel)		132 €
LFE070317	<input checked="" type="radio"/> Handrim mount, narrow (17 mm)		Std
LFE070327	<input type="radio"/> Handrim mount, wide (27 mm) (not with Ellipse 3R, Maxgrepp 24")		-
Handrim Cover			
LFE070250	<input type="checkbox"/> Handrim cover, Silicon (not with 25", Supergrip, Maxgrepp)		140 €
WHEEL CAMBER (choose only one option; Armrests only with 0° & 1°)			
LFE070600	<input checked="" type="radio"/> Camber 0°		Std
LFE070601	<input type="checkbox"/> Camber 1° LFE070603	<input type="checkbox"/> Camber 3°	41 €
LFE070602	<input type="checkbox"/> Camber 2° LFE070604	<input type="checkbox"/> Camber 4°	41 €
		For more information www.SunriseMedical.eu © Sunrise Medical 2025	

FORKS



Fork Aluminium,
black, orange & silver



Expansion kit, designed for
hemiplegic users

REAR WHEELS



Universal cross spaced



Design wheel



Lightweight wheel



Spinergy wheel, 18 spokes



Drum brake wheel

TYRES



Right Run Tyre



Solid Tyre



Marathon Plus
Evolution



Supergrip



Max Grepp



Ellipse 3R Handrim

WHEEL LOCKS



Knee lever wheel lock



Brake lever extension



Compact wheel lock



One arm wheel lock



Standard brake

ACCESSORIES



Caddy



Back pack



High pressure pump



Air pump



Mobile phone pocket

ANODIZED COLOUR AXLE PLATE



Axle plate
anodized, orange



Axle plate
anodized, silver



Axle plate
anodized, black

24 " LIGHT WEIGHT WHEEL COLOURS



Lightweight wheel, anodized orange



Lightweight wheel, anodized blue



Lightweight wheel, anodized red

Quickie Life				<input checked="" type="radio"/> = Standard <input type="radio"/> = No cost option <input type="checkbox"/> = Cost option		€	
SPOKE GUARDS (not with spinergy wheels, not with 25", not with mountainbike wheel)							
LFE070831	<input type="checkbox"/> Transparent					202 €	
		LFE070806	<input type="checkbox"/> Chrome			268 €	
LFE070822	<input type="checkbox"/> Pool	LFE070820	<input type="checkbox"/> Oriental			268 €	
LFE070811	<input type="checkbox"/> Feather	LFE070817	<input type="checkbox"/> Hot Wheels	LFE070802	<input type="checkbox"/> Basketball	268 €	
LFE070814	<input type="checkbox"/> Graffiti			LFE070827	<input type="checkbox"/> Silk	268 €	
LFE070828	<input type="checkbox"/> Space	LFE070812	<input type="checkbox"/> Flying Water	LFE070803	<input type="checkbox"/> Basketball 2	268 €	
WHEEL LOCKS (choose only one option)							
LFE060001	<input type="radio"/> Standard brake					-	
LFE060002	<input type="checkbox"/> Knee-lever brake					67 €	
LFE060010	<input type="checkbox"/> Brake lever extension, knee-lever brake, fold-down					32 €	
LFE060023	<input type="checkbox"/> Compact brake EVO					97 €	
LFE060015	<input type="checkbox"/> Push to lock brake, one-arm, left actuation					281 €	
LFE060016	<input type="checkbox"/> Push to lock brake, one-arm, right actuation					281 €	
ACCESSORIES							
LFE090000	<input checked="" type="radio"/> Tool kit					Std	
LFE090025	<input type="checkbox"/> Air pump					17 €	
LFE090024	<input type="checkbox"/> High pressure pump (0 - 11 bar)					92 €	
LFE090030	<input type="checkbox"/> Back pack					86 €	
LFE090026	<input type="checkbox"/> Mobile phone pocket					28 €	
COLOURS ⁽¹⁾ Please consult our colour catalogue to review your colour choice							
	Frame	Handrim	Hub/rim		Frame	Handrim	Hub/rim
	98 €	73 €	195 €		98 €	172 €	195 €
02 Brilliant Silver	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	27 Graphitemetallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
03 Silver Matt	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	29 Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
06 Sun Yellow	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	31 Matt Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
14 Blue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	32 Signal White	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
15 Brilliant Marine	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	33 Pink	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
16 Matt Marine	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	34 Dark Chocolate Metallic	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
19 Brilliant Purple	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	35 Skyblue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
21 Red Ruby	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	36 Cyan-Metallic	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22 Brilliant Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	37 Marsred (Orange)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
23 Matt Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	40 Kiwi	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24 Tutti Frutti Orange	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	41 Flamingo-Glitter	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
				42 Matt-Anthracite	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>
				46 Mamba Green (new)	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
STICKER COLOURS							
LFE100550	<input type="radio"/> Sticker colour - black						-
LFE100551	<input type="radio"/> Sticker colour - silver						-
ANODIZED COLOUR AXLE PLATE (not with one-hand operation)							
LFE100580	<input checked="" type="radio"/> Axle plate anodized, black						Std
LFE100581	<input type="checkbox"/> Axle plate anodized, silver						110 €
LFE100582	<input type="checkbox"/> Axle plate anodized, orange						110 €
24 " LIGHT WEIGHT WHEEL COLOURS							
	<input checked="" type="radio"/> Black						Std
LFE100336	<input type="checkbox"/> Anodized, orange						195 €
LFE100337	<input type="checkbox"/> Anodized, blue						195 €
LFE100338	<input type="checkbox"/> Anodized, red						195 €
(1) Handrim: only with handrim colour powder coated; Hub/rim: only for universal wheel 22" or 24"						For more information www.SunriseMedical.eu © Sunrise Medical 2025	