



Edition: 11.24  
Valid until: 11.25



Deliver to:

Name:

Street:

Town:

E-Mail:

Tel:

Postal code:

Fax:



Deliver to:

Name:

Street:

Town:

E-Mail:

Tel:

Postal code:

Fax:

Order date:

Purchase order:

Number required:

☐ Order ☐ Calculation ☐ Built 4 me (longer leadtime)

Marked for:

Sunrise order no:

☐ = Standard ☐ = Cost option  
☐ = Price reduction ☐ = No cost option

Max. user weight: 110 kg/140 Kg (Hybrid)

USER WEIGHT

kg

Mandatory information to ensure correct configuration

Art. No. DE: 76 800 000

Base Price

Xenon<sup>2</sup>

☐ = Standard ☐ = No cost option ☐ = Cost option

CLIENT MEASUREMENTS (all measurements in mm without cushion)

SEAT WIDTH SELECTION - measured from outside frame tube to outside frame tube (mm)

320	340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
									only Hybrid frame

SEAT DEPTH SELECTION - in mm

340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

FRONT SEAT HEIGHT SELECTION: in dependance of front castors, forks and frame versions

Fork (length)		Frame:		43	44	45	46	47	48	49	50	51	52	53	54	55
3"	111 mm	low		<input type="radio"/>	<input type="radio"/>	<input type="radio"/>										
	123 mm (1)					<input type="radio"/>	<input type="radio"/>									
	Frog Leg (1)							<input type="checkbox"/>	<input type="checkbox"/>							
	One Arm (1)			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>										
3"	111 mm	high						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>						
	123 mm (1)								<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	Frog Leg (1)										<input type="checkbox"/>	<input type="checkbox"/>				
	One Arm (1)							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>		
4"	111 mm	low			<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>								
	123 mm (1)					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
	Frog Leg (1)						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
	One Arm (1)			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
4"	111 mm	high						<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>					
	123 mm (1)								<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>				
	Frog Leg (1)										<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	One Arm (1)							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>	<input type="checkbox"/>	
5"	111 mm	low					<input type="radio"/>	<input type="radio"/>	<input type="radio"/>							
	123 mm (1)							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
	Frog Leg (1)								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
	One Arm (1)						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
5"	111 mm	high								<input type="radio"/>	<input type="radio"/>	<input type="radio"/>				
	123 mm (1)										<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Frog Leg (1)											<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	One Arm (1)									<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
6"	123 mm (1)								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
	Frog Leg (1)									<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	123 mm (1)										<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	Frog Leg (1)											<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

REAR SEAT HEIGHT SELECTION: in dependance of rear wheel size

Rear seat height in mm

Rear Wheels	370	380	390	400	410	420	430	440	450	460	470	480	490	500
24"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
25"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

(1) compatible with power add on devices.

(2) only with lightweight adapter

For more information  
visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 2025

# Xenon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## SEAT

### Seat sling

XFF020002	●	Seat sling, nylon, tension adjustable, 1 utility bag	
XFF020003	□	Seat sling, nylon, tension adjustable, 2 utility bags (only from SW 360 mm)	
XFF020004	□	Seat sling, nylon, tension adjustable, 1 utility bag, 1 catheter bag (only from SW 360 mm)	
XFF020000	○	Without seat sling	

### Seat cushion, cover with zipper, black

Foam (medium-soft)      Thickness:    XFF020101    □ 30 mm    XFF020102    □ 50 mm

### Seat accessories

XFF090018	□	Lap belt with buckle, metal	
XFF020050	□	Additional padding for seat upholstery, removable	

### JAY Cushions - for complete product range see Jay order forms

JAY020002	□	JAY Basic	
JAY020003	□	JAY Soft Combi P	
JAY020001	□	JAY Easy Visco, Curved Base	
JAY020004	□	JAY Lite	
JAY020007	□	JAY Easy Fluid, Curved Base	

## BACKREST

### Back height range

Choose back height factory setting (mm).

Height adjustments of in total 50 mm possible.

250 - 300	300 - 350	350 - 400	400 - 450	450 - 475
○ 250	○ 300	○ 350	○ 400	○ 450
○ 275	○ 325	○ 375	○ 425	○ 475
○ 300	○ 350	○ 400	○ 450	

### Backrest height, carbon fiber tubes (1)

Fixed height, non height adjustable

□ 250	□ 275	□ 300	□ 325
□ 350	□ 375	□ 400	□ 425

XFF030010	○	Backrest, fixed (90°)	
XFF030029	○	Desk sideguard preparation for later replacement (2)	
XFF030031	□	Backrest, fold-down (to the front) (4) (from back height 325 mm, fixed back height, not adjustable)	
XFF030012	□	Backrest, angle-adjustable (75° to 103°)	
XFF030029	○	Desk sideguard preparation for later replacement (3)	
XFF030031	□	Backrest, fold-down (to the front) (4) (from back height 325 mm, fixed back height, not adjustable)	

### Pushhandles

XFF030200	○	Without push handles	
XFF030201	○	Push handles long	
XFF030203	□	Push handles, fold-down (only available from back height 325 mm upwards)	
XFF030205	□	Push handles, height-adjustable for Jay backs (only available from back height 325 mm upwards)	
XFF030204	□	Push handles, height-adjustable (only available from back height 300 mm upwards)	

### Backrest upholstery

XFF030316	○	EXO Backrest sling, padded, adjustable, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black    ○ Orange    ○ Silver    ○ Red    ○ Blue	
XFF030317	□	EXO PRO Backrest, padded, adj., breathable spacer fabric, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black    ○ Orange    ○ Silver    ○ Red    ○ Blue	
XFF030300	○	Without backrest upholstery	


### Accessories for backrest

XFF030402	□	Stabilising bar, automatic folding (only available from SW 380mm; not with fixed height backrest)	
-----------	---	---	--

- (1) not with Jay J3 backs, not with Desk sideguard, not for use with add on bikes  
 (2) necessary for retrofitting Desk sideguard or Tubular armrest with Backrest, fixed 90°  
 (3) necessary for retrofitting Desk sideguard  
 (4) not with height adj; push handles, not with JAY backs; in combination with push handles fold-down only from backrest height 400 mm upwards; please see Folding Matrix on next side

For more information  
 visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

© Sunrise Medical 2025



XFF 020002 Seat sling with utility bag 1



XFF 020003 Seat sling with utility bags 2



XFF 030316 EXO backrest sling (EXO Evo version)



XFF 030317 EXO PRO backrest sling (EXO Evo version)



XFF 030201 Push handles long



XFF 030203 Push handles, fold-down



XFF 030204 Push handles, height-adjustable



XFF 090018 Lap belt with buckle, metal

### Folding-Matrix for XFF 030031 Backrest, fold-down to the front.



Seat width	Back height						
	325	350	375	400	425	450	475
320	✓	✓	✓	✓	✓	✓	✓
340	✓	✓	✓	✓	✓	✓	✓
360	✓	✓	✓	✓	✓	✓	✓
380		✓	✓	✓	✓	✓	✓
400			✓	✓	✓	✓	✓
420				✓	✓	✓	✓
440				✓	✓	✓	✓
460					✓	✓	✓
480						✓	✓
500							✓



=folds 90° to the front on folded chair



=can't be folded full 90° to the front on folded chair in combination with fold-down pushhandles



=can't be folded full 90° to the front on folded chair



=will fold less than 45° to the front on folded chair, folds 90° on unfolded chair

# Xenon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## Jay 3 Back - Shallow Contour (SC)

Mid Thoracic:	<input type="checkbox"/> JAY030010	Medium - 340 mm	<input type="checkbox"/> JAY030033	Tall - 380 mm		
Upper Thoracic:	<input type="checkbox"/> JAY030012	Short - 420 mm	<input type="checkbox"/> JAY030019	Medium - 460 mm	<input type="checkbox"/> JAY030025	Tall - 500 mm
Shoulder Height:	<input type="checkbox"/> JAY030015	Short - 530 mm	<input type="checkbox"/> JAY030022	Medium - 570 mm	<input type="checkbox"/> JAY030028	Tall - 610 mm

## Jay 3 Back - Mid Contour (MC)

Mid Thoracic:	<input type="checkbox"/> JAY030011	Medium - 340 mm	<input type="checkbox"/> JAY030018	Tall - 380 mm		
Upper Thoracic:	<input type="checkbox"/> JAY030013	Short - 420 mm	<input type="checkbox"/> JAY030020	Medium - 460 mm	<input type="checkbox"/> JAY030026	Tall - 500 mm
Shoulder Height:	<input type="checkbox"/> JAY030016	Short - 530 mm	<input type="checkbox"/> JAY030023	Medium - 570 mm	<input type="checkbox"/> JAY030029	Tall - 610 mm

## JAY J3 Back - Mid Deep Contour (MDC)

Mid Thoracic:	<input type="checkbox"/> JAY030044	Medium - 340 mm	<input type="checkbox"/> JAY030041	Tall - 380 mm	
Upper Thoracic:	<input type="checkbox"/> JAY030045	Short - 420 mm	<input type="checkbox"/> JAY030042	Medium - 460 mm	<input type="checkbox"/> JAY030046 Tall - 500 mm
Shoulder Height:	<input type="checkbox"/> JAY030043	Short - 530 mm	<input type="checkbox"/> JAY030047	Medium - 570 mm	<input type="checkbox"/> JAY030048 Tall - 610 mm

## Jay 3 Carbon Back - Shallow Contour (SC)

Lower Thoracic:	<input type="checkbox"/> JAY060001	Short - 170 mm	<input type="checkbox"/> JAY060002	Tall - 240 mm	
Mid Thoracic:	<input type="checkbox"/> JAY060003	Short - 300 mm	<input type="checkbox"/> JAY060004	Tall - 380 mm	

## Accessories for JAY J3 backs

WHT030415	<input type="checkbox"/> Whitmyer Headrest, COBRA XTRA, flip back	(only with JAY J3 backs)
WHT030416	<input type="checkbox"/> Whitmyer Headrest, AXYS, Standard	(only with JAY J3 backs)
JAY090058	<input type="checkbox"/> JAY Accessory bag	(only with Jay 3 Carbon backs)

## ARMREST AND SIDEGUARDS

### Armrest

XFF040005	<input type="checkbox"/> Desk sideguard, flip-up, removable	
XFF040001	<input type="checkbox"/> Desk sideguard, armpad (250 mm), fixed, flip-up, removable	
XFF040002	<input type="checkbox"/> Desk sideguard, armpad (250 mm), height-adjustable, flip-up, removable	
XFF040003	<input type="checkbox"/> Desk sideguard, armpad (310 mm), fixed, flip-up, removable	
XFF040004	<input type="checkbox"/> Desk sideguard, armpad (310 mm), height-adjustable, flip-up, removable	
XFF040052	<input type="checkbox"/> Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable	
XFF040117	<input type="checkbox"/> Tubular armrest, Aluminium, swing away (sideways), padded, removable	

### Accessories armrest

XFF040020	<input type="checkbox"/> Release assistance for tetraplegics	(only with desk sideguard)
-----------	--	----------------------------

### Sideguards

XFF040107	<input type="checkbox"/> Sideguard, composite, removable, 210 mm high (310 mm long), black	
XFF040101	<input type="checkbox"/> Sideguard, carbon fibre, cold weather protection, slim (black)	
XFF040102	<input type="checkbox"/> Sideguard, aluminium, cold weather protection (in frame colour)	
XFF040103	<input type="checkbox"/> Sideguard, aluminium, fender, cold weather protection (in frame colour)	
XFF040104	<input type="checkbox"/> Sideguard, carbon fibre, fender, cold weather protection (black)	

## FOOTBOARDS

XFF050023	<input type="checkbox"/> Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), heel loop, frame colour	
XFF050126	<input type="checkbox"/> Swing-away mechanism for platform footboard	(minimum LLL 30 mm longer)
XFF050059	<input type="checkbox"/> Footboard, platform, Performance, aluminium, angle-adjustable, flip-up (sideways), calf strap	(1)
XFF050131	<input type="checkbox"/> Footboard, platform, lightweight, plastic, angle & depth-adj. flip up (sideways), calf strap, black	
XFF050132	<input type="checkbox"/> Footboard, platform, lightweight, carbon fibre, angle & depth-adj. flip up (sideways), calf strap, black	
XFF050061	<input type="checkbox"/> Footboard, platform, lightweight, plastic, angle & depth-adj. flip up (sideways), high-mount black	(LLL 220-340 mm)

### Footboard, auto-folding

XFF050027	<input type="checkbox"/> Footboard, platform, lightweight, carbon fibre, auto-folding, angle & depth-adjustable, calf strap, black	
XFF050028	<input type="checkbox"/> Footboard, platform, lightweight, plastic, auto-folding, angle and depth-adjustable, calf strap, black	

### Accessories for footboards

XFF050127	<input type="checkbox"/> Adjustable side protection, to create a side-stop on the footboard	(2)
-----------	---	-----

## LOWER LEG LENGTH (LLL): has to be 25 mm less than front seat height; all measurements are in mm

<input type="radio"/> 220	<input type="radio"/> 230	<input type="radio"/> 240	<input type="radio"/> 250	<input type="radio"/> 260	<input type="radio"/> 270	<input type="radio"/> 280	<input type="radio"/> 290	<input type="radio"/> 300	<input type="radio"/> 310
<input type="radio"/> 320	<input type="radio"/> 330	<input type="radio"/> 340	<input type="radio"/> 350	<input type="radio"/> 360	<input type="radio"/> 370	<input type="radio"/> 380	<input type="radio"/> 390	<input type="radio"/> 400	<input type="radio"/> 410
<input type="radio"/> 420	<input type="radio"/> 430	<input type="radio"/> 440	<input type="radio"/> 450	<input type="radio"/> 460	<input type="radio"/> 470	<input type="radio"/> 480	<input type="radio"/> 490	<input type="radio"/> 500	

(1) not in SW320 with 10 mm inset

(2) only for footboard, platform lightweight or Performance

For more information  
visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 2025

JAY J3 Shallow Contour



JAY J3 Mid Contour



JAY J3 Mid Deep Contour



XFF040107  
Composite  
sideguard  
detachable



XFF040101  
carbon fibre, cold  
weather protection,  
slim



XFF040102  
Aluminium, cold weather  
protection



XFF040103  
Aluminium, cold weather  
protection, fender



XFF04000x  
Desk sideguard, armpad, flip-  
up, removable



XFF040104  
Carbon fibre, cold weather  
protection, fender



XFF050132  
Footboard platform,  
lightweight, carbon fibre



XFF050059  
Footboard platform,  
performance



XFF050027  
Footboard platform,  
lightweight, carbon fibre,  
auto folding



XFF050127  
Adjustable side protection

# Xenon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## FRAME VERSIONS

### Open frame (max. user weight 110kg)

XFF010012	○	Frame FF 100 ° (integrated footrest) (UK 80 °)	
		XFF010101    □ Inset: Slight (10 mm)	
XFF010013	□	Frame FF 92 ° (integrated footrest), performance (UK 88 °)	
		XFF010101    □ Inset: Slight (10 mm)	

### Hybrid frame (max user weight 140kg)

XFF010020	□	Frame FF 100 ° (integrated footrest), hybrid	
XFF010021	□	Frame FF 92 ° (integrated footrest), hybrid	

## CENTRE OF GRAVITY (COG)

XFF070703	Adapter, standard	○ 2    ○ 3    ○ 4    ○ 5    ○ 6 ○ 7    ○ 8    ○ 9    ○ 10	
XFF070701	Adapter, lightweight (1)	○ 5    ○ 6    ○ 7    ○ 9    ○ 10    ○ 11	
XFF010411	○	Axle plate, aluminium, long (full adjustment)	
XFF010412	□	Axle plate, aluminium, shorted (max. adjustment 10 mm)	

## FRAME ACCESSORIES

(choose either tip assist or anti- tip tubes per side)

XFF090016	□	Frame protection, black	
XFF090100	□	Tie down brackets for transportation	
XFF090010	□	Transit wheels	
XFF090014	□	Frame reinforcement kit for add on bike	
XFF090050	□	Anti-tip tubes, plug-in, pair	
<b>Anti-tip tubes, swing-away</b>			
XFF090004	□	Left    XFF090005    □ Right	
XFF090006	□	Pair	
<b>Tip assist</b>			
XFF090001	□	Left    XFF090002    □ Right	
XFF090003	□	Pair	
XFF090009	□	Crutch holder, loop (only in combination with tip assist)	

## PREPARATION FOR ATTITUDE ADD-ON BIKE

XFF090211	□	without ATTITUDE order	Please use & include Attitude OF
XFF090212	○	with ATTITUDE order	Please use & include Attitude OF

## EMPULSE F55 BUNDLING

EMO000101	○	With F55 order, not mounted to the wheelchair	(Please use separate OF to order the F55. Not with Lightweight Compact brake)
EMO000102	○	With F55 order, mounted to the wheelchair	(Please use separate OF to order the F55. Not with Lightweight Compact brake)

## EMPULSE F35 BUNDLING

EMO000301	○	With F35 order, not mounted to the wheelchair	(Please use separate OF to order the F35.)
EMO000302	○	With F35 order, mounted to the wheelchair	(Please use separate OF to order the F35.)

## FORKS

(compare with your choice on the first page)

XFF080013	●	Castor fork, Carbotecture ®    Colour: Black	
		Colour:    □ Silver    □ Orange    □ Red (2)    □ Blue (2)	
XFF080029	○	Castor fork, Aluminium, 123mm    Colour: Black	
		Colour:    □ Silver    □ Orange	
XFF080008	□	Castor fork, aluminium, one-arm    Colour: Black	
		Colour:    □ Silver    □ Orange    □ Red    □ Blue	
XFF080007	□	Castor fork, aluminium, suspension (Frog Legs ®), black	

## CASTORS

(compare with your choice on the first page)

XFF080301	○	Castor, solid (hard), rim: composite, black	
XFF080302	○	Castor, solid (soft), rim: composite, black	(only 5" & 6")
XFF080303	□	Castor, solid (hard), light-up, rim: composite, black	(only 4" & 5")
XFF080304	□	Castor, solid (hard), rim: aluminium    Colour, rim: Silver	(only 4" & 5")
XFF080306	□	Castor, solid (soft), rim: aluminium    Colour, rim: Silver	(only 4" & 5")
		Colour, rim:    □ Black    □ Orange    □ Red    □ Blue	

(1) only with 0° wheel camber

(2) only in combination with castor, soft (XFF080302 and XFF080306)

For more information  
visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
© Sunrise Medical 2025





XFF010012 / 13  
Open Frame



XFF090020 / 21  
Hybrid frame



XFF090001  
Tip assist



XFF090004  
Anti-tip tubes, swing away



XFF090050  
Anti-tip tubes, plug in, pair



XFF090010  
Transit wheels

#### EMPULSE F35 & F55 BUNDLING - NOGO MATRIX FOR 24" REAR WHEELS

COG	Seat Depth						
	320-380	400	420	440	460	480	500
20	✓	✓	✓	✓	✓	✓	✓
30	✓	✓	✓	✓	✓	✓	✓
40	✓	✓	✓	✓	✓	✓	✓
50	X	✓	✓	✓	✓	✓	✓
60	X	✓	✓	✓	✓	✓	✓
70	X	X	✓	✓	✓	✓	✓
80	X	X	✓	✓	✓	✓	✓
90	X	X	X	✓	✓	✓	✓
100	X	X	X	✓	✓	✓	✓
110	X	X	X	X	✓	✓	✓

#### EMPULSE F35 & F55 BUNDLING - NOGO MATRIX FOR 25" REAR WHEELS

COG	Seat Depth						
	320-380	400	420	440	460	480	500
20	✓	✓	✓	✓	✓	✓	✓
30	X	✓	✓	✓	✓	✓	✓
40	X	✓	✓	✓	✓	✓	✓
50	X	X	✓	✓	✓	✓	✓
60	X	X	✓	✓	✓	✓	✓
70	X	X	X	✓	✓	✓	✓
80	X	X	X	✓	✓	✓	✓
90	X	X	X	X	✓	✓	✓
100	X	X	X	X	✓	✓	✓
110	X	X	X	X	X	✓	✓

# Xenon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## REAR WHEELS

(compare with your choice on the first page; choose only one wheel type)

XFF070000	○	Without rear wheel	
XFF070001	○	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes	(only 24")
XFF070002	○	Rear wheel, Design, 36 spokes (silver), straight spokes	(only 24")
XFF070003	○	Rear wheel, Lightweight, 24 spokes (silver), straight spokes	(please see colour options on last page)
XFF070004	□	Rear wheel, Proton, 24 spokes (black), straight spokes	(user weight limit 100kg)
XFF070005	□	Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver	
XFF070007	□	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub black	
XFF070008	□	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes	(1) (only 24")
XFF070009	□	Rear wheel, Mountainbike, 36 spokes (silver), straight spokes	(only 24")
XFF070010	□	Mountainbike wheel, replacement kit incl. tyre (extra coarse profile)	(only 24")

## Axle

XFF070150	●	Quick-release axle	
XFF070151	□	Quick-release axle, designed for quads	(2)

## TYRES

(choose only one tyre type)

XFF070101	○	Tyre, pneumatic, fine profile, Schwalbe Right Run	(only 24")
XFF070102	○	Tyre, puncture proof, tread	
XFF070107	□	Tyre, pneumatic, slick, Schwalbe One	
XFF070104	□	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	
XFF090163	□	Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube	(not with solid tyres)

## HANDRIMS

XFF070204	○	Handrims, stainless steel, silver	
XFF070202	○	Handrims, aluminium, anodised, silver	
XFF070201	○	Handrims, aluminium, powder-coated, clear	
XFF070203	□	Handrims, aluminium, powder-coated, coloured	(select colour on last page)
XFF070205	□	Handrims, titanium, untreated	
XFF070206	□	Handrims, titanium, powder-coated, clear	
XFF070207	□	Handrims, Supergrip®	
XFF070208	□	Handrims, Maxgrepp®	(only 24")
XFF070212	□	Handrims, Ellipse 3R, black	(ergonomic handrims with rubber-strip)

## Hand rim mounting

XFF070307	□	Handrim mount, super narrow (7 mm)	(only with stainless steel)
XFF070317	●	Handrim mount, narrow (17 mm)	
XFF070327	○	Handrim mount, wide (27 mm)	

## Handrim accessories

XFF070250	□	Handrim cover, rubber	(only 24")
-----------	---	-----------------------	------------

## WHEEL CAMBER

(choose only one option)

XFF070600	○	0°	
XFF070602	○	2°	(not with Centre of Gravity Adapter, lightweight)
XFF070604	○	4°	(not with Centre of Gravity Adapter, lightweight)

## SPOKE GUARDS

(only 24")

XFF070831	□	Transparent	
	XFF070813	□	Galactica
XFF070823	□	Princess	XFF070818 □ Jungle Fever
XFF070816	□	Happy Lions	XFF070811 □ Feather
XFF070819	□	Monkeys	XFF070821 □ Pirates
XFF070820	□	Oriental	XFF070806 □ Chrome
XFF070829	□	Superhero	XFF070814 □ Graffiti
XFF070817	□	Hot Wheels	XFF070802 □ Basketball XFF070803 □ Basketball 2
XFF070822	□	Pool	XFF070812 □ Flying Water
XFF070827	□	Silk	XFF070828 □ Space
XFF070830	□	Tiger	XFF070833 □ Zebra

- (1) only with Adapter, standard and only with 0° & 4° camber; not with push handles, fold down & not without push handles.  
 (2) not with Rear wheels Spinergy/ Mountainbike/ drum brakes

For more information  
 visit our webpage  
[www.SunriseMedical.eu](http://www.SunriseMedical.eu)  
 © Sunrise Medical 2025





XFF070001  
Universal wheel



XFF070002  
Design wheel



XFF070003  
Lightweight wheel



XFF07000  
Proton wheel



XFF07010  
Mountainbike wheel



XFF070101  
Schwalbe Right Run



XFF070107  
Schwalbe One



XFF070104  
Schwalbe Marathon Plus  
Evolution



XFF070102  
Tyre, puncture proof



XFF070202  
Handrims aluminum,  
anodized



XFF070203  
Handrims aluminum,  
powder coated, coloured



XFF070208  
Handrims Maxgrepp



XFF070212  
Handrims Ellipse 3R

# Xenon<sup>2</sup>

● = Standard    ○ = No cost option    □ = Cost option

Price in €

## WHEEL LOCKS

(choose only one option)

XFF060001	● Standard brake	
XFF060002	□ Knee-lever brake	
XFF060010	□ Brake lever extension, knee-lever brake, fold-down	
XFF060023	□ Compact brake EVO	
XFF060024	□ Lightweight compact brake EVO (1)	

## ACCESSORIES

XFF090000	● Tool kit	
XFF090025	□ Air pump	
XFF090024	□ High pressure pump (0 - 11 bar)	
XFF090030	□ Back pack	
XFF090026	□ Mobile phone pocket	

## COLOURS<sup>(2)</sup>

Please consult our colour catalogue to review your colour choice

	Frame	Handrim	Hub/rim	X-brace		Frame	Handrim	Hub/rim	X-brace	
02 Brilliant Silver	○	□	□	□	27 Graphitemetallic	○	□	□	□	
03 Silver Matt	□	□	□	□	29 Black	○	□	□	□	
06 Sun Yellow	○	□	□	□	31 Matt Black	○	□	□	□	
14 Blue	○	□	□	□	32 Signal White	○	□	□	□	
15 Brilliant Marine	○	□	□	□	33 Pink	○	□	□	□	
16 Matt Marine	□	□	□	□	34 Dark Chocolate Metallic	□	□	□	□	
19 Brilliant Purple	○	□	□	□	35 Skyblue	□	□	□	□	
21 Red Ruby	○	□	□	□	36 Cyan-Metallic	□	□	□	□	
22 Brilliant Red	□	□	□	□	37 Marsred (Orange)	□	□	□	□	
23 Matt Red	□	□	□	□	40 Kiwi	□	□	□	□	
24 Tutti Frutti Orange	○	□	□	□	41 Flamingo-Glitter	□	□	□	□	
					42 Matt-Anthracite	○	□	□	□	
					46 Mamba Green (new)	□				
					Rear wheel, lightweight anodized, orange				□	
					Rear wheel, lightweight anodized, blue				□	
					Rear wheel, lightweight anodized, red				□	

## COLOURS - for stickers ("Xenon<sup>2</sup>" on the frame and "Quickie" on sideguard)

XFF100550	○ black	(Sticker "Quickie" on sideguard will be silver in combination with carbon sideguards)	
XFF100551	○ silver		
XFF100552	○ orange		
XFF100553	○ red		
XFF100554	○ blue		

(1) not with Hybrid frame; not available with 4° camber; not with seat depth 340 mm; not for use with add on bikes; adjustment of CoG +/- approx. 25mm in the field

(2) Handrim: only with handrim colour powder coated; Hub/rim: only for universal wheel 24"

(3) only with 24" Lightweight wheel, Hub and rim anodized, black spokes

Comments for special build - please apply for a quote; lead time conditions my change.



## Quickie Xenon<sup>2</sup> - Basic Price

Option price	..... €
Packaging and transport	..... €
Exclusive VAT	..... €
VAT.	..... €
<b>Total Price</b>	..... €

Sunrise Medical Customer service

Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg

Tel. +497253 980-0

Fax. +497253 980-222

E-mail: [export@sunrisemedical.de](mailto:export@sunrisemedical.de)

Website: [sunrisemedical.eu](http://sunrisemedical.eu)

For more information

visit our webpage

[www.SunriseMedical.eu](http://www.SunriseMedical.eu)

© Sunrise Medical 2025



XFF060001  
Standard brake



XFF060002  
Knee lever brake



XFF060023  
Compact brake EVO



XFF060024  
Lightweight compact  
brake EVO



XFF090000  
Tool kit



XFF090024  
High pressure pump



XFF090025  
Air pump



XFF090030  
Back pack



XFF090026  
Mobile phone pocket