TECHNICAL DATA

Max. user weight:	65 kg
Seat width:	22 – 36 cm (growth +2 cm without additional parts)
Seat depth:	24 – 38 cm (growth +4 cm without additional parts)
Front seat height:	36 – 48 cm
Rear seat height:	29 – 46 cm
Backrest height:	20 – 40 cm
Backrest angle:	-16° to +20°
Frame versions:	straight or abduction frame (3 cm each side)
	straight of abaaction name to one cach stacy
Frame angle:	95°
Frame angle:	95°
Frame angle: Front wheel	95° 3"/4"/5"/6"
Frame angle: Front wheel Rear wheel	95° 3"/4"/5"/6" 20"/22"/24"
Frame angle: Front wheel Rear wheel Wheelchair weight:	95° 3"/4"/5"/6" 20"/22"/24" from 8 kg





For further information on product options please check the order form on our website. You can view this document also in a larger font on a PDF-file at www.SunriseMedical.eu





Sunrise Medical GmbH & Co. KG

Kahlbachring 2-4 D-69254 Malsch/HD Tel.: +49 (0) 72 53 / 9 80-0 Fax: +49 (0) 72 53 / 9 80-222 E-Mail: kundenservice@sunrisemedical.de www.SunriseMedical.de



Generation 2015

A lightweight wheelchair with growth adjustability

ΤМ

f (





(

One thing we're certain of is that kids will grow! That's why the Simba features 4 cm of depth growth built right into the frame. The 2 cm width growth capability works also without additional parts. The Simba – your loyal friend!

6 upholstery colours

A DESIGN THAT'S FIT FOR KIDS

With independence comes confidence. Every centimeter of the Simba was designed to give kids the best possible wheel access for efficient propulsion. With features such as the abduction frame, flip-up angle adjustable footrest and safari wheel lock integrated into side guard, kids are able to independently access the wheels, transfer in and out of the wheelchair, and access the wheel locks without interference during propulsion.





IT'S ALL ABOUT ERGONOMICS

The push bar can be easily adjusted in angle and height for an enhanced ergonomic position. The new stabilizer reinforces the backrest especially when push handles are chosen.

A FULLY ADJUSTABLE AXLE ADAPTER

The optimised axle adapter offers full adjustability of the rear seat height, centre of gravity and camber without compromising on any performace.



THE CONVINCING SAFARI WHEEL LOCK

With its improved handling the updated safari side guard wheel lock let you stop quickly and safely.

COLOURFUL OPTIONS

We want the Simba to be an expression of your child's unique personality. Aside from its super stylish frame, the Simba also offers 28 bright frame colours, 6 different upholstery colours, coloured castor forks and a full range of stylish spoke guards.



SIMBA – AN INNOVATIVE DESIGN

The outstanding design of the new castor link tube, combined with the innovative carbotecture forks that are made from extremely lightweight stiff high-tech material, makes the wheelchair just that little bit lighter.

