

Patient Lifts

Sunlift




Sunlift



The SUNLIFT range consists of 4 models: MICRO, MINI, MIDI and MAJOR. All models feature exceptional rigidity thanks to the advanced design of all of the moving connections, eliminating play and offering a rigid feel with the hoist fully loaded. Actuators incorporate a double security button, light indicator of the battery charge level and an acoustic warning before total discharge happens. The battery capacity (2.9 Amph) provides a range of use at 50% higher than usual, allowing up to 60 double uses (up and down).

Sunlift Range

Measurements





		SUNLIFT Micro	SUNLIFT Mini
9	Safe working load	130 kg	130 kg
	Max total length with wheels facing outwards)	1.010 mm	1.130 mm
	Min total length with wheels facing inwards)	970 mm	1.065 mm
1	Max widht with legs closed	545 mm	545 mm
1	Max height (without load) (A)	1.725 mm	1.870 mm
1	Max height (without load) (B)	1.130 mm	1.280 mm
1	Max spreader bar height (without load) (C)	1.483 mm	1.628 mm
1	Min spreader bar height (without load)	460 mm	513 mm
1	Turning radius with legs open (D)	1.130 mm	1.240 mm
N	Legs open: (E) Max widht (external) Max widht (internal)	870 mm 740 mm	940 mm 810 mm
N	Legs closed: (F) Max widht (external) Max width (internal)	545 mm 465 mm	540 mm 470 mm
1	Total leg height	100 mm	100 mm
(Ground clearance	20 mm	20 mm
(Castor wheels	75 mm	75 mm
F	Rear wheels	100 mm	100 mm
	Mast, boom and electric components assembly weight	16,6 kg	17,54 kg
E	Base assembly weight	13,34 kg	13,6 kg
E	Batteries weight	2,76 kg	2,76 kg







Technical Features (common for the 4 models)

Batteries: 2 x 12v de gel Batteries capacity: 2,9 Amph Charger input: 230 v / 50 Hz 0.12 A Charger output: 28 v 0.5 Amph

SUNLIFT Midi	SUNLIFT Major
150 kg	175 kg
1.225 mm	1.340 mm
1.160 mm	1.275 mm
640 mm	650 mm
1.940 mm	2.006 mm
1.285 mm	1.435 mm
1.688 mm	1.813 mm
555 mm	575 mm
1.380 mm	1.500 mm
1.075 mm 980 mm	1.185 mm 1.105 mm
640 mm 540 mm	650 mm 520 mm
100 mm	110 mm
20 mm	15 mm
75 mm	75 mm
100 mm	100 mm
17,6 kg	22,75 kg
15,1 kg	18,1 kg
2,76 kg	2,76 kg

Common features



Legs can be opened and closed by pressing the ergnomic pedal.



Hand control classified with IPX4 (Water resistant) and ergonomic design.



Ergonomic push handles allowing elbows position at 90° for easier handling.



Slings

Slings



	Quick fit sling	Quick fit sling deluxe	Access sling
Recommended use	85 % of patients	95% of patients	For easier access to toilet
Material	Polyester, Padded or Net	Padded	Polyester
Sizes	Medium or large	Medium	Medium
Independece level *	1 and 2	1 and 2	2 and 3
Others	-	With or wihout head support	_

	Long seat sling	
Recommended use	Amputee patients	
Material	Polyester	
Sizes	Medium	
Independece level*	1 and 2	
Others	-	

Selection of the most appropriate sling for each case:

		_				
	Toilet	Commode chair to bath	From chair to chair	From bed to chair	From floor to bed	From chair to commode chair
Quick fit sling					-	
Quick fit sling deluxe						
Access sling						
Long seat sling				-	-	

*Level of independence

1 = Total independece

7 = Total independece



All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text

Sunrise Medical S.L

Polígono Bakiola 41 48498 Arrankudiaga Vizcaya, España

Tel: 902 14 24 34 Fax: 94 648 15 75 www.SunriseMedical.es