

Wheel Add-On

 | m90



Rev:01/2026
Valid:12/2026

ORDER FORM

Errors and omissions excepted.

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address
C.-ID:			C.-ID:		
Name:			Name:		
Street:			Street:		
Town/PC:			Town/PC:		
E-Mail:			E-Mail:		
Phone:			Phone:		
Order Information			User Weight		
Order	Calculation	Built-4-Me®	Max. user weight: 120 kg		
Purchase Order No.:			User Weight: kg		
Art. No.: 99200001			Reimbursement Code: xxx		
Order Date:	Day	Month	Year		

Base Price € X.XXX,-

Technical Data Overview

Max. user weight	120 kg (265 lbs)	Motor	150W
Rear Wheel Sizes	24" (540mm)	Operating Voltage	24V
Total Device Weight	12 kg (26.4 lbs)	Maximum Motor Noise	60db
Per Battery Weight (4.0Ah)	0.6 kg (1.3 lbs)	Battery Type	Lithium Ion (21.6V DC and 86.4Wh each battery)
Max. Range Per Charge	15 km (9.3 mi)	Device Charger Type	28V DC 5.7A 160W
Maximum Speed	6 kph (3.7 mph)	Docking Charger	28V DC 1.5A IP44
Maximum Slope	6°	Control Box	3.7V 1Ah 3.7Wh
Bluetooth Standard	BLE (125Kbps Coded PHY)		

Base

€

M90 unit	Incl. Dual Controller box, Joystick and Speed Controller with Handrim Steering, Batteries, long and short mounting hardware, Quick Release Axles, Cable Charger, Docking Charger	–
-----------------	--	---

Bundling

€

EMM090213	Without wheelchair order	–
EMM090217	With wheelchair order, not mounted Wheelchair order on separate order form. Seperate delivery of the wheelchair and the M90.	–
EMM090218	With wheelchair order, mounted to wheelchair Wheelchair order on separate order form.	–

Sunrise Medical Wheelchair Models - Specific mounting hardware will be provided for the chair selected

€

Specific wheelchair model the M90 Wheels will be used on

EMM000035	Parix ² Manual wheels with 12mm axle not retrofittable.	EMM000184	Neon ²
EMM000036	Basix ² Manual wheels with 12mm axle not retrofittable.	EMM000052	Ibis
EMM000037	Rubix ² Manual wheels with 12mm axle not retrofittable.	EMM000075	SR45
EMM000038	Helix ²	EMM000065	Xenon ²
EMM000055	Easy Life Eco	EMM000076	IRIS
EMM000056	Easy Life i	EMM000103	Argon ²
EMM000102	Easy Life T	EMM000174	QS5 X
EMM000111	Easy Life RT	EMM000101	Helium or Nitrum
EMM000059	Style X / Style X Ultra	EMM000182	Zippie IRIS
EMM000108	Easy Life	EMM000183	SIMBA
EMM000109	Easy Life R	EMM000149	Youngster 3
EMM000129	Easy 200/300		

RGK - Empulse M90 is compatible with RGK models. Contact your local RGK representative for more information.

Controller configuration: Only for wheelchair orders, factory pre-assembly

€

EMM110053	Joystick mounting location, left	–
EMM110052	Joystick mounting location, right	–
EMM110214	Mounting bracket, long	–
EMM110213	Mounting bracket, short	–

Handrims		€
EMM070096	Handrims, aluminium, powder coated, black Mounting position, wide (27mm)	–
EMM070207	Handrims, Supergrip R Mounting position, wide (27mm)	XXX,-
EMM070212	Handrims, Ellipse, 3R, black Mounting position, narrow (17mm)	XXX,-
FMG-MF06-3247-01	Handrims, Gekko Mounting position, narrow (17mm)	XXX,-
FMG-MF06-3247-03	Handrims, Quadro Mounting position, narrow (17mm)	XXX,-

Tyres		€
EMM070104	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	–
EMM070102	Solid tyre, puncture proof, tread	–

Options and Accessories - Batteries are interchangeable and can be used with either wheel		€
EMM090227	Leg Mount Controller Includes thigh strap and hardware	XXX,-
EMM160042	Additional Li Ion Battery Set 2 additional batteries (21.6V DC and 86.4Wh each battery)	XXX,-
EMM090233	Wrist Support, Joystick	XXX,-
EMM090232	Swing Away Mount, Joystick, right available soon	XXX,-
EMM090235	Swing Away Mount, Joystick, left available soon	XXX,-

Comments

M90

€

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day Month Year

**Sunrise Medical GmbH**

Kahlbachring 2-4
D-69254 Malsch/HD
Phone: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: Export@sunrisemedical.de
www.SunriseMedical.eu