



NEATECH™

—
JoB SERIES



Easy to transport



Foldable



Lightweight





NEATECH™

Free to be disabled with the JoB the holiday is re-discovered

JoB responds to the need of people with disabilities to have a specific product for the holidays and especially for bathing.

The Neatech has created and produced "J'amme 'o bagn" which is not only eye-catching invitation in Neapolitan to bathe in the sea, instead the JoB It has effectively solved the problem of the inability of people with disabilities or temporary limitations to enter the water, thus satisfying the natural desire of sea and atmosphere of fun and lightheartedness that characterize a vacation day. The vacation of a person with a disability should not be conceived as a simple moment of escape and entertainment, but also as an opportunity for growth.

Thanks to its waterproof wheels with variable section and controlled deformation, the JoB enables users and their companions to face pebbles, gravel, snow and dirt, transforming each new obstacle into a source of entertainment. To facilitate the transport has been studied a tightening system with no fixed restrictions, so the wheelchair is assembled and dismantled in a few minutes and without the help of any tools.

JoB It is the ideal solution for swimming pools and bathing facilities not yet equipped for the stay of the disabled, as it allows you to adapt the structure to the rules in force for the removal of architectural barriers.

FEATURES



Built with sea water resisting materials



Demountable with quick-release systems



Enters the water



Ergonomic and comfortable



Wide range of accessories




Easy to transport even by airplane

TECHNICAL SPECIFICATIONS

Total length	1250 mm
Total width	900 mm
Minimum transport length	1020 mm
Seat height	360 mm
Total weight	7 kg + 3kg x wheel
Maximum user weight	125 kg
Seat width	465 mm
Seat depth	420 mm
Backrest width	465 mm
Backrest height	710 mm
Backrest angle-adjustable	fixed
Legrest angle-adjustable	fixed
Warranty	2 years
EN12182 CE	yes

AVAILABLE COLOURS



www.neatech.it

NEATECH.it s.r.l.
Via A. De Curtis 4/a
80040 Cercola (NA) ITALIA
tel 081 555.19.46
fax 081 555.25.07
e-mail info@neatech.it