



JAY LitePro

JAY LitePro OF EU Template 21.01.2026

Invoice to:
Name:
Street:
Town: Postal code:
E-Mail:
Tel: Fax:

Deliver to:
Name:
Street:
Town: Postal code:
E-Mail:
Tel: Fax:



Order date:
Purchase order:
Number required:

Marked for:
Sunrise order no:

- ☐ Order
☐ Calculation

JAY® LitePro Cushion Standard Features

Multizone Elastomer Insert with Pressure Adjustable pods for superior stability and pressure management
Perforated multi-density foam base for improved moisture and airflow.
Washable microclimate outer cover and inserts
Low profile 50mm cushion height with only 0.9Kg cushion weight
2 Year guarantee, 6 months on cover

Max. user weight: 150 Kg

USER WEIGHT kg

Mandatory information to ensure correct configuration

Reimbursement code:

Application:

Client at high risk of skin breakdown.

Price -. XXX€

To order a JAY LitePro cushion, follow the steps below to build your part number.

fixed	Width	Depth	Insert	Outer Cover	fixed
JLP	40	42	S	M	EU
JAY LitePro Cushion	STEP 1	STEP 2	STEP 3		

Example: JLP4042SMEU

CUSHION WIDTH AND DEPTH - STEP 1 + STEP 2


Standard Cushion (width: 340 - 500 mm; max user weight 150 kg)

Width →

mm	340	360	380	400	420	440	460	480	500
340									
360	◇	○							
380	○	○	○	○	◇				
400	○	○	○	○	○	○	○	○	
420	○	○	○	○	○	○	○	○	◇
440		○	○	○	○	○	○	○	◇
460			○	○	○	○	○	○	○
480						○	○	○	○
500							○	○	○
520									
600									

Insert Size

JAY LitePro Measuring Guide



Width

Xmm = Width

Measurement made from outside edge to outside edge of cushion. (Including cover).

Depth

Ymm = Depth

Measurement made from back to front topside. (or bottom of contracture cut) of cushion. (Including Cover)

WIDTH

STEP 1

DEPTH

STEP 2

Lead Time

○ 5 Days

◇ 15 Days



For more information
visit our webpage

www.sunrisemedical.com

© Sunrise Medical 2026



Lead Time Overview

Standard Insert -5 Day Lead Time

Seat Width	Seat Depth	JAY LitePro Part No.
340	380	JLP3438SMEU
	400	JLP3440SMEU
	420	JLP3442SMEU
360	360	JLP3636SMEU
	380	JLP3638SMEU
	400	JLP3640SMEU
	420	JLP3642SMEU
380	440	JLP3644SMEU
	380	JLP3838SMEU
	400	JLP3840SMEU
	420	JLP3842SMEU
400	440	JLP3844SMEU
	460	JLP3846SMEU
	380	JLP4038SMEU
	400	JLP4040SMEU
420	420	JLP4042SMEU
	440	JLP4044SMEU
	460	JLP4046SMEU
	400	JLP4240SMEU
440	420	JLP4242SMEU
	440	JLP4244SMEU
	460	JLP4246SMEU
	400	JLP4440SMEU
460	420	JLP4442SMEU
	440	JLP4444SMEU
	460	JLP4446SMEU
	480	JLP4448SMEU
480	400	JLP4640SMEU
	420	JLP4642SMEU
	440	JLP4644SMEU
	460	JLP4646SMEU
500	480	JLP4648SMEU
	500	JLP4650SMEU
	400	JLP4840SMEU
	420	JLP4842SMEU
500	440	JLP4844SMEU
	460	JLP4846SMEU
	480	JLP4848SMEU
	500	JLP4850SMEU
500	460	JLP5046SMEU
	480	JLP5048SMEU
	500	JLP5050SMEU

Standard Insert -15 Day Lead Time

Seat Width	Seat Depth	JAY LitePro Part No.
340	360	JLP3436SMEU
420	380	JLP4238SMEU
500	420	JLP5042SMEU
	440	JLP5044SMEU





JAY[®] LitePro

JAY LitePro OF EU Template 21.01.2026

PELVIC LOADING AREA (PLA) INSERT* - STEP 3

	Select	Part #	INSERT
Standard Insert Option*	☐	S	S STEP 3

* The insert pad sizes supplied are based on PLA size and automatically determined when you select the cushion (see matrix)

* Includes as standard configuration a pair of Soft Insert Pods

OUTER COVER OPTION - STEP 4

	Select	Part #	OUTER COVER
Microclimatic	☐	M	M FIXED

JAY LitePro Smart Part Number

When ordering a JAY LitePro cushion, please follow the next steps

	Width	Depth	Insert	Outer Cover	
JLP	40	42	S	M	EU
fixed	STEP 1	STEP 2	STEP 3	fixed	fixed

Notes

The part # above (JLP4042SMEU) describes a JAY LitePro cushion with 400mm width and 420mm depth, Standard insert with two soft pods and microclimatic outer cover.

Additional explanations:

If any modifications of your cushion is needed please contact your Sunrise Medical authorized Customer Service.

JAY LitePro Costs	Unit price	Total price
Cushion Price		
Total Cost		



For more information
visit our webpage
www.sunrisemedical.com
© Sunrise Medical 2026



JAY® LitePro

JAY LitePro_OF_EU Template_21.01.2026

Foam Base

Charge

Cushion base (ww = width; dd = depth)

- ☐ BJLPwwddEU Base Assembly



€

Covers

Charge

Spare Outer Cover (ww = width; dd = depth)

- ☐ CJLPwwddMEU Microclimatic outer cover
☐ CJLPwwddIEU Incontinence Outer Cover

€

€

Inserts

Charge

PLA B

- ☐ JLPINSERTBEU Standard Multizone Elastomer Insert for PLA size B

€

PLA C

- ☐ JLPINSERTC Standard Multizone Elastomer Insert for PLA size C

€

PODS

- ☐ JLPSTAND Pressure Adjustable Pod - Soft
☐ JLPFIRM Pressure Adjustable Pod - Firm

€

€

Solid Seat

Charge

- ☐ JSIwwdd Solid Seat insert (wood plate)

€

JAY YOUR WAY Modifications

- ☐ NP Notch for Posts (Back Canes)
1" from the side of the cushion is removed. No modification is made to the cover.
- ☐ IW Increase Width
Available in incrementals of 1" or 2"
- ☐ IL Increase Length
Available in incrementals of 1" or 2"
- ☐ LLCL Left Leg Length
Adjustments of -1" or -2"
- ☐ LLCR Right Leg Length
Adjustments of -1" or -2"

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the regular item number construction



For more information
visit our webpage
www.sunrisemedical.com
© Sunrise Medical 2026