



SUNRISE
MEDICAL

Edition: 11.24
Valid until: 11.25

QUICKIE

NEON²

Deliver to:

Name:

Street:

Town:

Postal code:

E-Mail:

Tel:

Fax:



Invoice to:

Name:

Street:

Town:

Postal code:

E-Mail:

Tel:

Fax:

Order date:

Purchase order:

Number required:

☐ Order ☐ Calculation ☐ Built 4 me

Reimbursement code:

Art. no.: 05300001 IC

Marked for:

Sunrise order no:

- ☐ = Standard ☐ = Cost option
☐ = No cost option

Customer service fax number: +497253 980-222

Max. user weight: 140 kg ⁽¹⁾

USER WEIGHT

kg

Mandatory information to ensure correct configuration

Base Price

CLIENT MEASUREMENTS (all measurements without cushion)

SEAT WIDTH SELECTION - measured from outside frame tube to outside frame tube (mm)

300	320	340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Includes option NSW016602:

Userweight 140 kg, upcharge

SEAT DEPTH SELECTION - in mm

340	360	380	400	420	440	460	480	500
<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

(1) Max user weight 140kg only with option NSW010602 or with SW 440-500 mm;
Standard maximum user weight: 120kg.

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2025

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair
(approx. 10 - 20 mm altogether)
- b) On the wheelchair: From outside frame tube to outside frame tube.

Seat depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between
front edge of the seat upholstery and back of knee.
(approx. 30 - 60 mm). Care must be taken with adj. upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

- a) For the user: When sitting with 90° kneeangle:
From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front of seat tube to the ground.
(min. 30 mm longer than user measurement)
Remember to include the cushion height.

Rear seat height:

- b) On the wheelchair: From back of seat tube to the ground.
Note: Max. seat height difference rear to front is 110 mm

Back height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From back of seat tube to upper edge of backrest upholstery.
Note: The user's ability to sit upright has a major influence on the
back height and it can therefore differ considerably

Back angle:

Important: The back height depends on the degree of disability.
0° is the angle 90° to the ground. Negative angles are <90°
and positive angles are >90°.

Technical Data Overview:

Max. overall height:	1020 mm
Max. overall length:	1150 mm
Max. pivot width:	700 mm
Max. safe slope:	7° (for use of parking brakes)
Wheelchair weight:	from 11,5 kg
Dismantled weight:	from 9,1 kg
Mass of heaviest part:	2,1 kg
Max. overall width:	Camber
	For 0° camber: SW + 160 mm
	For 2° camber: SW + 220 mm
	For 4° camber: SW + 270 mm

Intended use and environment - for people with limited mobility; self-propelling or being pushed by an assistant, indoor and outdoor use.

Product is **crash tested according to ISO 7176-19** and therefore approved for transportation in a vehicle in combination with option "Tie down brackets for transportation". For transport without a user fix the chair safe in the vehicle (e.g. car or plane)

All changes in the configuration might change the seat heights.
Manufacturing tolerances of +/- 10 mm

Neon²

● = Standard ○ = No cost option □ = Cost option

Price in €

SEAT

Seat sling

NSW020002	●	Seat sling, nylon, tension adjustable, 1 utility bag	
NSW020003	□	Seat sling, nylon, tension adjustable, 2 utility bags (only from SW 380 mm)	
NSW020004	□	Seat sling, nylon, tension adjustable, 1 utility bag, 1 catheter bag (only from SW 340 mm)	
NSW020000	○	Without seat sling	

Seat cushion, cover with zipper, black

Foam (medium-soft)	Thickness:	NSW020101	□	30 mm	NSW020102	□	50 mm	
--------------------	------------	-----------	---	-------	-----------	---	-------	--

Seat accessories

NSW020050	□	Additional padding for seat upholstery, removable	
NSW090018	□	Lap belt with buckle, metal	

JAY Cushions - for complete product range see Jay order forms

JAY020002	□	JAY Basic	
JAY020003	□	JAY Soft Combi P	
JAY020001	□	JAY Easy Visco, Curved Base	
JAY020004	□	JAY Lite	
JAY020007	□	JAY Easy Fluid, Curved Base	

BACKREST

Back height range (mm)

Choose back height factory setting.

Height adjustments of in total 50 mm possible.

250 - 300	300 - 350	350 - 400	400 - 450	450 - 475
○ 250	○ 300	○ 350	○ 400	○ 450
○ 275	○ 325	○ 375	○ 425	○ 475
○ 300	○ 350	○ 400	○ 450	

NSW030010	○	Backrest, fixed (90°)	
NSW030031	□	Backrest, fold-down (to the front) (1)	
NSW030012	□	Backrest, angle-adjustable (75° to 103°)	
NSW030031	□	Backrest, fold-down (to the front) (1)	

Pushhandles

NSW030200	○	Without push handles	
NSW030201	○	Push handles long	
NSW030203	□	Push handles, fold-down (available from back height 325 mm)	
NSW030205	□	Push handles, height-adjustable for Jay backs (available from back height 300 mm)	
NSW030204	□	Push handles, height-adjustable (available from back height 300 mm)	

Backrest upholstery

NSW030316	○	EXO Backrest sling, padded, adjustable, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue	
NSW030317	□	EXO PRO Backrest, padded, adj., breathable spacer fabric, 1 accessory bag, black (EXO Evo version) Colour, binding tape: ● Black ○ Orange ○ Silver ○ Red ○ Blue	
NSW030300	○	Without backrest upholstery	


Accessories for backrest

NSW030402	□	Stabilising bar, automatic folding (only available from SW 360- 480 mm)	
NSW030401	□	Stabilising bar, swing-away	
NSW030412	□	Headrest, medium, height, depth and angle-adjustable	
NSW030411	□	Headrest, tall, height, depth and angle-adjustable, 200x110 mm	

(1)

Not with height adjustable push handles, JAY backs, stabilising bar or headrest. Please note folding-Matrix on next page. If push handles long chosen only from back height 300 mm available; if push handles fold-down chosen only from back height 350 mm available

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2025



NSW020002
Seat sling with 1 utility bag



NSW020003
Seat sling with 2 utility bags



NSW090018
Lap belt with buckle, metal



NSW030012
Backrest, angle adjustable



NSW030201
Push handles long



NSW030203
Push handles, fold-down



NSW030204
Push handles, height-adjustable



NSW030316
EXO backrest sling



NSW030317
EXO PRO backrest sling




NSW030412
Headrest, medium

Folding-Matrix for NSW030031 Backrest, fold-down to the front.



NSW030031
Backrest, fold-down



		Back height						
		325	350	375	400	425	450	475
Seat width	320	✓	✓	✓	✓	✓	✓	✓
	340	✓	✓	✓	✓	✓	✓	✓
	360	✓	✓	✓	✓	✓	✓	✓
	380		✓	✓	✓	✓	✓	✓
	400			✓	✓	✓	✓	✓
	420				✓	✓	✓	✓
	440				✓	✓	✓	✓
	460					✓	✓	✓
	480						✓	✓
500							✓	



=folds 90° to the front on folded chair



=can't be folded full 90° to the front on folded chair in combination with fold-down pushhandles



=can't be folded full 90° to the front on folded chair



=will fold less than 45° to the front on folded chair, folds 90° on unfolded chair

Neon²

● = Standard ○ = No cost option □ = Cost option

Price in €

JAY J3 BACKS

Jay J3 Back - Shallow Contour

Mid Thoracic:	<input type="checkbox"/>	JAY030010	Medium - 340 mm	<input type="checkbox"/>	JAY030033	Tall - 380 mm				
Upper Thoracic:	<input type="checkbox"/>	JAY030012	Short - 420 mm	<input type="checkbox"/>	JAY030019	Medium - 460 mm	<input type="checkbox"/>	JAY030025	Tall - 500 mm	
Shoulder Height:	<input type="checkbox"/>	JAY030015	Short - 530 mm	<input type="checkbox"/>	JAY030022	Medium - 570 mm	<input type="checkbox"/>	JAY030028	Tall - 610 mm	

Jay J3 Back - Mid Contour

Mid Thoracic:	<input type="checkbox"/>	JAY030011	Medium - 340 mm	<input type="checkbox"/>	JAY030018	Tall - 380 mm				
Upper Thoracic:	<input type="checkbox"/>	JAY030013	Short - 420 mm	<input type="checkbox"/>	JAY030020	Medium - 460 mm	<input type="checkbox"/>	JAY030026	Tall - 500 mm	
Shoulder Height:	<input type="checkbox"/>	JAY030016	Short - 530 mm	<input type="checkbox"/>	JAY030023	Medium - 570 mm	<input type="checkbox"/>	JAY030029	Tall - 610 mm	

JAY 3 Back - Mid Deep Contour (150 mm Contour Depth)

Mid Thoracic:							<input type="checkbox"/>	JAY030044	Medium - 340 mm	<input type="checkbox"/>	JAY030041	Tall - 380 mm				
Upper Thoracic:							<input type="checkbox"/>	JAY030045	Short - 420 mm	<input type="checkbox"/>	JAY030042	Medium - 460 mm	<input type="checkbox"/>	JAY030046	Tall - 500 mm	
Shoulder Height:							<input type="checkbox"/>	JAY030043	Short - 530 mm	<input type="checkbox"/>	JAY030047	Medium - 570 mm	<input type="checkbox"/>	JAY030048	Tall - 610 mm	

Jay J3 Carbon Back - Shallow Contour

Lower Thoracic:	<input type="checkbox"/> JAY060001	Short - 170 mm	<input type="checkbox"/> JAY060002	Tall - 240 mm	
Mid Thoracic:	<input type="checkbox"/> JAY060003	Short - 300 mm	<input type="checkbox"/> JAY060004	Tall - 380 mm	

Accessories for JAY J3 backs

WHT030415	<input type="checkbox"/> Whitmyer Headrest, COBRA XTRA, flip back	(only with JAY J3 backs)
WHT030416	<input type="checkbox"/> Whitmyer Headrest, AXYS, Standard	(only with JAY J3 backs)
JAY090058	<input type="checkbox"/> JAY Accessory bag	(only with Jay J3 Carbon backs)

ARMREST AND SIDEGUARDS

Armrest

NSW040005	<input type="checkbox"/> Desk sideguard, flip-up, removable	
NSW040001	<input type="checkbox"/> Desk sideguard, armpad (250 mm), fixed, flip-up, removable	
NSW040002	<input type="checkbox"/> Desk sideguard, armpad (250 mm), height-adjustable, flip-up, removable	
NSW040003	<input type="checkbox"/> Desk sideguard, armpad (310 mm), fixed, flip-up, removable	
NSW040004	<input type="checkbox"/> Desk sideguard, armpad (310 mm), height-adjustable, flip-up, removable	
NSW040050	<input type="checkbox"/> Single-post sideguard, armpad (250 mm), height-adjustable, removable (1)	
NSW040051	<input type="checkbox"/> Single-post sideguard, armpad (360 mm), height-adjustable, removable (1)	
NSW040052	<input type="checkbox"/> Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable (1)	
NSW040117	<input type="checkbox"/> Tubular armrest, Aluminium, swing away (sideways), padded, removable	(only with backrest, angle adj.)

Accessories and options for desk sideguards

NSW040020	<input type="checkbox"/> Release assistance for tetraplegics	(only with desk sideguard)
NSW090012	<input type="checkbox"/> Tray, swing-away (to the side), SW 340 - 500 mm, left-mount	(only with desk sideguard)
NSW090013	<input type="checkbox"/> Tray, swing-away (to the side), SW 340 - 500 mm, right-mount	(only with desk sideguard)
NSW090011	<input type="checkbox"/> Tray, push-fit, SW 380 - 480 mm	(only with desk sideguard)

Sideguards

NSW040107	<input type="checkbox"/> Sideguard, composite, removable, 210 mm high (310 mm long), black	
NSW040102	<input type="checkbox"/> Sideguard, aluminium, cold weather protection (in frame colour)	
NSW040103	<input type="checkbox"/> Sideguard, aluminium, fender, cold weather protection (in frame colour)	
NSW040104	<input type="checkbox"/> Sideguard, carbon fibre, fender, cold weather protection (black)	

HANGER

NSW050001	<input type="checkbox"/> Hanger angle 80°	
NSW050002	<input type="checkbox"/> Hanger angle 70° (2)	

Hanger, aluminium, swing-away

NSW050006	<input type="checkbox"/> Hanger, aluminium, swing-away, LH	
NSW050007	<input type="checkbox"/> Hanger, aluminium, swing-away, RH	
NSW050008	<input type="checkbox"/> Hanger, aluminium, swing-away, pair	

Hanger, aluminium, flip-up (0-90°), swing-away

NSW050009	<input type="checkbox"/> Hanger, aluminium, flip-up (0° - 90°), swing-away, LH	
NSW050010	<input type="checkbox"/> Hanger, aluminium, flip-up (0° - 90°), swing-away, RH	
NSW050011	<input type="checkbox"/> Hanger, aluminium, flip-up (0° - 90°), swing-away, pair	(from SW 380 mm)

- (1) Only with COG Adapter Standard
(2) In combination with NSW010024 "Frame, Hemi" only 70° hangers available

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2025

JAY J3 Back
Shallow Contour



JAY J3 Back
Mid Contour



JAY J3 Back
Mid Deep Contour



Jay 3 Carbon Back
Shallow Contour



NSW04000x
Desk sideguard, armpad, flip-
up, removable



NSW040050
Single-post sideguard



NSW040052
Single-post sideguard, armpad
(330 mm)



NSW040117
Tubular armrest, swing away,
removable, padded



NSW040107
Sideguard, composite,
removable (black)



NSW040103
Aluminium, cold weather
protection, fender



NSW040104
Carbon fibre, cold weather
protection, fender



NSW050011
Hanger, aluminum, flip-up (0-
90°, Swing-away

Neon²

● = Standard ○ = No cost option □ = Cost option

Price in €

FOOTBOARDS

NSW050104	□ Footboard, divided, composit, flip-up, heel loop, black (from SW 380 mm)	
NSW050020	□ Footboard, divided, composit, angle-adjustable, flip-up, heel loop, black (from SW 380 mm)	
NSW050022	□ Footboard, divided, aluminum, angle-adjustable, flip-up, heel loop, frame colour	
NSW050051	□ Footboard, divided, aluminum, angle-adjustable, flip-up, high-mount, heel loop, frame colour (1)	
NSW050023	□ Footboard, platform, aluminum, angle-adjustable, flip-up, heel loop, frame colour	
NSW050126	□ Swing-away mechanism for platform footboard	
NSW050059	□ Footboard, platform, Performance, aluminum, angle-adjustable, flip-up, calf strap (from SW 320 mm)	
NSW050131	□ Footboard, platform, lightweight, plastic, angle & depth-adjustable flip up, calf strap, black	
NSW050132	□ Footboard, platform, lightweight, carbon fibre, angle & depth-adjustable, flip up, calf strap, black (1)	

Accessories for footboards

NSW050127	□ Adjustable side protection, to create a side-stop on the footboard (2)	
-----------	--	--

LOWER LEG LENGTH (LLL)

(LLL has to be 25 mm less than front seat height; measurements in mm)

○ 220	○ 230	○ 240	○ 250	○ 260	○ 270	○ 280	○ 290	○ 300	○ 310
○ 320	○ 330	○ 340	○ 350	○ 360	○ 370	○ 380	○ 390	○ 400	○ 410
○ 420	○ 430	○ 440	○ 450	○ 460	○ 470	○ 480	○ 490	○ 500	○ 510

LLL below approximately 320 mm require high-mount footboards- depends on configuration, please check configurator for details

FRAME VERSIONS

NSW010015	○ Frame SA	
NSW010602	□ Frame reinforcement for max. user weight 140 kg	
NSW010024	○ Frame, Hemi for lower Seat heights (compare with your choice on the first page)	

CENTRE OF GRAVITY (COG)

(20 mm= very passive; 120 mm= very active)

NSW070703	Adapter standard	○ 20	○ 30	○ 40	○ 50	○ 60	○ 70	○ 80	○ 90	
		○ 100	○ 110	○ 120						
NSW070701	Adapter lightweight	○ 50	○ 60	○ 70	○ 80	○ 90	○ 100	○ 110	○ 120	

FRAME ACCESSORIES

(choose either tip assist or anti- tip tube per side)

NSW090016	□ Frame protection, black	
NSW090100	□ Tie down brackets for transportation	
NSW090010	□ Transit wheels (not in combination with anti-tip tubes or tip assist)	

Anti-tip tubes, swing-away

NSW090004	□ Left NSW090005	□ Right	
NSW090006	□ Pair	(no additional tip assist possible)	

Tip assist

NSW090001	□ Left NSW090002	□ Right	
NSW090003	□ Tip assist, pair	(no additional anti-tip tube possible)	
NSW090009	□ Crutch holder, loop	(only in combination with tip assist)	

EMPULSE F55 BUNDLING

* Comming soon

EMO000101	○ With F55 order, not mounted to the wheelchair (Please use separate OF to order the F55.)	
EMO000102	○ With F55 order, mounted to the wheelchair (Please use separate OF to order the F55.)	

FORKS

(compare with your choice on the first page)

NSW080013	● Castor fork, Carbotecture®, 111mm Colour: Black	
	Colour: □ Silver □ Orange □ Red (3) □ Blue (3)	
NSW080029	○ Castor fork, Aluminium Colour: Black (only for 123mm, 138mm and 174mm forks)	
	Colour: □ Silver □ Orange	
NSW080008	□ Castor fork, aluminium, one-arm Colour: Black	
	Colour: □ Silver □ Orange □ Red □ Blue	

STANDARD CASTORS

(compare with your choice on the first page)

NSW080301	○ Castor, solid (hard), rim: composite, black	
NSW080302	○ Castor, solid (soft), rim: composite, black (only 5" & 6")	
NSW080303	□ Castor, solid (hard), light-up, rim: composite, black (only in 4" and 5")	
NSW080304	□ Castor, solid (hard), rim: aluminium Colour, rim: Silver (only 4" & 5")	
NSW080305	○ Castor, pneumatic, rim: composite, black (only 7")	
NSW080306	□ Castor, solid (soft), rim: aluminium Colour, rim: Silver (only 4" & 5")	
	Colour, rim: □ Black □ Orange □ Red □ Blue	

- (1) From SW 340 mm available
 (2) Only for footboard, platform lightweight or Performance
 (3) Only available in combination with castor, soft (NSW080302 or NSW080306)

For more information
 visit our webpage
www.SunriseMedical.eu
 © Sunrise Medical 2025



NSW050022
Footboard divided,
aluminium, frame colour



NSW050059
Footboard platform,
Performance



NSW050132
Footboard platform,
lightweight, carbon fibre



NSW050127
Adjustable side protection



NSW090010
Transit wheels



NSW090006
Anti-tip tubes, swing away



NSW090001
Tip assist



NSW080013 & NSW080306
Castor fork, Carbotecture & Castor solid (soft) aluminum



NSW080008
Castor fork, aluminum, one-arm, red

F55 BUNDLING - NOGO MATRIX FOR 24" BACKWHEELS

COG	Seat Depth						
	380	400	420	440	460	480	500
40	X	X	X	✓	✓	✓	✓
60	X	X	X	X	✓	✓	✓
70	X	X	X	X	✓	✓	✓
80	X	X	X	X	X	✓	✓
90	X	X	X	X	X	✓	✓
100	X	X	X	X	X	X	✓
110	X	X	X	X	X	X	✓
120	X	X	X	X	X	X	X
130	X	X	X	X	X	X	X

Neon²

● = Standard ○ = No cost option □ = Cost option

Price in €

REAR WHEELS

(compare with your choice on the first page)

NSW070000	○	Without rear wheel	-
NSW070001	○	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes (not in 25")	-
NSW070002	○	Rear wheel, Design, 36 spokes (silver), straight spokes (only 24" and 26")	-
NSW070003	○	Rear wheel, Lightweight, 24 spokes (silver), straight spokes (not in 26"; see colour options on last page)	-
NSW070004	□	Rear wheel, Proton, 24 spokes (black), straight spokes (not in 26"; user weight limit 100kg)	376 €
NSW070005	□	Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver (not in 22")	678 €
NSW070007	□	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub black (only 24" & 25")	951 €
NSW070008	□	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes (1) (only 22" & 24")	379 €
NSW070009	□	Rear wheel, Mountainbike, 36 spokes (silver), straight spokes (only 24")	344 €
NSW070010	□	Mountainbike wheel, replacement kit incl. tyre (extra coarse profile) (only 24")	522 €

Axle

NSW070150	●	Quick-release axle	
NSW070151	□	Quick-release axle, designed for quads (2)	

TYRES

NSW070101	○	Tyre, pneumatic, fine profile, Schwalbe Right Run (not in 25")	
NSW070102	○	Tyre, puncture proof, tread (not in 26")	
NSW070107	□	Tyre, pneumatic, slick, Schwalbe One (only 24" and 25")	
NSW070104	□	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution	
NSW090163	□	Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube (not with solid tyres)	

HANDRIMS

NSW070204	○	Handrims, stainless steel, silver	
NSW070202	○	Handrims, aluminium, anodised, silver	
NSW070201	○	Handrims, aluminium, powder-coated, clear	
NSW070203	□	Handrims, aluminium, powder-coated, coloured (select colour on last page)	
NSW070205	□	Handrims, titanium, untreated (not in 22")	
NSW070206	□	Handrims, titanium, powder-coated, clear (not in 22")	
NSW070207	□	Handrims, Supergrip ®	
NSW070208	□	Handrims, Maxgripp ® (only 22" & 24")	
NSW070212	□	Handrims, Ellipse 3R, black (ergonomic handrims with rubber-strip, 24" and 25" only)	

Hand rim mounting

NSW070307	□	Handrim mount, super narrow (7 mm) (only with stainless steel)	
NSW070317	●	Handrim mount, narrow (17 mm)	
NSW070327	○	Handrim mount, wide (27 mm)	

Handrim accessories

NSW070250	□	Handrim cover, rubber (not with 25" and 26")	
-----------	---	---	--

WHEEL CAMBER

NSW070600	○	0°	
NSW070602	○	2° (3)	
NSW070604	○	4° (3)	

SPOKE GUARDS

(coloured: only 22 & 24"; transparent: only 22", 24" & 26")

STONE COLOURS						
NSW070831	<input type="checkbox"/>	Transparent	NSW070813	<input type="checkbox"/>	Glaktica	
NSW070823	<input type="checkbox"/>	Princess	NSW070818	<input type="checkbox"/>	Jungle Fever	
NSW070816	<input type="checkbox"/>	Happy Lions	NSW070811	<input type="checkbox"/>	Feather	
NSW070819	<input type="checkbox"/>	Monkeys	NSW070821	<input type="checkbox"/>	Pirates	
NSW070820	<input type="checkbox"/>	Oriental	NSW070811	<input type="checkbox"/>	Feather	
NSW070829	<input type="checkbox"/>	Superhero	NSW070814	<input type="checkbox"/>	Graffiti	
NSW070817	<input type="checkbox"/>	Hot Wheels	NSW070802	<input type="checkbox"/>	Basketball	NSW070803 <input type="checkbox"/> Basketball 2
NSW070822	<input type="checkbox"/>	Pool	NSW070812	<input type="checkbox"/>	Flying Water	
NSW070827	<input type="checkbox"/>	Silk	NSW070828	<input type="checkbox"/>	Space	
NSW070830	<input type="checkbox"/>	Tiger	NSW070833	<input type="checkbox"/>	Zebra	

- (1) Not with push handles fold down & not without push handles
only with adapter, standard and only with 0° and 4° wheel camber
- (2) Not with rear wheels Spinergy / mountainbike/ drum brakes
- (3) Not with COG Adapter, lightweight

For more information
visit our webpage
www.SunriseMedical.eu
© Sunrise Medical 2025



NSW070001
Universal wheel



NSW070002
Design wheel



NSW070003
Lightweight wheel



NSW070004
Proton wheel



NSW070010
Mountainbike wheel



NSW070101
Schwalbe Right Run



NSW070107
Schwalbe One



NSW070104
Schwalbe Marathon Plus
Evolution



NSW070102
Tyre, puncture proof



NSW070208
Handrims Maxgrepp



NSW070212
Handrims Ellipse 3R

Neon²

● = Standard ○ = No cost option □ = Cost option

Price in €

WHEEL LOCKS

(choose only one option)

NSW060001	● Standard brake		
NSW060002	□ Knee-lever brake		
NSW060010	□ Brake lever extension, knee-lever brake, fold-down		
NSW060015	□ Push to lock brake, one-arm, left actuation		
NSW060016	□ Push to lock brake, one-arm, right actuation		
NSW060003	□ Compact brake		

ACCESSORIES

NSW090000	● Tool kit		
NSW090025	□ Air pump		
NSW090024	□ High pressure pump (0 - 11 bar)		
NSW090030	□ Back pack		
NSW090026	□ Mobile phone pocket		

COLOURS

(Please consult our colour catalogue to review your colour choice)

Frame	Handrim	Hub/rim (1)	Frame	Handrim	Hub/rim (1)
02 Brilliant Silver	○	□	27 Graphitemetallic	○	□
03 Silver Matt	□	□	29 Black	○	□
06 Sun Yellow	○	□	31 Matt Black	○	□
14 Blue	○	□	32 Signal White	○	□
15 Brilliant Marine	○	□	33 Pink	○	□
16 Matt Marine	□	□	34 Dark Chocolate Metallic	□	□
19 Brilliant Purple	○	□	35 Skyblue	□	□
21 Red Ruby	○	□	36 Cyan-Metallic	□	□
22 Brilliant Red	□	□	37 Marsred (Orange)	□	□
23 Matt Red	□	□	40 Kiwi	□	□
24 Tutti Frutti Orange	○	□	41 Flamingo-Glitter	□	□
			42 Matt-Anthracite	○	□
			46 Mamba Green (new)	□	□

(1) only with rear wheel, Universal 22" and 24"

COLOURS FOR REAR WHEEL, LIGHTWEIGHT, 24"

(hub and rim anodized, black spokes)

NSW100336	□ Anodized, orange	(only with Rear wheel, lightweight in 24")	
NSW100337	□ Anodized, blue	(only with Rear wheel, lightweight in 24")	
NSW100338	□ Anodized, red	(only with Rear wheel, lightweight in 24")	

COLOURS - for stickers ("Neon²" on the frame and "Quickie" on sideguard)

NSW100550	○ black		
NSW100551	○ silver		
NSW100552	○ orange		
NSW100553	○ red		
NSW100554	○ blue		

Comments for special build - please apply for a quote; lead time conditions my change.



Quickie Neon² - Basic Price

Option price €
Packaging and transport €
Exclusive VAT €
VAT. €
Total Price €

Sunrise Medical Customer service

Address: Kahlbachring 2-4, D-69254 Malsch/Heidelberg

Tel. +497253 980-0

Fax. +497253 980-222

E-mail: export@sunrisemedical.de

Website: sunrisemedical.eu

For more information

visit our webpage

www.SunriseMedical.eu

© Sunrise Medical 2025



NSW060001
Standard brake



NSW060002
Knee lever brake



NSW060003
Compact brake



NSW090000
Tool kit



NSW090024
High pressure pump



NSW090025
Air pump



NSW090030
Back pack



NSW090026
Mobile phone pocket



NSW100336
Lightweight wheel, anodized orange



NSW100337
Lightweight wheel, anodized blue



NSW100338
Lightweight wheel, anodized red