

QS7 X

SA



ORDER FORM

Errors and omissions excepted.

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address
C.-ID:			C.-ID:		
Name:			Name:		
Street:			Street:		
Town/PC:			Town/PC:		
E-Mail:			E-Mail:		
Phone:			Phone:		
Order Information			User Weight		
Order	Calculation	Built-4-Me®	Max. user weight Art. Nr. 773000006 / 773000007:		125 kg
Purchase Order No.:			Art. No.: 773000006 DE / 773000007 IC		
Reimbursement code: XXX			User Weight:		kg
Order Date:	Day	Month	Year	Base Price QS7 X SA € X.XXX,-	

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether).
- b) On the wheelchair: From outside frame tube to outside frame tube.

Seat depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee (approx. 30 - 60 mm).
Care must be taken with adjustable upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

- a) For the user: When sitting with 90° kneeangle: From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front edge of the seat-upholstry, next to the seat-tube to the ground (min. 30 mm longer than user measurement).
Remember to include the cushion height.

Rear seat height:

On the wheelchair: From rear edge of the seat-upholstery, next to the seat-tube to the ground.

Note: Maximum seat height difference rear to front is 110 mm.

Back height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From top of the seat-upholstery to upper edge of backrest upholstery.

Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably.

Back angle:

Backrest-angle is measured in relation to the seat. 90°= perpendicular to the seat.

Frame			€
QSS010015	Frame SA (removable footrest)		XXX,-

Technical Data

Max. overall height:	1050 mm	Max. overall width with:	24" Rear wheel	25" Rear wheel
Max. overall length:	1050 mm	For 0° camber:	SW + 160 mm	SW + 160 mm
Max turning circle:	700 mm	For 2° camber:	SW + 210 mm	SW + 220 mm
Overall weight:	10,4 kg	For 4° camber:	SW + 260 mm	SW + 270 mm
Weight of heaviest detachable part:	2,1 kg			
Transport weight (without rear wheels):		from 7,8 kg		
Max. safe slope (with wheels locked):		10°		

Please note that options such as i.e. armrest-sideguards or drum brakes can increase the overall width of the wheelchair. If this is the case, under certain circumstances, it may not be possible to use some or all of the available escape routes from a building. Please note that in certain configurations the wheelchair may exceed a width of 700 mm. It may be more difficult or impossible to travel on public transport.

Product is **crash tested according to ISO 7176-19** and therefore approved for transportation in a vehicle. In combination with QSS090100 (Tie down brackets for transportation) For transport without a user the chair must be safely clamped in the vehicle.

Please note that the technical data is based on a std configuration:

Standard wheelchair 400 x 400 mm; 6 COG; Seat angle of 8°

Manufacturing tolerances can be +/- 10 mm on all configurations

Lightweight-Guide for ordering the QS7 X

Sunrise Medical offers a wide variety of options for the QS7 X to suit every user's needs. However, most options add significant weight to the wheelchair. To make this more transparent when filling in the order-form, and to assist you in ordering the lightest possible wheelchair, we indicate the approximate weight-increase of the options on the order form:

(L) Indicates the lightest possible option within the range of options

(+) Indicates that this option adds more than 0,5 kg to the wheelchair's weight compared to the lightest option within the range

Client Measurements all measurements without cushion**Seat Width** measured from outside frame tube to outside frame tube

(mm)	320	340	360	380	400	420	440	460
	480	500						

Seat Depth

(mm)	340	360	380	400	420	440	460	480	500
	340/ 360/ 380 mm = 380 frame length			each frame in a separate length (20 mm increments)					

Front Seat Height Selection in dependence of front castors and frame versions

Fork Type (mm)	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570
3"															
Aluminium (90 mm)			–	–	–	–	–	–	–	–	–	–	–	–	–
Carbotecture (105 mm)	–	–			–	–	–	–	–	–	–	–	–	–	–
Aluminium (115 mm)	–	–	–			–	–	–	–	–	–	–	–	–	–
Aluminium (131 mm)	–	–	–	–	–			–	–	–	–	–	–	–	–
Aluminium (166 mm)	–	–	–	–	–	–	–	–			–	–	–	–	–
Aluminium One Arm (105 mm)	–	–			–	–	–	–	–	–	–	–	–	–	–
Aluminium Frog Leg	–	–	–	–	–			–	–	–	–	–	–	–	–
4"															
Aluminium (90 mm)	–	–			–	–	–	–	–	–	–	–	–	–	–
Carbotecture (105 mm)	–	–				–	–	–	–	–	–	–	–	–	–
Aluminium (115 mm)	–	–	–			–		–	–	–	–	–	–	–	–
Aluminium (131 mm)	–	–	–	–	–				–	–	–	–	–	–	–
Aluminium (166 mm)	–	–	–	–	–	–	–	–			–		–	–	–
Aluminium One Arm (105 mm)	–	–				–	–	–	–	–	–	–	–	–	–
Aluminium Frog Leg	–	–	–	–	–				–	–	–	–	–	–	–
5"															
Aluminium (90 mm)	–	–	–	–		–	–	–	–	–	–	–	–	–	–
Carbotecture (105 mm)	–	–	–	–		–		–	–	–	–	–	–	–	–
Aluminium (115 mm)	–	–	–	–	–				–	–	–	–	–	–	–
Aluminium (131 mm)	–	–	–	–	–					–	–	–	–	–	–
Aluminium (166 mm)	–	–	–	–	–	–	–	–		–				–	–
Aluminium One Arm (105 mm)	–	–	–	–		–		–	–	–	–	–	–	–	–
Aluminium Frog Leg	–	–	–	–	–	–	–		–		–	–	–	–	–
6"															
Aluminium (115 mm)	–	–	–	–	–	–	–	–			–	–	–	–	–
Aluminium (131 mm)	–	–	–	–	–	–	–		–			–	–	–	–
Aluminium (166 mm)	–	–	–	–	–	–	–	–	–	–					–
Aluminium (Frog Leg)	–	–	–	–	–	–	–	–	–	–		–	–	–	–
7"															
Aluminium (131 mm)	–	–	–	–	–	–	–	–	–	–			–	–	–
Aluminium (166 mm)	–	–	–	–	–	–	–	–	–	–	–				

Carbotecture (105 mm) not approved for use with powered add-ons (powered wheels or push-devices). In general, always verify whether a combination agreement is in place before fitting powered add-ons.

Rear Seat Height Selection in dependence of rear wheel size

Rear seat height in mm	370	380	390	400	410	420	430	440	450	460	470	480	490	500
24"													1	–
25"	–													1

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

(1) Only with lightweight axle adapter

Seat Upholstery and Cushions

€

Seat sling

QSS020153	Seat sling, nylon, tension adjustable	–
QSS020154	1 utility bag	–
QSS020155	2 utility bags	XXX,-
QSS020000	Without seat sling	–

Seat cushion, cover with zipper, black

Foam (medium-soft) Thickness:

QSS020101	30 mm	XXX,-
QSS020102	50 mm	XXX,-

JAY Cushions – for complete product range see JAY order forms

JAY020002	JAY Basic	XXX,-
JAY020003	JAY Soft Combi P	XXX,-
JAY020001	JAY Easy Visco C	XXX,-
JAY020007	JAY Easy Fluid	XXX,-
JAY020067	JAY LitePro	XXX,-

Seat accessories

QSS090018	Lap belt with buckle, metal	XXX,-
-----------	-----------------------------	-------

Backrest

€

	250 - 300	300 - 350	350 - 400	400 - 450	450 - 475
Backrest height range Choose back height factory setting (mm). Height adjustments of in total 50 mm possible.	250	300	350	400	450
	275	325	375	425	475
	300	350	400	450	

Push handles

QSS030200	Without push handles (L)	–
QSS030201	Push handles, long	–
QSS030203	Push handles, fold-down Only available from back height 325 mm upwards	XXX,-
QSS030204	Push handles, height-adjustable (+) Only available from back height 300 mm upwards	XXX,-
QSS030205	Push handles, height-adjustable for JAY backs (+) Only available from back height 325 mm upwards	XXX,-

Back Design

QSS030010	Backrest, fixed (90°)	–
QSS030029	Desk sideguard preparation for later replacement necessary for retrofitting Desk sideguard or Tubular armrest with Backrest, fixed 90°	–
QSS030031	Backrest, fold-down (to the front) not with height adj. push handles; not with JAY backs; only from back height 350 mm; fixed non-adjustable back height, only from backrest height 400 mm upwards in combination with push handles fold-down; please see Folding Matrix on next side; not with stabilising bar options	XXX,-
QSS030012	Backrest, angle-adjustable 75° to 103°	XXX,-
QSS030029	Desk sideguard preparation for later replacement necessary for retrofitting Desk sideguard	–
QSS030031	Backrest, fold-down (to the front) not with height adj. push handles; not with JAY backs; only from back height 350 mm; fixed non-adjustable back height, only from backrest height 400 mm upwards in combination with push handles fold-down; please see Folding Matrix on next side; not with stabilising bar options	XXX,-

Backrest upholstery (black)

QSS030316	EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black	–
	Colour binding tape: Black Orange Silver Red Blue	–
QSS030317	EXO PRO Backrest sling, padded, adjustable, breathable, spacer fabric, 1 accessory bag, black	XXX,-
	Colour binding tape: Black Orange Silver Red Blue	–
QSS030300	Without backrest upholstery	–

Seat Upholsteries & Cushions




// QSS020101-QSS020102
Seat cushion, cover with zipper, black



// QSS090018
Lap belt with buckle, metal



// JAY020002
JAY Basic (pictured without cover)



// JAY020003
JAY Soft Combi P (pictured without cover)



// JAY020001
JAY Easy Visco (pictured without cover)



// JAY020007
JAY Easy Fluid (pictured without cover)



// JAY020067
JAY Lite Pro (pictured without cover)

Backrest



// QSF030201
Push handles, long



// QSF030203
Push handles, fold-down



// QSF030204
Push handles, height-adjustable

Overview of Foldability – 030031 Backrest Foldable Forward

Backheight	350	375	400	425	450	475
Seatwidth						
320						
340						
360						
380						
400						
420						
440						
460						
480						
mm	500					



- The backrest can also be folded 90° forward when the wheelchair is folded.
- In combination with fold-down push handles, the backrest cannot be folded a full 90° forward when the wheelchair is folded.
- When the wheelchair is folded, the backrest cannot be folded a full 90° forward.
- When the wheelchair is folded, the backrest can be folded less than 45° forward. When the wheelchair is unfolded, the backrest can be folded 90° forward.

Backrest accessories

QSS030401	Stabilising bar, swing-away only with height adjustable pushhandles and from SW 380 mm	XXX,-
QSS030411	Headrest, height-, depth- and angle-adjustable medium (200 x 110 mm)	XXX,-
QSS030410	Headrest, height-, depth- and angle-adjustable tall (230 x 110 mm)	XXX,-
QSS030402	Stabilising bar, automatic folding only available from SW 380mm; not with fold-down backrest	XXX,-

JAY BACKS plus ACCESSORIES max. user weight 136 kg

€

JAY J3 Back, Shallow Contour		Short (S)		Medium (M)		Tall (T)		
Mid Thoracic				JAY030010	340 mm	JAY030033	380 mm	XXX,-
Upper Thoracic		JAY030012	420 mm	JAY030019	460 mm	JAY030025	500 mm	XXX,-
JAY J3 Back, Mid Contour		Short (S)		Medium (M)		Tall (T)		
Mid Thoracic				JAY030011	340 mm	JAY030018	380 mm	XXX,-
Upper Thoracic		JAY030013	420 mm	JAY030020	460 mm	JAY030026	500 mm	XXX,-
JAY J3 Back - Mid Deep Contour		Short (S)		Medium (M)		Tall (T)		
Mid Thoracic				JAY030044	340 mm	JAY030041	380 mm	XXX,-
Upper Thoracic		JAY030045	420 mm	JAY030042	460 mm	JAY030046	500 mm	XXX,-
JAY J3 Carbon Back - Shallow Contour		Short (S)		Medium (M)		Tall (T)		
Low Thoracic		JAY060001	170 mm			JAY060002	240 mm	XXX,-
Mid Thoracic		JAY060003	300 mm			JAY060004	380 mm	XXX,-

Accessories for JAY J3 backs

WHT030415	Whitmyer Headrest, COBRA XTRA, flip back only with JAY J3 backs; not with J3 carbon	XXX,-
WHT030416	Whitmyer Headrest, AXYS, Standard only with JAY J3 backs; not with J3 carbon	XXX,-
JAY090058	JAY Accessory bag only with JAY J3 Carbon backs	XXX,-

Armrests & Sideguards

€

Sideguard

QSS040101	Sideguard, carbon fibre, cold weather protection, slim (black) (L)	XXX,-
QSS040102	Sideguard, aluminium , cold weather protection (in frame colour)	XXX,-
QSS040103	Sideguard, aluminium , fender, cold weather protection (in frame colour)	XXX,-
QSS040104	Sideguard, carbon fibre (EVO), fender, cold weather protection (black)	XXX,-
QSS040107	Sideguard, composite, removable, 210 mm high (310 mm long), black	XXX,-
QSS040120	Sideguard, composite, fender (black)	XXX,-

Armrest

QSS040005	Desk sideguard, flip-up, removable	XXX,-
QSS040001	Desk sideguard, armpad (250 mm), fixed, flip-up, removable	XXX,-
QSS040002	Desk sideguard, armpad (250 mm), height-adjustable, flip-up, removable	XXX,-
QSS040003	Desk sideguard, armpad (310 mm), fixed, flip-up, removable	XXX,-
QSS040004	Desk sideguard, armpad (310 mm), height-adjustable, flip-up, removable	XXX,-
QSS040052	Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable (+)	XXX,-
QSS040117	Tubular armrest. Alu, swing away, padded, removable (+)	XXX,-

Accessories armrest

QSS090011	Tray, push-fit, SW 380 - 480 mm	XXX,-
-----------	---------------------------------	-------

JAY BACKS plus ACCESSORIES



// JAY030012
JAY J3 Back, Shallow Contour
Short (S) Upper Thoracic



// JAY030013
JAY J3 Back, Mid Contour Short (S)
Upper Thoracic



// JAY030045
JAY J3 Back - Mid Deep Contour
Short (S) Upper Thoracic

Armrests & Sideguards



// QSS040101
Sideguard, carbon fibre, cold wea-
ther protection, slim (black)



// QSS040102
Sideguard, aluminium, cold wea-
ther protection (in frame colour)



// QSS040103
Sideguard, aluminium, fender, cold
weather protection (in frame colour)



// QSS040104
Sideguard, carbon fibre EVO, fender,
cold weather protection (black)



// QSS040107
Sideguard, composite, removable,
210 mm high (310 mm long), black



// QSS040120
Sideguard, composite, fender
(black)



// QSS040005
Desk sideguard, flip-up, removable



// QSS040001
Desk sideguard, armpad (250 mm),
fixed, flip-up, removable



// QSS040002
Desk sideguard, armpad (250 mm),
height-adjustable, flip-up, removab-
le

Hanger Angle Angle measured between seat and hanger tube		
QSS050001	Hanger angle 80° (EU: 100°)	–
QSS050002	Hanger angle 70° (EU: 110°)	–
QSS050041	Hanger angle 90°	–
Hanger		
QSS050006	Hanger, aluminium, swing-away, LH	XXX,-
QSS050007	Hanger, aluminium, swing-away, RH	XXX,-
QSS050008	Hanger, aluminium, swing-away, pair	XXX,-
	Hanger, aluminium, flip-up (0° - 90°), swing-away	
	QSS050009 left	XXX,-
	QSS050010 right	XXX,-
	QSS050011 pair	XXX,-

Footboards			€
QSS050020	Footboard, divided, composite, angle-adjustable, flip-up (sideways), heel loop, black from SW 380 mm		–
QSS050050	Footboard, divided, aluminium, angle-adjustable, flip-up (sideways), heel loop, black		XXX,-
QSS050102	Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), heel loop, black		XXX,-
	QSS050126	Swing-away mechanism for platform footboard minimum LLL 30 mm longer	XXX,-
	QSS050062	Locking system, platform footplate, aluminium	XXX,-
QSS050059	Footboard, platform, Performance, aluminium, angle-adjustable, flip-up (sideways), calf strap		XXX,-
QSS050131	Footboard, platform, lightweight, plastic, angle & depth-adj. flip up (sideways), calf strap, black		XXX,-
QSS050132	Footboard, platform, lightweight, carbon fibre, angle & depth-adj. flip up (sideways), calf strap, black (L)		XXX,-
Accessories for footboards			
QSS050127	Adjustable side protection, to create a side-stop on the footboard only for footboard, platform lightweight or Performance		XXX,-

Lower Leg Length LLL has to be 25 mm less than front seat height; all measurements are in mm										
220	230	240	250	260	270	280	290	300	310	320
330	340	350	360	370	380	390	400	410	420	430
440	450	460	470	480	490	500				

Centre of Gravity (COG) measured from front of back post; choose only one option (0 = very passive 13 = very active)									€
The COG scale reflects the distance in cm from the backrest tube to the middle of the axle stem in a 6° configuration. The exact cm measurement will differ slightly depending on the seat height configuration of the chair.									
QSS070703	Axle plate, standard	1 ¹	2	4	5	6	8		–
		9	10						
	¹ comes installed with passive adapter COG 9 only from SD 40 COG 10 only from SD 42								
QSS070701	Axle plate, lightweight (L) ²	4 ¹	5	7	8	9	11		–
	¹ comes installed with passive adapter ² only with 0° wheel camber and not compatible with R90 folding axle adapters COG 9 only from SD 40 COG 11 only from SD 42								
QSS010411	Axle plate, aluminium, long (full adjustment)								–
QSS010412	Axle plate, aluminium, shorted (max. adjustment rear seat height +10mm) (L)								XXX,-

Armrests & Sideguards



// QSS040052
Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable



// QSS040117
Tubular armrest. Alu, swing away, padded, removable



// QSS050001
Hanger angle 80° (100°),
QSS050006-8 Hanger aluminium,
Swing-Away



// QSS050002
Hanger angle 70° (110°)
QSS050006-8 Hanger aluminium,
Swing-Away



// QSS050041
Hanger angle 90°
QSS050006-8 Hanger aluminium,
Swing-Away



// QSS050009 - QSS050010
Hanger, aluminium, flip-up (0° -
90°), swing-away left & right

Footboards



// QSS050020
Footboard, divided, composite, angle-adjustable, flip-up (sideways), heel loop, black



// QSS050050
Footboard, divided, aluminium, angle-adjustable, flip-up (sideways), heel loop, black



// QSS050102
Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), heel loop, black



// QSS050126
Swing-away mechanism for platform footboard



// QSS050062
Locking system, platform footplate, aluminium



// QSS050059
Footboard, platform, Performance, aluminium, angle-adjustable, flip-up (sideways), calf strap



// QSS050131
Footboard, platform, lightweight, plastic, angle & depth-adjust. flip up (sideways), calf strap, black



// QSS050132
Footboard, platform, lightweight, carbon fibre, angle & depth-adjust. flip up (sideways), calf strap, black



// QSS050127
Adjustable side protection, to create a side-stop on the footboard

Centre of Gravity (CoG)



// QSS070703
Axle plate Adapter, standard



// QSS070701
Axle plate Adapter, lightweight

Frame Accessories (choose either tip assist or anti- tip tubes per side) €

QSS090100	Tie down brackets for transportation according to ISO 7176-19 (max. user-weight 136 kg)	XXX,-
QSS110156	Twin-LED light Mounted on both sides of frame	XXX,-
QSS010069	Transport Protection, black Not compatible with desk sideguards (including preparation for desk sideguards), tubular armrest, anti-tips, tip assist, transit wheels, or tie down brackets	XXX,-

Tip assist

QSS090001	Left	QSS090002	Right	XXX,-
QSS090003	Pair			XXX,-
QSS090009	Crutch holder, loop only in combination with tip assist			XXX,-
QSS090050	Anti-tip tubes, plug-in, pair			XXX,-

Anti-tip tubes, swing-away (+)

QSS090004	Left	QSS090005	Right	XXX,-
QSS090006	Pair			XXX,-
QSS090010	Transit wheels			XXX,-

Empulse F55 Bundling €

EM0000101	With F55 order, not mounted to the wheelchair Please use separate OF to order the F55. Only with compact brake options (QSS060023 & QSS060031)	XXX,-
EM0000102	With F55 order, mounted to the wheelchair Please use separate OF to order the F55. Only with Lightweight Compact brake	XXX,-

Empulse F35 Bundling €

EM0000301	With F35 order, not mounted to the wheelchair Please use separate OF to order the F35. Only with Lightweight Compact brake	XXX,-
EM0000302	With F35 order, mounted to the wheelchair Please use separate OF to order the F35. Only with Lightweight Compact brake	XXX,-

Empulse R90 Bundling €

EM0000401	With R90 order, not mounted to the wheelchair Please use separate OF to order the R90, R90 must only be used in combination with anti-tips; R90 only with QSS070703	XXX,-
EM0000402	With R90 order, mounted to the wheelchair Please use separate OF to order the R90, R90 must only be used in combination with anti-tips; R90 only with QSS070703	XXX,-

Empulse M90 Bundling €

EM0000501	With M90 order, not mounted to the wheelchair Please use separate OF to order the M90, M90 must only be used in combination with anti-tips; M90 only with 0° wheel camber	XXX,-
EM0000502	With M90 order, mounted to the wheelchair Please use separate OF to order the M90, M90 must only be used in combination with anti-tips; M90 only with 0° wheel camber	XXX,-

Forks compare with your choice on the first page €

QSS080013	Castor fork, Carbotechnology ®, (105 mm) Colour: Black (L) not compatible with power add on devices	–
	Orange Silver Red Blue Red & Blue only in combination with castor, soft (QSS080302 and QSS080306)	XXX,-
QSS080029	Castor fork, Aluminium, Black (anodized) only 90 mm, 115 mm, 131 mm or 166 mm fork; 166 mm not compatible with power add on devices	–
	Colour Cast. Fork Silver (anodized)	XXX,-
QSS080064	Castor fork, Aluminium, Powder coated only 90 mm, 115 mm, 131 mm or 166 mm fork; 166 mm not compatible with power add on devices	–
	Colour Cast. Fork Starlight White (04)	XXX,-
	Colour Cast. Fork Royal Blue (12)	XXX,-
	Colour Cast. Fork Red Pepper (15)	XXX,-
	Colour Cast. Fork Solar Orange (21)	XXX,-
QSS080008	Castor fork, aluminium, one-arm (105 mm) Colour: Black	XXX,-
	Orange Silver Red Blue	XXX,-
QSS080007	Castor fork, aluminium, suspension (Frog Legs ®), black (+)	XXX,-

Forks



// QSS090100
Tie down brackets for transportation according to ISO 7176-19 (max. user-weight 136 kg)



// QSS090009
Crutch holder, loop



// QSS090001
Tip assist



// QSS090010
Transit wheels



// QSS090004
Anti-tip tubes, swing-away



// QSS090050
Anti-tip tubes, plug-in, pair



// QSS110156
Twin-LED light

Bundling - QS7 X SA / 24" & 25" Wheels (orange is no go with 25" wheels) / F55 & F35

COG	Seath Depth						
	380	400	420	440	460	480	500
1	✓ 24", x 25"	✓	✓	✓	✓	✓	✓
2	x	✓ 24", x 25"	✓	✓	✓	✓	✓
4	x	x	✓	✓	✓	✓	✓
5	x	x	✓ 24", x 25"	✓	✓	✓	✓
6	x	x	x	✓ 24", x 25"	✓	✓	✓
7	x	x	x	✓ 24", x 25"	✓	✓	✓
8	x	x	x	x	✓	✓	✓
9	x	x	x	x	✓ 24", x 25"	✓	✓
10	x	x	x	x	x	✓ 24", x 25"	✓
11	x	x	x	x	x	✓ 24", x 25"	✓

Forks



// QSS080013
Castor fork, Carbotecture®, (105 mm) Colour: Black



// QSS080064
Castor fork, Aluminium, Powder coated (115 mm)



// QSS080008
Castor fork, aluminium, one-arm (105 mm) Colour: Black

Castor compare with your choice on the first page						€
QSS080305	Castor, pneumatic, rim: composite, black only for 7" castors					–
QSS080301	Castor, solid (hard), rim: composite, black					–
QSS080302	Castor, solid (soft), rim: composite, black only for 5" & 6" castors					–
QSS080303	Castor, solid (hard), light-up, rim: composite, black only for 4" & 5" castors					–
QSS080304	Castor, solid (hard), rim: aluminium Colour, rim: Silver only for 4" & 5" castors					XXX,-
QSS080306	Castor, solid (soft), rim: aluminium Colour, rim: Silver only for 4" & 5" castors					XXX,-
	Colour rim:	Black	Orange	Red	Blue	XXX,-
QSS080060	Castor, solid (soft), rim: carbon-composite, ultralight (L) only for 4" castors					XXX,-

Rear Wheels compare with your choice on the first page		€
QSS070000	Without rear wheel	–
QSS070001	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes (+) only 24"	–
QSS070002	Rear wheel, Design, 36 spokes (silver), straight spokes only 24"	–
QSS070003	Rear wheel, Lightweight, 24 spokes (silver), straight spokes please see colour options on last page, max. user weight: 125 kg	–
QSS070004	Rear wheel, Proton, 24 spokes (black), straight spokes (L) user weight limit 100kg	XXX,-
QSS070005	Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver	XXX,-
QSS070007	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub: black	XXX,-
QSS070008	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes only 24", only with Adapter, standard and only with 0° & 4° camber; not with push handles, fold down & not without push handles	XXX,-
QSS070009	Rear wheel, Mountainbike, 36 spokes (silver), straight spokes only 24"	XXX,-
QSS070010	Mountainbike wheel, replacement kit incl. tyre (extra coarse profile) only 24"	XXX,-
Axle		
QSS070150	Quick-release axle	–
QSS070153	Quick-release axle superlight, Titan (L) not with Universal or Mountainbike wheel	XXX,-
QSS070151	Quick-release axle, designed for quads not with Rear wheels Spinergy/ Mountainbike/ drum brakes	XXX,-

Tyres choose only one tyre type					€
QSS070101	Tyre, pneumatic, fine profile, Schwalbe Right Run only with 24"				–
QSS070107	Tyre, pneumatic, slick, Schwalbe One (L)				XXX,-
QSS070102	Tyre, puncture proof, tread				–
QSS070104	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution (+)				XXX,-
QSS090163	Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube not with solid tyres				XXX,-

Castor



// QSS080301
Castor, solid (hard), rim: composite,
black



// QSS080302
Castor, solid (soft), rim: composite,
black



// QSS080303
Castor, solid (hard), light-up, rim:
composite, black



// QSS080304
Castor, solid (hard), rim: aluminium
Colour, rim: Silver



// QSS080306
Castor, solid (soft), rim: aluminium
Colour, rim: Silver



// QSS080060
Castor, solid (soft), rim: carbon-com-
posite, ultralight



// QSS080305
Castor, pneumatic, rim: composite,
black

Rear Wheels



// QSS070001
Rear wheel, Universal, 36 spokes (silver), cross-wise spokes



// QSS070002
Rear wheel, Design, 36 spokes (silver), straight spokes



// QSS070003
Rear wheel, Lightweight, 24 spokes (silver), straight spokes



// QSS070004
Rear wheel, Proton, 24 spokes (black), straight spokes



// QSS070005
Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver



// QSS070007
Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub: black



// QSS070008
Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes



// QSS070009
Rear wheel, Mountainbike, 36 spokes (silver), straight spokes



// QSS070153
Quick-release axle superlight, Titan



// QSS070151
Quick-release axle, designed for quads

Tyres



// QSS0701011
Tyre, pneumatic, fine profile, Schwalbe Right Run



// QSS070107
Tyre, pneumatic, slick, Schwalbe One



// QSS070102
Tyre, puncture proof, tread (+)



// QSS070104
Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution

Handrims		€
QSS070201	Handrims, aluminium, powder-coated, clear (L)	–
QSS070202	Handrims, aluminium, anodised, silver (L)	–
QSS070203	Handrims, aluminium, powder-coated, coloured (L) select colour on last page	XXX,-
QSS070204	Handrims, stainless steel, silver	–
QSS070205	Handrims, titanium, untreated	XXX,-
QSS070206	Handrims, titanium, powder-coated, clear	XXX,-
QSS070207	Handrims, Supergrip® (+)	XXX,-
QSS070212	Handrims, Ellipse 3R, black (+) ergonomic handrims with rubber-strip	XXX,-
QSS070105	Handrims, Airy (Omobic), black (+) ergonomic handrims with rubber-strip	XXX,-
QSS070112	Handrims Flow (Omobic), Lizard grip, black (+) tetraplegic handrims with anti-slip coating	XXX,-
Handrim accessories		
QSS070250	Handrim cover, rubber only with 24"	XXX,-
Handrim mounting		
QSS070307	Handrim mount, super narrow (7 mm) only with stainless steel, Ellipse, Flow or Airy	XXX,-
QSS070317	Handrim mount, narrow (17 mm)	–
QSS070327	Handrim mount, wide (27 mm) not with Ellipse, Flow or Airy	–

Wheel Camber choose only one option

QSS070600	0°	–
QSS070602	2° not with Centre of Gravity Adapter, lightweight	–
QSS070604	4° not with Centre of Gravity Adapter, lightweight	–

Spoke Guards only 24"

QSS070831	Transparent					XXX,-	
QSS070823	Princess		QSS070821	Pirates	QSS070816	Happy Lions	XXX,-
QSS070814	Graffiti		QSS070819	Monkeys	QSS070802	Basketball	XXX,-
QSS070820	Oriental		QSS070812	Flying Water	QSS070829	Superhero	XXX,-
QSS070833	Zebra		QSS070817	Hot Wheels	QSS070818	Jungle Fever	XXX,-
QSS070822	Pool		QSS070811	Feather	QSS070827	Silk	XXX,-
QSS070806	Chrome		QSS070830	Tiger	QSS070803	Basketball 2	XXX,-
QSS070813	Galactica		QSS070828	Space			XXX,-

Wheel Locks choose only one option

Wheel Locks		€
QSS060001	Standard brake	–
QSS060002	Knee-lever brake from seat depth 340 mm	XXX,-
	QSS060010 Brake lever extension, knee-lever brake, fold-down	XXX,-
QSS060023	Compact brake EVO (L)	XXX,-
QSS060009	Brake, pull to lock, brake lever, long	XXX,-
QSS060031	Aluminium compact brake (Omobic), black	XXX,-

Handrims



// QSS070202
Handrims, aluminium, anodised, silver



// QSS070203
Handrims, aluminium, powder-coated, coloured



// QSS070207
Handrims, Supergrip®



// QSS070212
Handrims, Ellipse 3R, black



// QSS070105
Handrims, Airy (Omobic), black



// QSS070112
Handrims Flow (Omobic), Lizard grip, black

Wheel Locks



// QSS060001
Standard brake



// QSS060002
Knee-lever brake



// QSS060023
Compact brake EVO



// QSS060009
Brake, pull to lock, brake lever, long + extension (automatically included)



// QSS060031
Aluminium compact brake (Omobic), black

Accessories			€
QSS090000	Tool kit		–
QSS090024	High pressure pump (0 - 11 bar)		XXX,-
QSS090025	Air pump		XXX,-
QSS090030	Back pack		XXX,-
QSS090026	Mobile phone pocket		XXX,-

Frame Colours Handrim: only with handrim colour powder coated; Hub/rim: only for universal wheel

	Frame	Handrim	Hub / Rim		Frame	Handrim	Hub / Rim
01 Matte Black				12 Royal Blue			
02 Gloss Black				13 Brilliant Ocean			
03 Metallic Shadow				14 Scarlet Noir			
04 Starlight White				15 Red Pepper			
05 Silver Ore				16 Cherry Lustre			
06 Storm Gray				17 Turbo Green			
07 Desert Tan				18 Pink Power			
08 Desert Sage				19 Satin Rose			
09 Salute Green				20 Mystic Purple			
10 Velvet Sapphire				21 Solar Orange			
11 Midnight Blue				22 Velocity Yellow			
	–	–	XXX,-		–	–	XXX,-

Sticker Colours QS7 X" on the frame and "Quickie / SOPUR" on sideguard

		€
QSS100550	Black Sticker "Quickie" on sideguard will be silver in combination with carbon sideguards	–
QSS100551	Silver	–
QSS100552	Orange	–
QSS100553	Red	–
QSS100554	Blue	–

Component Colours

		€
QSS100588	Footboard in frame colour option only with aluminium footboard QSS or QSS050050	XXX,-

24 " Light Weight Wheel Colours Hub and rim

		€
	Black	–
QSS100336	Anodized, orange only with 24" Lightweight wheel, Hub and rim anodized, black spokes	XXX,-
QSS100337	Anodized, blue only with 24" Lightweight wheel, Hub and rim anodized, black spokes	XXX,-
QSS100338	Anodized, red only with 24" Lightweight wheel, Hub and rim anodized, black spokes	XXX,-

Accessories

// QSS090000
Tool kit



// QSS090024
High pressure pump (0 - 11 bar)



// QSS090025
Air pump



// QSS090030
Back pack



// QSS090026
Mobile phone pocket

Comments



Quickie QS7 X SA

€

Base Price	00,000.00
Option Price	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day Month Year



Sunrise Medical GmbH

Kahlbachring 2-4

D-69254 Malsch/HD

Phone: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: kundenservice@sunrisemedical.de

www.SunriseMedical.eu