

QS7 X

FF & Hybrid



ORDER FORM

Errors and omissions excepted.

Billing Address Or Customer-ID			Shipping Address	Identical To Billing Address
C.-ID:			C.-ID:	
Name:			Name:	
Street:			Street:	
Town/PC:			Town/PC:	
E-Mail:			E-Mail:	
Phone:			Phone:	
Order Information			User Weight	
Order	Calculation	Built-4-Me®	Max. user weight Art. Nr. 77300000 / 77300001:	110 kg
			Max. user weight Art. Nr. 77300003 / 77300004:	140 kg
Purchase Order No.:			Art. No.: 7730000 DE / 77300001 IC 77300003 DE / 77300004 IC	
Reimbursement code: XXX XXX			User Weight:	kg
Order Date:	Day	Month	Year	Base Price QS7 X
				€ X.XXX,-

Please follow this measuring method:

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.

Seat width (SW):

- a) For the user: Allow space between sideguard and hip on each side of the chair (approx. 10 - 20 mm altogether).
- b) On the wheelchair: From outside frame tube to outside frame tube.

Seat depth (SD):

- a) For the user: When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee (approx. 30 - 60 mm).
Care must be taken with adjustable upholstery.
- b) On the wheelchair: From back tube to front of seat upholstery.

Front seat height:

- a) For the user: When sitting with 90° kneeangle: From back of knee to bottom of shoe sole.
- b) On the wheelchair: From front edge of the seat-upholstry, next to the seat-tube to the ground (min. 30 mm longer than user measurement).
Remember to include the cushion height.

Rear seat height:

On the wheelchair: From rear edge of the seat-upholstry, next to the seat-tube to the ground.

Note: Maximum seat height difference rear to front is 110 mm.

Back height:

- a) For the user: From the seat, including cushion, to be no more than 10 - 20 mm below the lower shoulder blade.
- b) On the wheelchair: From top of the seat-upholstry to upper edge of backrest upholstery.

Note: The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably.

Back angle:

Backrest-angle is measured in relation to the seat. 90°= perpendicular to the seat

Technical Data

Max. overall height:	1050 mm	Max. overall width with:	24" Rear wheel	25" Rear wheel
Max. overall length:	1050 mm	For 0° camber:	SW + 160 mm	SW + 160 mm
Max turning circle:	700 mm	For 2° camber:	SW + 210 mm	SW + 220 mm
Overall weight:	from 9,3 kg (Hybrid frame from 9,8 kg)	For 4° camber:	SW + 260 mm	SW + 270 mm
Weight of heaviest detachable part:	2,1 kg			
Transport weight (without rear wheels):	from 6,7 kg (Hybrid frame from 7,2 kg)			
Max. safe slope (with wheels locked):	10°			

Please note that options such as i.e. armrest-sideguards or drum brakes can increase the overall width of the wheelchair. If this is the case, under certain circumstances, it may not be possible to use some or all of the available escape routes from a building. Please note that in certain configurations the wheelchair may exceed a width of 700 mm. It may be more difficult or impossible to travel on public transport.

Product is **crash tested according to ISO 7176-19** and therefore approved for transportation in a vehicle. In combination with QSF090100 (Tie down brackets for transportation) For transport without a user the chair must be safe clamped in the vehicle.

Please note that the technical data is based on a std configuration:

Standard wheelchair 400 x 400 mm; 6 COG; Seat angle of 8°

Manufacturing tolerances can be +/- 10 mm on all configurations

Lightweight-Guide for ordering the QS7 X

Sunrise Medical offers a wide variety of options for the QS7 X to suit every user's needs. However, most options add significant weight to the wheelchair. To make this more transparent when filling in the order-form, and to assist you in ordering the lightest possible wheelchair, we indicate the approximate weight-increase of the options on the order form:

(L) Indicates the lightest possible option within the range of options

(+) Indicates that this option adds more than 0,5 kg to the wheelchair's weight compared to the lightest option within the range

Frame		€
QSF010067	Open frame, backrest-angle and CoG adjustable, 110 kg max user-weight	XXX,-
QSF010035	Hybrid frame, backrest-angle and CoG adjustable, 140 kg max user-weight not possible with QSF010068 PRO frame	XXX,-

Fixed welded Pro Options

Important Note:

Please check your measurements carefully as a QS7 X is non returnable.


± 10 mm center of gravity adjustment with QSF070703 adapter, standard. Only -10 mm center of gravity adjustment with QSF070701 adapter, lightweight by using additional passive adapter (76807018-ET).

QSF010068	PRO: open frame, fixed welded CoG, 110 kg max. user-weight	XXX,-
-----------	--	-------

Backrest-style for PRO versions

QSF010044	fixed welded, non angle-adjustable, non foldable (choose angle and backrest-height)	–
-----------	---	---

Angle (to the ground)	-4° (86°)	-2° (88°)	0° (90°)	+2° (92°)	+4° (94°)				
Backrest Height (mm)	250	275	300	325	350	375	400	425	450



Client Measurements all measurements without cushion

Seat Width measured from outside frame tube to outside frame tube

(mm)	320	340	360	380	400	420	440	460
	480 ¹	500 ¹						

¹ 480 mm & 500 mm only Hybrid frame QSF010035

Seat Depth

(mm)	340	360	380	400	420	440	460	480	500
	340/ 360/ 380 mm = 380 frame length			each frame in a separate length (20 mm increments)					


Front Seat Height Selection in dependence of front castors and frame versions

Fork Type (mm)		400	410	420	430	440	450	460	470	480	490	500	510	520	530	540	550	560
3"																		
Carbotecture (105 mm)		L	L	–	M	M	H	–	H	XH	XH	–	–	–	–	–	–	–
Aluminium (115 mm)		–	–	L	L	M	M	–	H	H	XH	XH	–	–	–	–	–	–
Aluminium One Arm (105 mm)		L	L	–	M	M	H	–	H	XH	XH	–	–	–	–	–	–	–
Aluminium Frog Leg		–	–	–	L	L	–	M	M	H	H	–	XH	XH	–	–	–	–
4"																		
Carbotecture (105 mm)		L	L	L	M	M	H	–	H	XH	XH	XH	–	–	–	–	–	–
Aluminium (115 mm)		–	–	L	L	M	M	M	H	H	H	XH	XH	XH	–	–	–	–
Aluminium One Arm (105 mm)		L	L	L	M	M	H	–	H	XH	XH	XH	–	–	–	–	–	–
Aluminium Frog Leg		–	–	–	L	L	L	M	M	H	–	H	XH	XH	XH	–	–	–
5"																		
Carbotecture (105 mm)		–	–	–	L	L	M	–	M	H	H	–	XH	XH	–	–	–	–
Aluminium (115 mm)		–	–	–	L	L	L	M	M	H	H	H	XH	XH	XH	–	–	–
Aluminium One Arm (105 mm)		–	–	–	L	L	M	–	M	H	H	–	XH	XH	–	–	–	–
Aluminium Frog Leg		–	–	–	–	–	–	L	L	M	–	M	H	H	XH	XH	–	–
6"																		
Aluminium (115 mm)		–	–	–	–	–	–	L	–	L	M	M	–	H	H	XH	XH	–
Aluminium Frog Leg		–	–	–	–	–	–	–	–	L	–	–	M	–	H	–	–	XH

Carbotecture (105 mm) not approved for use with powered add-ons (powered wheels or push-devices). In general, always verify whether a combination agreement is in place before fitting powered add-ons.

Frame height legend

Low (L)	10" / 255 mm
Med (M)	11" / 280 mm
High (H)	12" / 305 mm
X high (XH)	13" / 330 mm



Rear Seat Height Selection in dependence of rear wheel size

Rear seat height in mm	370	380	390	400	410	420	430	440	450	460	470	480	490	500
24"													²	—
25"	—													²

The rear seat height can only be between 0 mm and 110 mm below the selected front seat height.

² Only with lightweight adapter

Seat Upholstery and Cushions

€

Seat sling

QSF020153	Seat sling, nylon, tension adjustable	—
QSF020154	1 utility bag	—
QSF020155	2 utility bags	XXX,-
QSF020000	Without seat sling	—

Seat cushion, cover with zipper, black

Foam (medium-soft) Thickness:

QSF020101	30 mm	XXX,-
QSF020102	50 mm	XXX,-


JAY Cushions – for complete product range see JAY order forms

JAY020002	JAY Basic	XXX,-
JAY020003	JAY Soft Combi P	XXX,-
JAY020001	JAY Easy Visco C	XXX,-
JAY020007	JAY Easy Fluid	XXX,-
JAY020067	JAY LitePro	XXX,-

Seat accessories

QSF090018	Lap belt with buckle, metal	XXX,-
-----------	-----------------------------	-------

Seat Upholsteries & Cushions




// QSF020101-QSF020102
Seat cushion, cover with zipper,
black



// QSF090018
Lap belt with buckle, metal



// JAY020002
JAY Basic (pictured without cover)



// JAY020003
JAY Soft Combi P (pictured without
cover)



// JAY020001
JAY Easy Visco (pictured without co-
ver)



// JAY020007
JAY Easy Fluid (pictured without
cover)



// JAY020067
JAY Lite Pro (pictured without cover)

Backrest

€

Backrest height range Choose back height factory setting (mm). Height adjustments of in total 50 mm possible.	250 - 300	300 - 350	350 - 400	400 - 450	450 - 475
	250	300	350	400	450
	275	325	375	425	475
	300	350	400	450	

Push handles

QSF030200	Without push handles (L)	–
QSF030201	Push handles, long	–
QSF030203	Push handles, fold-down Only available from back height 325 mm upwards	XXX,-
QSF030204	Push handles, height-adjustable (+) Only available from back height 300 mm upwards	XXX,-
QSF030205	Push handles, height-adjustable for JAY backs (+) Only available from back height 325 mm upwards	XXX,-

Back Design

QSF030010	Backrest, fixed (90°)	–
QSF030029	Desk sideguard preparation for later replacement necessary for retrofitting Desk sideguard or Tubular armrest with Backrest, fixed 90°	–
QSF030031	Backrest, fold-down (to the front) not with height adj. push handles; not with JAY backs; only from back height 350 mm; fixed non-adjustable back height, only from backrest height 400 mm upwards in combination with push handles fold-down; please see Folding Matrix on next side; not with stabilising bar	XXX,-
QSF030012	Backrest, angle-adjustable 75° to 103°	XXX,-
QSF030029	Desk sideguard preparation for later replacement necessary for retrofitting Desk sideguard	–
QSF030031	Backrest, fold-down (to the front) not with height adj. push handles; not with JAY backs; only from back height 350 mm; fixed non-adjustable back height, only from backrest height 400 mm upwards in combination with push handles fold-down; please see Folding Matrix on next side; not with stabilising bar	XXX,-

Backrest upholstery (black)

QSF030316	EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black	–
	Colour binding tape: Black Orange Silver Red Blue	–
QSF030317	EXO PRO Backrest sling, padded, adjustable, breathable, spacer fabric, 1 accessory bag, black	XXX,-
	Colour binding tape: Black Orange Silver Red Blue	–
QSF030300	Without backrest upholstery	–

Backrest accessories

QSF030402	Stabilising bar, automatic folding only available from SW 380 mm; only available from SW 380mm; not with fold-down backrest	XXX,-
-----------	--	-------

Backrest



// QSF030201
Push handles, long



// QSF030203
Push handles, fold-down



// QSF030204
Push handles, height-adjustable



// QSF030031
Backrest, fold-down (to the front)



// QSF030316
EXO Backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag, black



// QSF030317
EXO PRO Backrest sling, padded, adjustable, breathable, spacer fabric, 1 accessory bag, black

Overview of Foldability – 030031 Backrest Foldable Forward

Backheight	350	375	400	425	450	475
Seatwidth						
320						
340						
360						
380						
400						
420						
440						
460						
480						
mm	500					



- The backrest can also be folded 90° forward when the wheelchair is folded.
- In combination with fold-down push handles, the backrest cannot be folded a full 90° forward when the wheelchair is folded.
- When the wheelchair is folded, the backrest cannot be folded a full 90° forward.
- When the wheelchair is folded, the backrest can be folded less than 45° forward. When the wheelchair is unfolded, the backrest can be folded 90° forward.

JAY BACKS plus ACCESSORIES max. user weight 136 kg

€

JAY J3 Back, Shallow Contour	Short (S)		Medium (M)		Tall (T)		
Mid Thoracic			JAY030010	340 mm	JAY030033	380 mm	XXX,-
Upper Thoracic	JAY030012	420 mm	JAY030019	460 mm	JAY030025	500 mm	XXX,-
JAY J3 Back, Mid Contour	Short (S)		Medium (M)		Tall (T)		
Mid Thoracic			JAY030011	340 mm	JAY030018	380 mm	XXX,-
Upper Thoracic	JAY030013	420 mm	JAY030020	460 mm	JAY030026	500 mm	XXX,-
JAY J3 Back - Mid Deep Contour	Short (S)		Medium (M)		Tall (T)		
Mid Thoracic			JAY030044	340 mm	JAY030041	380 mm	XXX,-
Upper Thoracic	JAY030045	420 mm	JAY030042	460 mm	JAY030046	500 mm	XXX,-
JAY J3 Carbon Back - Shallow Contour	Short (S)		Medium (M)		Tall (T)		
Low Thoracic	JAY060001	170 mm			JAY060002	240 mm	XXX,-
Mid Thoracic	JAY060003	300 mm			JAY060004	380 mm	XXX,-

Accessories for JAY J3 backs

WHT030415	Whitmyer Headrest, COBRA XTRA, flip back only with JAY J3 backs; not with J3 Carbon	XXX,-
WHT030416	Whitmyer Headrest, AXYS, Standard only with JAY J3 backs; not with J3 Carbon	XXX,-
JAY090058	JAY Accessory bag only with JAY J3 Carbon backs	XXX,-

Armrests & Sideguards

€

Sideguard

QSF040101	Sideguard, carbon fibre, cold weather protection, slim (black) (L)	XXX,-
QSF040102	Sideguard, aluminium, cold weather protection (in frame colour)	XXX,-
QSF040103	Sideguard, aluminium, fender, cold weather protection (in frame colour)	XXX,-
QSF040104	Sideguard, carbon fibre (EVO), fender, cold weather protection (black)	XXX,-
QSF040107	Sideguard, composite, removable, 210 mm high (310 mm long), black	XXX,-
QSF040120	Sideguard, composite, fender (black)	XXX,-

Armrest

QSF040005	Desk sideguard, flip-up, removable	XXX,-
QSF040001	Desk sideguard, armpad (250 mm), fixed, flip-up, removable	XXX,-
QSF040002	Desk sideguard, armpad (250 mm), height-adjustable, flip-up, removable	XXX,-
QSF040003	Desk sideguard, armpad (310 mm), fixed, flip-up, removable	XXX,-
QSF040004	Desk sideguard, armpad (310 mm), height-adjustable, flip-up, removable	XXX,-
QSF040052	Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable (+)	XXX,-
QSF040117	Tubular armrest. Alu, swing away, padded, removable (+)	XXX,-

Footboards

€

QSF050102	Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), heel loop, black	XXX,-
QSF050126	Swing-away mechanism for platform footboard minimum LLL 30 mm longer	XXX,-
QSF050062	Locking system, platform footplate, aluminium	XXX,-
QSF050059	Footboard, platform, Performance, aluminium, angle-adjustable, flip-up (sideways), calf strap not in SW 320 with 10 mm inset	XXX,-
QSF050131	Footboard, platform, lightweight, plastic, angle & depth-adj. flip up (sideways), calf strap, black	–
QSF050132	Footboard, platform, lightweight, carbon fibre, angle & depth-adj. flip up (sideways), calf strap, black (L)	XXX,-
QSF050061	Footboard, platform, lightw., plastic, angle & depth-adj. flip up (sideways), high-mount L frame = LLL 220 - 320mm; M frame = LLL 220 - 340 mm; H frame = LLL 220 - 360 mm; XH frame = LLL 220 - 380 mm (see front seat height matrix for frame height reference)	XXX,-

Footboard, auto-folding

QSF050027	Footboard, platform, lightweight, carbon fibre, auto-folding, angle & depth-adjustable, calf strap, black	XXX,-
-----------	---	-------

Accessories for footboards

QSF050127	Adjustable side protection, to create a side-stop on the footboard only for footboard, platform lightweight or Performance	XXX,-
-----------	---	-------

JAY BACKS plus ACCESSORIES



// JAY030012
JAY J3 Back, Shallow Contour
Short (S) Upper Thoracic



// JAY030013
JAY J3 Back, Mid Contour Short (S)
Upper Thoracic



// JAY030045
JAY J3 Back - Mid Deep Contour
Short (S) Upper Thoracic

Armrests & Sideguards



// QSF040101
Sideguard, carbon fibre, cold wea-
ther protection, slim (black)



// QSF040102
Sideguard, aluminium, cold wea-
ther protection (in frame colour)



// QSF040103
Sideguard, aluminium, fender, cold
weather protection (in frame colour)



// QSF040104
Sideguard, carbon fibre EVO, fender,
cold weather protection (black)



// QSF040107
Sideguard, composite, removable,
210 mm high (310 mm long), black



// QSF040120
Sideguard, composite, fender
(black)



// QSF040005
Desk sideguard, flip-up, removable



// QSF040001
Desk sideguard, armpad (250 mm),
fixed, flip-up, removable



// QSF040002
Desk sideguard, armpad (250 mm),
height-adjustable, flip-up, removab-
le

Armrests & Sideguards



// QSF040052
Single-post sideguard, armpad (330 mm), depth adjustable, height-adjustable, removable



// QSF040117
Tubular armrest. Alu, swing away, padded, removable

Footboards



// QSF050102
Footboard, platform, aluminium, angle-adjustable, flip-up (sideways), heel loop, black



// QSF050126
Swing-away mechanism for platform footboard



// QSF050062
Locking system, platform footplate, aluminium



// QSF050059
Footboard, platform, Performance, aluminium, angle-adjustable, flip-up (sideways), calf strap



// QSF050131
Footboard, platform, lightweight, plastic, angle & depth-adjustable, flip up (sideways), calf strap, black



// QSF050132
Footboard, platform, lightweight, carbon fibre, angle & depth-adjustable, flip up (sideways), calf strap, black



// QSF050061
Footboard, platform, lightw., plastic, angle & depth-adjustable, flip up (sideways), high-mount



// QSF050027
Footboard, platform, lightweight, carbon fibre, auto-folding, angle & depth-adjustable, calf strap, black



// QSF050127
Adjustable side protection, to create a side-stop on the footboard

Lower Leg Length LLL has to be 25 mm less than front seat height; all measurements are in mm

220	230	240	250	260	270	280	290	300	310	320
330	340	350	360	370	380	390	400	410	420	430
440	450	460	470	480	490	500				

Frame

€

Frame angle

QSF010012	Frame 100 °	–
QSF010013	Frame 92 ° , Performance	XXX,-
QSF010009	Frame 88 ° , Performance Compact	XXX,-

Frame inset

QSF010100	Inset: Standard (0 mm)	–
QSF010101	Inset: Slight (10 mm) not in combination with 88° frame angle	XXX,-

Accessories

QSF110156	Twin-LED light (mounted on caster stems, battery and controller on calf-strap) not in combination with Attitude	XXX,-
QSF010040	Design Frame Protection	XXX,-
QSF010069	Transport Protection, Black Not compatible with desk sideguards (including preparation for desk sideguards), tubular armrest, anti-tips, tip assist, transit wheels, tie down brackets, or frame reinforcement kit	XXX,-
QSF010070	Folding front frame Only with H or XH QSF010067 Open Frame & QSF010100 0 mm inset, not compatible with QSF010012 Frame 100° or 6" castors	XXX,-

Centre of Gravity (CoG) measured from front of back post; choose only one option (0 = very passive 13 = very active)

The CoG scale reflects the distance in cm from the backrest tube to the middle of the axle stem in a 6° configuration. The exact cm measurement will differ slightly depending on the seat height configuration of the chair.

€

QSF070703	Axle plate, standard	0 ¹	1 ²	2	4	5	6	—
		8	9	10	11			
	¹ comes installed with passive adapter							
	² Hybrid frame comes with passive adapter							
QSF070701	Axle plate, lightweight (L) ³	3	4 ²	5	7	8	9	—
		11	12	13				
	² Hybrid Frame comes with passive adapter							
	³ only with 0° wheel camber and not compatible with R90 folding axle adapters COG 3, 12, & 13 not with Hybrid Frame							
QSF010411	Axle plate, aluminium, long (full adjustment)							—
QSF010412	Axle plate, aluminium, shorted (max. adjustment rear seat height +10mm) (L)							XXX,-

Frame Accessories (choose either tip assist or anti- tip tubes per side)

€

QSF090100	Tie down brackets for transportation according to ISO 7176-19 (max. user-weight 136 kg)	XXX,-
-----------	---	-------

Tip assist

QSF090001	Left	QSF090002	Right	XXX,-
QSF090003	Pair			XXX,-
QSF090009	Crutch holder, loop only in combination with tip assist			XXX,-

Anti-tip tubes (+)

QSF090050	Anti-tip tubes, plug-in, pair	XXX,-
-----------	-------------------------------	-------

Anti-tip tubes, swing-away (+)

QSF090004	Left	QSF090005	Right	XXX,-
QSF090006	Pair			XXX,-
QSF090010	Transit wheels			XXX,-

Frame



// QSF010012
Frame 100 °



// QSF010013
Frame 92 °, Performance



// QSF010009
Frame 88 °, Performance Compact



// QSF010070
Folding front frame



// QSF110156
Twin-LED light (mounted on caster stems, battery and controller on calf-strap, details on backpage)



Integrated Li-Ion Battery (3,6V, 2600mAh), on-off-button on controller-box, 3 modes: full power, low power, flash-mode, incl. USB to micro-USB cable for charging



// QSF010069
Transport Protection, Black (light-weight axle plate adapter pictured, also selectable with axle plate adapter, standard)

Centre of Gravity (CoG)



// QSF070703
Axle plate Adapter, standard



// QSF070701
Axle plate Adapter, lightweight

Frame Accessories



// QSF090100
Tie down brackets for transportation according to ISO 7176-19 (max. user-weight 136 kg)



// QSF090009
Crutch holder, loop



// QSF090001
Tip assist



// QSF090010
Transit wheels



// QSF090004
Anti-tip tubes, swing-away



// QSF090050
Anti-tip tubes, plug-in, pair

Preparation Add-on Bike Attitude only with FF-Frame €

QSF090211	Preparation for ATTITUDE, w/o ATTITUDE order docking pins assembled	XXX,-
QSF090212	Preparation for ATTITUDE, with ATTITUDE order docking pins assembled; ATTITUDE order separately	–

Empulse F55 Bundling €

EMO000101	With F55 order, not mounted to the wheelchair Please use separate OF to order the F55	XXX,-
EMO000102	With F55 order, mounted to the wheelchair Please use separate OF to order the F55	XXX,-

Empulse F35 Bundling €

EMO000301	With F35 order, not mounted to the wheelchair Please use separate OF to order the F35	XXX,-
EMO000302	With F35 order, mounted to the wheelchair Please use separate OF to order the F35	XXX,-

Empulse R90 Bundling €

EMO000401	With R90 order, not mounted to the wheelchair Please use separate OF to order the R90, R90 must only be used in combination with anti-tips; R90 only with QSF070703	XXX,-
EMO000402	With R90 order, mounted to the wheelchair Please use separate OF to order the R90, R90 must only be used in combination with anti-tips; R90 only with QSF070703	XXX,-

Empulse M90 Bundling €

EMO000501	With M90 order, not mounted to the wheelchair Please use separate OF to order the M90, M90 must only be used in combination with anti-tips; M90 only with 0° wheel camber	XXX,-
EMO000502	With M90 order, mounted to the wheelchair Please use separate OF to order the M90, M90 must only be used in combination with anti-tips; M90 only with 0° wheel camber	XXX,-

Forks compare with your choice on the first page €

QSF080013	Castor fork, Carbotecture ®, (105 mm) Colour: Black (L) not compatible with power add on devices					–
	Orange	Silver	Red	Blue	Red & Blue only in combination with castor, soft (QSF080302 and QSF080306)	XXX,-
QSF080029	Castor fork, Aluminium, Black (115 mm) (anodized)					–
	Colour Cast. Fork Silver (anodized)					XXX,-
QSF080064	Castor fork, Aluminium, Powder coated (115 mm)					–
	Colour Cast. Fork Starlight White (04)					XXX,-
	Colour Cast. Fork Royal Blue (12)					XXX,-
	Colour Cast. Fork Red Pepper (15)					XXX,-
	Colour Cast. Fork Solar Orange (21)					XXX,-
QSF080008	Castor fork, aluminium, one-arm (105 mm) Colour: Black					XXX,-
	Orange	Silver	Red	Blue		XXX,-
QSF080007	Castor fork, aluminium, suspension (Frog Legs ®), black (+)					XXX,-

Bundling - QS7 X FF/Hybrid: 92° or 100° Frames / 24" & 25" wheels (orange is no go with 25" wheels) / F55 & F35

COG	Seath Depth						
	380	400	420	440	460	480	500
0	✓	✓	✓	✓	✓	✓	✓
1	✓	✓	✓	✓	✓	✓	✓
2	✓	✓	✓	✓	✓	✓	✓
3	✓	✓	✓	✓	✓	✓	✓
4	✓	✓	✓	✓	✓	✓	✓
5	✓	✓	✓	✓	✓	✓	✓
6	✓ 24", x 25"	✓	✓	✓	✓	✓	✓
7	✓ 24", x 25"	✓	✓	✓	✓	✓	✓
8	x	✓ 24", x 25"	✓	✓	✓	✓	✓
9	x	x	✓	✓	✓	✓	✓
10	x	x	✓ 24", x 25"	✓	✓	✓	✓
11	x	x	x	✓	✓	✓	✓
12	x	x	x	✓ 24", x 25"	✓	✓	✓
13	x	x	x	x	✓	✓	✓

Bundling - QS7 X FF/Hybrid: 88° Frames / 24" & 25" Wheels (orange is no go with 25" wheels) / F55 & F35

COG	Seath Depth						
	380	400	420	440	460	480	500
0	✓	✓	✓	✓	✓	✓	✓
1	✓	✓	✓	✓	✓	✓	✓
2	✓	✓	✓	✓	✓	✓	✓
3	✓	✓	✓	✓	✓	✓	✓
4	✓	✓	✓	✓	✓	✓	✓
5	✓	✓	✓	✓	✓	✓	✓
6	✓ 24", x 25"	✓	✓	✓	✓	✓	✓
7	✓ 24", x 25"	✓	✓	✓	✓	✓	✓
8	x	✓ 24", x 25"	✓	✓	✓	✓	✓
9	x	x	✓	✓	✓	✓	✓
10	x	x	✓ 24", x 25"	✓	✓	✓	✓
11	x	x	x	✓ 24", x 25"	✓	✓	✓
12	x	x	x	✓ 24", x 25"	✓	✓	✓
13	x	x	x	x	✓	✓	✓

Forks



// QSF080013
Castor fork, Carbotecture ®, (105 mm) Colour: Black



// QSF080064
Castor fork, Aluminium, Powder coated (115 mm)



// QSF080008
Castor fork, aluminium, one-arm (105 mm) Colour: Black

Castor compare with your choice on the first page

€

QSF080301	Castor, solid (hard), rim: composite, black	–
QSF080302	Castor, solid (soft), rim: composite, black only for 5" & 6" castors	–
QSF080303	Castor, solid (hard), light-up, rim: composite, black only for 4" & 5" castors	–
QSF080304	Castor, solid (hard), rim: aluminium Colour, rim: Silver only for 4" & 5" castors	XXX,-
QSF080306	Castor, solid (soft), rim: aluminium Colour, rim: Silver only for 4" & 5" castors	XXX,-
	Colour rim: Black Orange Red Blue	XXX,-
QSF080060	Castor, solid (soft), rim: carbon-composite, ultralight (L) only for 4" castors	XXX,-

Rear Wheels compare with your choice on the first page

€

QSF070000	Without rear wheel	–
QSF070001	Rear wheel, Universal, 36 spokes (silver), cross-wise spokes (+) only in 24"	–
QSF070002	Rear wheel, Design, 36 spokes (silver), straight spokes only 24"	–
QSF070003	Rear wheel, Lightweight, 24 spokes (silver), straight spokes please see colour options on last page, max. user weight: 125 kg	–
QSF070004	Rear wheel, Proton, 24 spokes (black), straight spokes (L) user weight limit 100 kg	XXX,-
QSF070005	Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver	XXX,-
QSF070007	Rear wheel, Spinergy, 12 spokes (black), straight spokes, hub: black	XXX,-
QSF070008	Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes only 24", only with Adapter, standard and only with 0° & 4° camber; not with push handles, fold down & not without push handles	XXX,-
QSF070009	Rear wheel, Mountainbike, 36 spokes (silver), straight spokes only 24"	XXX,-
QSF070010	Mountainbike wheel, replacement kit incl. tyre (extra coarse profile) only 24"	XXX,-

Axle

QSF070150	Quick-release axle	–
QSF070153	Quick-release axle superlight, Titan (L) not with Universal or Mountainbike wheel	XXX,-
QSF070151	Quick-release axle, designed for quads not with Rear wheels Spinergy/ Mountainbike/ drum brakes	XXX,-

Tyres choose only one tyre type

€

QSF070101	Tyre, pneumatic, fine profile, Schwalbe Right Run only with 24"	–
QSF070107	Tyre, pneumatic, slick, Schwalbe One (L)	XXX,-
QSF070102	Tyre, puncture proof, tread (+)	–
QSF070104	Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution (+)	XXX,-
QSF090163	Spares-kit: 1 extra tyre (same as on the chair) and 1 extra inner tube not with solid tyres	XXX,-

Castor



// QSF080301
Castor, solid (hard), rim: composite,
black



// QSF080302
Castor, solid (soft), rim: composite,
black



// QSF080303
Castor, solid (hard), light-up, rim:
composite, black



// QSF080304
Castor, solid (hard), rim: aluminium
Colour, rim: Silver



// QSF080306
Castor, solid (soft), rim: aluminium
Colour, rim: Silver



// QSF080060
Castor, solid (soft), rim: carbon-com-
posite, ultralight

Rear Wheels



// QSF070001
Rear wheel, Universal, 36 spokes (silver), cross-wise spokes



// QSF070002
Rear wheel, Design, 36 spokes (silver), straight spokes



// QSF070003
Rear wheel, Lightweight, 24 spokes (silver), straight spokes



// QSF070004
Rear wheel, Proton, 24 spokes (black), straight spokes



// QSF070005
Rear wheel, Spinergy, 18 spokes (black), straight spokes; hub: silver



// QSF070007
Rear wheel, Spinergy, 12 spokes (black), straight spokes; hub: black



// QSF070008
Rear wheel, drum brakes, 36 spokes (silver), cross-wise spokes



// QSF070009
Rear wheel, Mountainbike, 36 spokes (silver), straight spokes



// QSF070153
Quick-release axle superlight, Titan



// QSF070151
Quick-release axle, designed for quads

Tyres



// QSF070101
Tyre, pneumatic, fine profile, Schwalbe Right Run



// QSF070107
Tyre, pneumatic, slick, Schwalbe One



// QSF070102
Tyre, puncture proof, tread (+)



// QSF070104
Tyre, pneumatic, puncture proof, tread, Schwalbe Marathon Plus Evolution

Handrims						€
QSF070201	Handrims, aluminium, powder-coated, clear (L)					–
QSF070202	Handrims, aluminium, anodised, silver (L)					–
QSF070203	Handrims, aluminium, powder-coated, coloured (L) select colour on last page					XXX,-
QSF070204	Handrims, stainless steel, silver					–
QSF070205	Handrims, titanium, untreated					XXX,-
QSF070206	Handrims, titanium, powder-coated, clear					XXX,-
QSF070207	Handrims, Supergrip® (+)					XXX,-
QSF070212	Handrims, Ellipse 3R, black (+) ergonomic handrims with rubber-strip					XXX,-
QSF070105	Handrims, Airy (Omobic), black (+) ergonomic handrims with rubber-strip					XXX,-
QSF070112	Handrims Flow (Omobic), Lizard grip, black (+) tetraplegic handrims with anti-slip coating					XXX,-
Handrim accessories						
QSF070250	Handrim cover, rubber only with 24"					XXX,-
Handrim mounting						
QSF070307	Handrim mount, super narrow (7 mm) only with stainless steel, Ellipse, Flow or Airy					XXX,-
QSF070317	Handrim mount, narrow (17 mm)					–
QSF070327	Handrim mount, wide (27 mm) not with Ellipse, Flow or Airy					–
Wheel Camber choose only one option						
QSF070600	0°					–
QSF070602	2° not with Centre of Gravity Adapter, lightweight					–
QSF070604	4° not with Centre of Gravity Adapter, lightweight					–
Spoke Guards only 24" €						
QSF070831	Transparent					XXX,-
QSF070823	Princess	QSF070821	Pirates	QSF070816	Happy Lions	XXX,-
QSF070814	Graffiti	QSF070819	Monkeys	QSF070802	Basketball	XXX,-
QSF070820	Oriental	QSF070812	Flying Water	QSF070829	Superhero	XXX,-
QSF070833	Zebra	QSF070817	Hot Wheels	QSF070818	Jungle Fever	XXX,-
QSF070822	Pool	QSF070811	Feather	QSF070827	Silk	XXX,-
QSF070806	Chrome	QSF070830	Tiger	QSF070803	Basketball 2	XXX,-
QSF070813	Galactica	QSF070828	Space			XXX,-
Wheel Locks choose only one option €						
QSF060001	Standard brake					–
QSF060002	Knee-lever brake from seat depth 340 mm					XXX,-
	QSF060010	Brake lever extension, knee-lever brake, fold-down				XXX,-
QSF060023	Compact brake EVO					XXX,-
QSF060024	Compact brake lightweight EVO (L) not available with 4° camber; not with seat depth 340 mm; not for use with add on bikes; adjustment of CoG +/- approx. 25 mm in the field					XXX,-
QSF060009	Brake, pull to lock, brake lever, long					XXX,-
QSF060031	Aluminium compact brake (Omobic), black					XXX,-

Handrims



// QS F070202
Handrims, aluminium, anodised, silver



// QS F070203
Handrims, aluminium, powder-coated, coloured



// QS F070207
Handrims, Supergrip®



// QS F070212
Handrims, Ellipse 3R, black



// QS F070105
Handrims, Airy (Omobic), black



// QS F070112
Handrims Flow (Omobic), Lizard grip, black

Wheel Locks



// QS F060001
Standard brake



// QS F060002
Knee-lever brake



// QS F060023
Compact brake EVO



// QS F060009
Brake, pull to lock, brake lever, long + extension (automatically included)



// QS F060031
Aluminium compact brake (Omobic), black

Accessories			€
QSF090000	Tool kit		XXX,-
QSF090024	High pressure pump (0 - 11 bar)		XXX,-
QSF090025	Air pump		XXX,-
QSF090030	Back pack		XXX,-
QSF090014	Frame reinforcement kit for add on bike		XXX,-
QSF090026	Mobile phone pocket		XXX,-

Frame Colours Handrim: only with handrim colour powder coated; Hub/rim: only for universal wheel							
	Frame	Handrim	Hub / Rim		Frame	Handrim	Hub / Rim
01 Matte Black				12 Royal Blue			
02 Gloss Black				13 Brilliant Ocean			
03 Metallic Shadow				14 Scarlet Noir			
04 Starlight White				15 Red Pepper			
05 Silver Ore				16 Cherry Lustre			
06 Storm Gray				17 Turbo Green			
07 Desert Tan				18 Pink Power			
08 Desert Sage				19 Satin Rose			
09 Salute Green				20 Mystic Purple			
10 Velvet Sapphire				21 Solar Orange			
11 Midnight Blue				22 Velocity Yellow			
	–	–	XXX,-		–	–	XXX,-

Sticker Colours QS7 X" on the frame and "Quickie / SOPUR" on sideguard			€
QSF100550	Black Sticker "Quickie" on sideguard will be silver in combination with carbon sideguards		–
QSF100551	Silver		–
QSF100552	Orange		–
QSF100553	Red		–
QSF100554	Blue		–

Frame Graphic "Floating Dots" Colour No choice of sticker colours (QSF100550) in combination with frame graphic "floating dots" possible			€
QSF100575	Matte Black		XXX,-
QSF100571	Silver (Grey)		XXX,-
QSF100570	Black		XXX,-
QSF100572	Orange		XXX,-
QSF100573	Red		XXX,-
QSF100574	Blue		XXX,-

Component Colours			€
QSF100588	Footboard in frame colour option only with aluminium footboard QSF050102		XXX,-

24 " Light Weight Wheel Colours Hub and rim			€
	Black		–
QSF100336	Anodized, orange only with 24" Lightweight wheel, Hub and rim anodized, black spokes		XXX,-
QSF100337	Anodized, blue only with 24" Lightweight wheel, Hub and rim anodized, black spokes		XXX,-
QSF100338	Anodized, red only with 24" Lightweight wheel, Hub and rim anodized, black spokes		XXX,-

Accessories



// QSF090000
Tool kit



// QSF090024
High pressure pump (0 - 11 bar)



// QSF090025
Air pump



// QSF090030
Back pack



// QSF090026
Mobile phone pocket

Comments



Quickie QS7 X FF & Hybrid

€

Base Price	00,000.00
Option Price	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day Month Year



Sunrise Medical GmbH

Kahlbachring 2-4
D-69254 Malsch/HD
Phone: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: kundenservice@sunrisemedical.de
www.SunriseMedical.eu