

QUICKIE[®]
Live without limits.



QS 3^x

PREMIUM FOLD. EFFORTLESS
RIDE. SMART VALUE.

QS
SERIES



QUICKIE QS3 X

THE STARTING POINT FOR ACTIVE FOLDING MOBILITY

The **QUICKIE QS3 X** brings the proven performance and modern design of the QS-range into an accessible wheelchair.

Its smooth **FreeFold™** folding makes folding, lifting, and transport easy, while the **lightweight** frame delivers smooth, responsive driving with less effort.

Built-in adjustability, clean styling, and a choice of frame colours create a comfortable, personalised experience that moves beyond the traditional clinical look. The QS3 X combines **contemporary design** with everyday practicality.

With a focused range of options, the QS3 X **simplifies setup**, supports faster delivery, and helps riders access mobility sooner.



FOLDING

SMOOTH FOLDING, EVERYDAY CONVENIENCE

The QS3 X brings proven **FreeFold** performance to a simplified 2-arm cross brace, delivering folding and unfolding performance comparable to premium folding wheelchairs. Effortless to operate, it can be **folded with just one finger** for smooth, simple handling.

LIGHTWEIGHT

LIGHT TO HANDLE, EASY TO TAKE ALONG

With a **transport weight from just 9 kg (SA) and 8,3 kg (FF)**, the QS3 X feels about the same to lift as a fully packed cabin suitcase, making it easy to push, lift, and transport while helping to reduce strain throughout the day. Whether you're heading out for a walk, visiting friends, or travelling by car, its lightweight design makes handling simple - so you can pack up quickly and go.



DUAL ARM CROSS BRACE DESIGN

The dual-arm cross brace design gives the QS3 X a clean, modern look while helping reduce overall weight. The result is a wheelchair that feels lighter to handle, less bulky in everyday use, and stable with every push.



FOLDING FIXATION

The optional folding fixation snap keeps the QS3 X securely folded during transport and storage, making handling easier and more convenient when it is not in use.





PERFORMANCE

EASY TO PROPEL

Designed for smooth, efficient movement, the QS3 X delivers straight-line tracking and a responsive, rigid-feel ride. Its stiff frame and stable wheel alignment help reduce pushing effort for easier, more controlled propulsion throughout the day.

BUILT FOR EVERYDAY RELIABILITY

Built from 6000 series aluminium, the QS3 X combines durability with long-term reliability. Its corrosion-resistant frame is designed to handle everyday life while maintaining a clean, modern look over time.

FORGED CASTOR LINK

The integrated castor design reduces weight, vibration, and service complexity while improving propulsion efficiency and creating a sleek, modern look.



CONFIGURABILITY

CHOOSE THE SETUP THAT FITS YOUR NEEDS

Swing-Away

The Swing-Away frame allows the footrests to move easily to the side, supporting more comfortable transfers. It also helps position the wheelchair closer to beds, chairs, or vehicles and reduces the overall transport size for easier storage and handling.

Fixed-Front

The Fixed-Front frame provides added stability and support, helping maintain a consistent leg position throughout the day. It's ideal for riders who benefit from increased positioning support and a more integrated feel while driving.





HIDDEN ADJUSTABILITY

DESIGNED TO FIT CHANGING NEEDS

Life changes, and your needs can change with it. Alongside seat height, centre of gravity, backrest height, backrest angle and camber adjustment, the seat depth can be adjusted within a 20 mm* range. This helps the QS3 X adapt to different needs and offers practical flexibility over time - now and in the future.

The new axle plate is lighter, offers a cleaner look, and provides full rear seat height adjustment with support for multiple camber positions.



*starting from seat depth 400mm

DESIGN

DESIGNED TO FEEL AS GOOD AS IT LOOKS

Confidence matters. With its modern appearance, the QS3 X feels right at home wherever you go. It's designed to look less clinical and more like something that fits your lifestyle, so you can move through your day feeling comfortable and self-assured.

Choose from **7 frame colors** to personalize your QS3 X. From subtle to more expressive options, the color range lets you create a look that suits your style.



Matte Black




Metallic Shadow



Solar Orange



Storm Grey



Midnight Blue



Red Pepper



Silver Ore





OPTION HIGHLIGHTS

DESK SIDEGUARD

Flip up and removable, armpad available in two size



**NEW COMPOSITE
SIDEGUARD WITH FENDER**

**ANTI-TIPS SA OR
PLUG-IN**



**TIP ASSIST WITH CRUTCH
HOLDER OPTION**



**DESIGN ELEMENT TO
CLOSE THE REAR FRAME**



ANTI-TIPS SA OR PLUG-IN



NEW UPHOLSTERY

The new upholstery offers an improved combination of comfort and breathability, with increased padding for enhanced support.



QUICK DELIVERY

READY WHEN YOU NEED IT

Staying mobile shouldn't depend on long waiting times or complicated processes. Designed with streamlined configurations and reliable availability, the QS3 X helps simplify selection and supporting quick delivery so you can get moving sooner.

TECHNICAL DATA	QS3 X
Max. user weight	130 kg
Seat width	360 - 540 mm
Seat depth	380 - 500 mm
Front seat height	400 - 540 mm
Rear seat height	370 - 500 mm
Backrest height	300 - 500 mm
Backrest angle	from -5 to +15
Total width	SW +240/+260/+270 at 0°, 1.5°, 3° and 24" wheels The total width can exceed 800 mm in certain configurations
Total length	950 mm (FF), 1070 mm (SA)
Total height	1080 mm
Max. safe slope (with wheels locked)	10° (7°)
Frame versions	Swing Away (SA) or Fixed Front (FF)
Frame angle	70° / 80° (SA); 80° (FF)
Front wheel	4", 5", 6", 7"
Rear wheel	22", 24"
Wheelchair weight	From 11,6 kg (SA); 10,9 kg (FF)
Dismantled weight	From 9 kg (SA); 8,3 kg (FF)
Mass of heaviest part	2,3 kg
Removable parts	Rear wheels, SA hangers
Frame colours	7 colours
Crash tested	Yes (approved for transport in a vehicle according ISO 7176-19:2022)
Transport without user	Fix the wheelchair firmly and securely in the vehicle (e.g. car or plane)
Intended use and environment	For people unable to walk and with limited mobility; inside and outside. Operation by user or assistant.



Sunrise Medical GmbH

Kahlbachring 2-4
69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de

For further information on the full specification, options, and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.



Rev. A © 06.2026



www.SunriseMedical.eu

