

100% ROLLATORS



GEMINO

GEMINO – 100% ROLLATOR FOR YOU

Are you looking for a rollator that fits your active lifestyle and allows you to walk comfortably and confidently, indoors and out? Then a Gemino rollator is the perfect choice for you. Our award-winning lightweight rollators set new standards in quality, comfort, safety and design. And let you be yourself!

The Gemino comes in various models, sizes and colours. Naturally, there's also a special/medical Gemino rollator for you, if you need something extra to walk safely. So there's always a Gemino that fits your needs - and personal style. Whatever your choice, you'll always get a 100% rollator. 100% comfort, safety and design.

Great when you can rely on something 100%





100% COMFORT



ENJOYABLE WALKING EXPERIENCE

Our Gemini rollators have been designed with your comfort and enjoyment in mind. A lightweight aluminium frame, high-quality materials and maximum stability ensure a great walking experience, indoors and out. Due to the height adjustment with memory function the ergonomic handles are always in the right position. And if it's time for a break, you can sit safely on the spacious seat. Naturally, you can also take your rollator with you when travelling further – it easily folds to a compact size.

Try the Gemini of your choice, and discover how easy and enjoyable walking can be!

// CHOOSE YOUR FRAME SIZE

Selecting the right size will not only help you walk comfortably and safely. Each frame size also has a different seat height - important for maintaining a correct posture when sitting (For more information, please see Technical Data on page 22).

// SMART PUSH HANDLE ADJUSTMENT


Find out how simple it is to adjust the ergonomic push handles to fit your height. They safely lock into place and due to the special memory function, they're always in the right position for you.

// EASY FOLDING WITH SELF-LOCKING SYSTEM

The Gemini's unique 'fold-click system' makes it very easy to fold and compact to store. Simply put the push handles in the lowest position to reduce the size further. Practical for small spaces and portability by car. By the way, did you know you can also safely manoeuvre in tight places when your rollator is (partially) folded?

// PRACTICAL EVERYDAY USE

With your Gemini you'll enjoy comfortable walking - and comfortable seating, if it's time for a break. Combine with a comfortable back support for additional support. Lots of things to carry? Simply use the practical basket – or choose a stylish bag.





100% SAFETY

WHEREVER YOU GO

Wherever you go, you'll feel confident with a Gemino rollator. Its sturdy, stable frame design ensures maximum safety and security in all settings. And the ergonomically designed push handles provide excellent support. But there's more. Many clever features are proof that we always go the extra mile for your safety. Gemino – enjoy your active live, in total safety!

// EXCELLENT, ULTRA-LIGHT BRAKING

Ergonomically designed, you'll find that operating the brake handles requires very little effort. The smart brake design is easy on the wheels and ensures optimal grip, even when the wheels show signs of wear.

// EASY NEGOTIATION OF OBSTACLES

Kick tubes make it easy to climb kerbs or thresholds, while the protective side bumpers help you avoid obstacles.

// WELL VISIBLE IN THE DARK

Going out at night? 6 integrated reflectors on all sides make sure you're well visible in the dark. For even more safety, simply add the optional light.





GEMINO 20

GEMINO 20

GREAT FUNCTIONALITY AT AN ATTRACTIVE PRICE

Are you looking for a high-quality rollator that offers great value for money? Then the Gemini 20 is just right for you. Simple and safe to use, and it looks good too. With this lightweight rollator you'll enjoy easy and stable walking, while the 'fold-click system' allows for compact storage and easy transport. It's also fitted with a practical basket and spacious seat – great if you want to take a break.

The 'basic' Gemini 20 comes in 3 frame sizes (colour Grey). Choose the model that fits your height and enjoy 100% quality and functionality – at an attractive price.

// 100% GEMINO - SMART CHOICE FOR MOBILITY

Expect a real Gemini walking experience from your Gemini 20 rollator. It has all the technical features of the award-winning Gemini 30, and can be combined with lots of practical accessories.

// EASY STORAGE AND PORTABILITY


Your Gemini 20 simply folds to a compact package due to the unique 'fold-click system'. For practical storage and easy portability by car.






GEMINO 30


THE ORIGINAL 100% ROLLATOR




Silver grey




Champagne




Sapphire blue




Pink



Apple green



Scarlett red



Black

// LIFE IS COLOURFUL

You are looking for individuality? Express your personality with a selection of 7 attractive anodized frame colours.

// ADVANCED BAGS TO MIX AND MATCH

Compliment your Gemini 30 with a stylish and practical bag with shoulder strap – in your favourite colour!



Graphite grey



Blue



Mint green



Pink



Alcantara black



Green



Black



Red



GEMINO 30

ENJOY LIFE AND EXPRESS YOURSELF

Looking for a lightweight rollator that stands out? Gemini 30 - the original 100% rollator – will exceed your expectations when it comes to comfort, safety and design. Its modern design styling and high finishing express quality and class. Naturally, it allows for full personalization to suit your active life – and personal style!

With the award-winning Gemini 30, you'll get 'the best of the best'. Balancing high manoeuvrability with maximum stability, this design rollator will truly enhance your independence. Personalize with your choice of colour and accessories – and enjoy life to the full!

// AWARD-WINNING DESIGN

Our original 100% rollator is recognized by international design institutes like red dot and iF and German consumer magazine Öko-Test.



GEMINO 30 carbon



GEMINO 30 CARBON

THE PREMIUM CARBON ROLLATOR FOR AN ACTIVE LIFESTYLE

May it be something special? Equipped with all the advantages of the Gemini 30, like Gemini's unique `fold-click-system` and the huge range of accessories, the Gemini 30 Carbon is the choice for those who are looking for something special. Particularly light, particularly elegant and particularly maneuverable are the features that characterize our premium rollator and make it really stand out.

// SUPER LIGHTWEIGHT

The low weight of 5.4 (M)/5.6 kg makes the Gemini 30 Carbon easy to lift, easy to transport and store, and also super maneuverable in tight spaces.

// FREEDOM OF LEG MOVEMENT

Parallel front tubes enable an increased seating space because the legs do not hit the tubes.

// IN FRAME INTEGRATED STEPTUBE

The steptube was unobtrusively integrated into the carbon frame and allows to easily overcome curbs.





100% LIGHTNESS



// ANATOMICAL SHAPED HANDGRIPS

Ergonomically designed hand grips provide extra comfort and a secure feel and support while walking due to the extra hand support.



// CARBON FRAME

Carbon is characterized by extreme strength and stiffness combined with low weight and this makes the Gemino 30 Carbon the perfect companion for an active lifestyle



// ADJUSTABLE BACK SUPPORT

The adjustable back support allows to lean back and relax. The built-in seat together with the back support allow for a comfortable rest when needed.



// EASY HANDLING OF PURCHASES

A lot to carry? The advanced net bag can be easily opened with one hand, offers plenty of space for purchases and carries up to 5 kg.



GEMINO 30

Comfort

GEMINO 30 COMFORT

HIGH INDOOR MANOEUVRABILITY WITH OUTDOOR PERFORMANCE

Are you looking for a lightweight rollator that allows you to walk comfortably and confidently, in- and outdoors? Then a Gemino 30 Comfort rollator is the perfect choice for you. Compact with large soft wheels and superb rolling characteristics.

// WHEELS THAT MAKE THE DIFFERENCE

Strolling the cobblestone streets or walking down the forest path. The large, soft, shock absorbing wheels allows you to walk comfortable – on any surface.

// INDOOR AND OUTDOOR USE

The compact frame gives you high indoor manoeuvrability and the large soft wheels allows you to negotiate obstacles with ease whether it is indoor thresholds or outdoor curbs.

// EASILY HANDLING EXTRA WEIGHT

Lots of things to carry? The extra strong basket carries up to 10 kg.



GEMINO 60



// WHEELS THAT MAKE THE DIFFERENCE

The extra-large, soft wheels allow you to negotiate virtually any obstacle with ease. And ensure ultimate walking comfort – on any surface.



// DESIGNED FOR SUPERIOR STABILITY

The Gemini 60's sturdy, extra wide and long frame together with the large wheels ensure excellent stability and safety, also on rough terrain.



// ENJOY MORE SPACE

The large wheelbase makes for maximum stability, whilst offering extra walking space and a very spacious seat – also available with adjustable comfort back support. Great if you want to take a break.



// EASILY HANDLING EXTRA WEIGHT

All Gemini 60 models are suitable for heavy-duty use and take a maximum user weight of 150 kg. Lots of things to carry? The extra strong basket carries up to 10 kg.



GEMINO 60

ULTIMATE COMFORT AND SAFETY – ON ANY TERRAIN

Let's go! Gemini 60 is the perfect companion for your active lifestyle – lightweight yet rock-solid. With its extra-large, soft wheels this design rollator allows you to walk comfortably and securely, also on uneven terrain. Strolling the cobblestone streets in the old town, enjoying a walk in the forest or simply negotiating thresholds in your home? It'll all be so much easier with the robust Gemini 60!

The Gemini 60 is designed for active outdoor use and easily copes with extra weight. At the same time, it quickly folds together – for easy storage and portability by car. Gemini 60 - stay active and enjoy ultimate comfort, wherever you go!



GEMINO 30 Walker



EXTRA SUPPORT AND STABILITY. ALSO GREAT FOR GAIT TRAINING

If you need extra support when walking, the Gemino 30 Walker is the perfect choice for you. This lightweight rollator is fitted with fully adjustable forearm supports, enabling you to walk comfortably and securely. The integrated push handles and brakes are easy to reach and very light to operate. Combining the benefits of the Gemino 30 with enhanced support and stability makes the Gemino 30 Walker also ideal for gait training.

// FOREARM SUPPORTS FOR ENHANCED STABILITY AND SAFETY

You can easily adjust the height, depth and angle of the forearm supports to suit your needs. The brakes and push handles are fully integrated – easy to reach and light to operate. For comfortable and safe walking.

// ADVANCED BRAKING OPTIONS FOR EVEN MORE SECURITY

Looking for even more security? Complement your Gemino 30 Walker with a slowdown brake or innovative SpeedControl braking system. With SpeedControl your rollator will brake automatically when you start to accelerate. You'll keep an upright posture and walk confidently and comfortably close to your rollator.



GEMINO 30 Walker



OPTIMAL SUPPORT AND STABILITY – INDOORS AND OUT

Going out a lot? With the Gemino 60 Walker you get the best of both worlds: The comfortable rolling characteristics and stability of the Gemino 60 and the additional support and safety you'd expect from a 'walker'. Adjust the forearm supports and push handles to fit your needs, and enjoy more freedom indoors - and outdoors. Naturally, you can also use this lightweight rollator for extended gait training!

// ADJUSTABLE FOREARM SUPPORTS FOR ENHANCED STABILITY

The height and angle adjustable forearm supports with integrated push handles and brakes are designed for optimal support and stability, and easy use. Adjust to suit your height and enjoy a comfortable walking position.

// COMFORTABLE AND SAFE OUTDOOR USE

The sturdy frame design coupled with extra-large soft wheels allow you to go out in total comfort and safety. Cobble streets or thresholds in your home? With the Gemino 60 Walker you'll overcome any obstacle with ease.



GEMINO 30

SpeedControl



// ALSO A PERFECT COMPANION FOR CHILDREN

The GEMINO 30 S SpeedControl, the smallest model, offers children every opportunity to discover the world. Thanks to the dynamic braking system they enjoy total safety when walking – making gait training fun!



GEMINO 30 SPEEDCONTROL

THE ROLLATOR INNOVATION


You're looking for a rollator that offers maximum security when walking? Then the GEMINO 30 SpeedControl with its innovative dynamic braking system is just what you need. You'll feel completely secure and enjoy a smooth gait and better, upright position. SpeedControl helps prevent falls and is also great for gait training – for adults and children.

// 3-LEVEL DYNAMIC BRAKING FOR MAXIMUM SECURITY

The innovative SpeedControl braking system features 3-level breaking and can be easily adjusted to your needs. When you start to accelerate, the centrifugal brakes in the rear wheels will automatically be activated. Naturally, you'll experience no resistance when walking at low speed or up slopes. SpeedControl promotes a steady pace, ensures that you walk comfortably, close to the rollator and prevents the rollator from 'running away' from you.

// GREAT AID FOR GAIT TRAINING

It's the GEMINO 30 SpeedControl's unique ability to smooth out your gait through offering controlled resistance that also makes it a great training partner if you need to learn to walk again. You can concentrate fully on the sequence of movements, while your rollator promotes an upright posture and a steady pace. And use the different braking levels for adjusting the rollator to the progress you are making in your rehabilitation.



GEMINO 30

Parkinson

SPECIALLY DESIGNED FOR PEOPLE WITH PARKINSON'S

Some people need a little extra help to be able to walk safely. If you've Parkinson's disease, you may be one of those people. Our Gemino 30 Parkinson rollator has been developed with you in mind and offers features that are fully tailored to your particular needs. Get out and about safely – and with confidence!

// SLOWDOWN BRAKE PROMOTES A STEADY PACE
The Gemino 30 Parkinson is fitted with a slowdown brake, which provides additional resistance and helps you walk at a steady pace. Adjust the level of resistance to your specific needs.

// INCREASED SAFETY WITH REVERSE BRAKE
Due to the reverse braking system the rollator immediately stops moving once the brake is released. This increases safety if you experience 'freezing' while walking - when it feels like your feet are glued to the ground.

// LASER DEVICE HELPS OVERCOME FREEZING
When 'freezing' occurs, you can activate a laser line with a button on the push handles. Through concentrating on the line that appears on the ground, you can then step over it and get moving again.



Left or right
mount



100% CHOICE

PRACTICAL ACCESSORIES FOR YOU

Gemino rollators are not just easily adjustable. They also come with a full range of attractive accessories to make everyday life easier.

A stylish bag, a practical crutch holder, a comfortable back support or a light for extra safety? Choose the options that fit your needs and lifestyle!

// COLOURFUL ADVANCED BAGS

Compliment your Gemino with a stylish, closable bag in your favourite colour. With handles and shoulder strap for easy carrying.

Not available for Gemino 60 M



Blue



Graphite grey



Black



Pink



Green



Red



Alcantara black



Mint green

// BAGS FOR GEMINO 30 S



Blue, black print



Grey, black print



// UMBRELLA****

Rainy weather conditions? With this height-adjustable umbrella you can go out safely, as you can keep both hands on the push handles.



// BACK SUPPORT*

Provides you support to sit comfortably and safely when it's time to take a break. No need to remove when folding the rollator.



// ADJUSTABLE COMFORT BACK SUPPORT**

For additional comfort, choose this large, softly padded back support. Quick to adjust, easily folds with the rollator.



// TABLE TRAY

Practical tray with cup holder allows you to transport food and beverages safely and with ease.



// CRUTCH HOLDER

Keeps your crutch or walking cane securely fastened when walking. Easy to mount, easy to access.



// BASKET FOR OXYGEN CYLINDER***/*

Carrier bag with shoulder strap for safe transport of oxygen cylinders.



// SLOWDOWN BRAKE****

Provides additional resistance to help you walk at a steady pace. Easy to adjust to the desired level – and easy to turn on/off.



// UNIVERSAL BRACKET

Suitable for fitting your choice of personal accessories.




// LIGHT

Ensures you're well visible in the dark – for additional safety.



// INFUSION HOLDER****

Height and angle adjustable drip holder.



Automatic braking at times of sudden acceleration


// SPEEDCONTROL*

With the innovative SpeedControl braking system – available with all Gemino 30 models – you'll feel completely secure when walking. Your rollator will brake automatically when you start to accelerate, ensuring smooth walking and a comfortable, upright position at all times. Naturally, you'll experience no resistance when rolling at normal pace or up slopes. SpeedControl helps prevent falls and is also great for gait training - for adults and children!



// ONE HAND BRAKE

Comfortable and safe usage with one hand.



// REVERSE BRAKE KIT

The rollator immediately stops moving when the brake is released.



// STABILIZER BAR

More safe feeling, easy removable for transport. Only available for Gemino 30 and Gemino 60 Walker

** check availability per model*

*** not available for Gemino 30 / 30 M Walker*

**** not available for Gemino 60 M*

***** not available for Gemino 30 Carbon / 30 M Carbon*



GEMINO – SELECTION GUIDE



GEMINO 20

Light-weight design rollator. Great functionality at an attractive price



Indoor and outdoor use
Everyday use



150 kg
130 kg (M)
125 kg (S)



150-200 cm
135-170 cm (M)
125-165 cm (S)
Also suitable for children



Small turning radius
for great indoor
manoeuvrability



GEMINO 30 Gemino 30 Carbon

Our original 100% rollator. Award-winning design combining high indoor maneuverability with outdoor performance. Lots of colour options and accessories for easy personalization.



Use indoors and out,
mainly in the city
Everyday use



150 kg
130 kg (M)
125 kg (S)



150-200 cm
135-170 cm (M)
125-165 cm (S)
Also suitable for children



Small turning radius
for great indoor
manoeuvrability



Easy personalization with
7 frame colours and
8 bag colours.
For Gemino 30 Carbon only
Carbon black available.



Super lightweight
5.4 (M)/5.6 kg Carbon
version available for easy
transportation and lifting



GEMINO 30 Comfort

Designed for comfortable walking in- and outdoors, also on uneven terrain. Large, soft wheels allow negotiating of obstacles with ease.



Indoor and outdoor use,
also on uneven terrain
Everyday use / Active use



150 kg
130 kg (M)



150-200 cm
135-170 cm (M)



Large, soft wheels ensure
ultimate comfort and easy
negotiation of obstacles



Extra strong bag
(10 kg load capacity)



GEMINO 60

Designed for any kind of terrain and heavy-duty use. Lightweight, yet rock-solid. All the benefits of the award-winning Gemino 30 and more.



Mostly outdoor use, also on
rugged terrain, occasionally
indoors

Active / heavy-duty use



150 kg
150 kg (M)



150-200 cm
135-170 cm (M)



Extra wide and long frame
for ultimate stability



Extra-large and soft wheels
ensure ultimate comfort
and easy negotiation of any
obstacle



Extra wide seat



Extra strong bag
(10 kg load capacity)



GEMINO 30 Walker

Lightweight walker with adjustable forearm supports. The perfect companion if you need extra support and stability indoors. Also great for gait training.



Indoor use
Everyday use /
gait training



150 kg
130 kg (M)



150-200 cm
135-170 cm (M)



Height-adjustable forearm
supports for extra support
and stability



Small turning radius
for great indoor
manoeuvrability



GEMINO 60 Walker

Optimal support and stability – indoors and out. Combines the comfortable rolling characteristics and stability of the Gemino 60 with the additional support and safety you'd expect from a walker.



Indoor and outdoor use,
also on uneven terrain

Everyday use / gait training



130 kg
130 kg (M)



150-200 cm
135-170 cm (M)



Height-adjustable forearm
supports for extra support
and stability



Extra wide and long frame
for ultimate stability



Extra-large and soft wheels
ensure ultimate comfort
and easy negotiation of any
obstacle



Extra strong bag (10 kg load
capacity)



GEMINO 30 SpeedControl

For extra security the brakes in the rear wheels will automatically come into action in case of sudden acceleration. Great for gait training.



Use indoors and out,
mainly in the city

Everyday use



100 kg



150-200 cm
135-170 cm (M)
125-165 cm (S)
Also suitable for children



Small turning radius
for great indoor
manoeuvrability



In case of sudden
acceleration, the
SpeedControl wheels
will automatically give
resistance



GEMINO 30 Parkinson

Specially designed for people with Parkinson's. Featuring a slowdown brake promoting a steady pace, and a reverse braking system and laser light to help with 'freezing'.



Use indoors and out,
mainly in the city

Everyday use



130 kg
130 kg (M)



150-200 cm
135-170 cm (M)



Small turning radius
for great indoor
manoeuvrability



Reverse braking system:
the rollator immediately
stops moving once the
brake is released



Slowdown brake provides
extra resistance and
promotes a steady pace



Laser light can be activated
to help overcome 'freezing'

CHOOSE THE GEMINO ROLLATOR THAT SUITS YOU



GEMINO TECHNICAL DATA

		GEMINO 20			GEMINO 30 GEMINO 30 Carbon			GEMINO 30 Comfort	
		20	20M	20S	30	30M	30S	30	30M
DIMENSIONS							not Carbon		
Maximum user weight	kg	150	130	125	150	130	125	150	130
Recommended user height	cm	150-200	135-170	125-165	150-200	135-170	125-165	150-200	135-170
Maximum load capacity in basket	kg	5	5	5	5	5	5	10	10
Total length A	cm	65	65	68	65	65	68	68	68
Total width B	cm	60	60	55	60	60	55	62	62
Total weight without basket (with basket)	kg	7.1 (7.3)	7.0 (7.2)	6.6 (6.8)	7.1 (7.4) 5.6 (5.9) Carbon	7.0 (7.3) 5.4 (5.7) Carbon	6.6 (6.8)	8 (7.7)	7.8 (7.5)
Width between centre of push handles C	cm	47	47	42	47	47	42	47	47
Seat height D	cm	62	55	48	62	55	48	64	57
Total height E	cm	78 - 100	69 - 88	65 - 77	78 - 100	69 - 88 73 - 88 Carbon	65 - 77	80 - 102	71 - 90
Total height folded	cm	80	76	69	80	76	69	82	77
Total lenght folded	cm	65	65	68	65	65	68	68	68
Total width folded	cm	23	23	23	23	23	23	23	23
Diameter front wheels F	cm	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	23 x 4.2	23 x 4.2
Diameter rear wheels G	cm	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	23 x 4.2	23 x 4.2
PERFORMANCE									
Turning radius	cm	84	84	79	84	84	79	88	88
COLOUR OPTIONS									
Black					✓				
Silver grey					✓ not Carbon	✓ not Carbon	✓	✓	✓
Champagne					✓ not Carbon	✓ not Carbon			
Sapphire blue					✓ not Carbon	✓ not Carbon	✓	✓	✓
Pink					✓ not Carbon		✓		
Apple green					✓ not Carbon	✓ not Carbon			
Scarlett red					✓ not Carbon			✓	✓
CERTIFICATES									
CE		✓	✓	✓	✓	✓	✓	✓	✓
EN 12182		✓	✓	✓	✓	✓	✓	✓	✓
ISO 11199-2		✓	✓	✓	✓	✓	✓	✓	✓
ISO 11199-3									

Gemino Walker rollators: EN 12182, ISO 11199-3 - All other Gemino rollators: EN 12182, ISO 11199-2

		GEMINO 60			GEMINO 30 Walker			GEMINO 60 Walker			GEMINO 30 SpeedControl			GEMINO 30 Parkinson	
		60	60M		30	30M		60	60M		30	30M	30S	30	30M
DIMENSIONS															
Maximum user weight	kg	150	150		150	130		130	130		100	100	100	130	130
Recommended user height	cm	150-200	135-170		150-200	135-170		150-200	135-170		150-200	135-170	125-165	150-200	135-170
Maximum load capacity in basket	kg	10	10		5	5		10	10		5	5	5	5	5
Total length A	cm	74	74		65	65		74	74		65	65	68	65	65
Total width B	cm	64	64		60	60		64	64		60	60	55	60	60
Total weight without basket (with basket)	kg	8.5 (8.8)	8.2 (8.5)		10.6 (10.9)	10.4 (10.7)		11.8 (12.1)	11.6 (11.9)		8.7 (8.9)	8.6 (8.8)	8.2 (8.4)	7.5 (7.8)	7.4 (7.7)
Width between centre of push handles C	cm	51	51		41	41		45	45		47	47	42	47	47
Seat height D	cm	62	53		62	55		62	53		62	55	48	62	55
Total height E	cm	76-95	73 - 83		98 - 115	84 - 100		98 - 115	89 - 98		78 - 100	69 - 88	65 - 77	78 - 100	69 - 88
Total height folded	cm	84	80		97	89		109	101		80	76	69	80	76
Total lenght folded	cm	74	74		65	65		74	74		65	65	68	65	65
Total width folded	cm	23	23		35.5	35.5		35.5	35.5		23	23	23	23	23
Diameter front wheels F	cm	28 x 4.2	28 x 4.2		20 x 3.6	20 x 3.6		28 x 4.2	28 x 4.2		20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6
Diameter rear wheels G	cm	23 x 4.2	23 x 4.2		20 x 3.6	20 x 3.6		23 x 4.2	23 x 4.2		20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6	20 x 3.6
PERFORMANCE															
Turning radius	cm	95	95		84	84		95	95		84	84	79	84	84
COLOUR OPTIONS															
Black															
Silver grey		✓	✓		✓	✓		✓	✓		✓	✓	✓	✓	✓
Champagne															
Sapphire blue													✓	✓	
Pink													✓		
Apple green															
Scarlett red															
CERTIFICATES															
CE		✓	✓		✓	✓		✓	✓		✓	✓	✓	✓	✓
EN 12182		✓	✓		✓	✓		✓	✓		✓	✓	✓	✓	✓
ISO 11199-2		✓	✓		✓	✓		✓	✓		✓	✓	✓	✓	✓
ISO 11199-3															

For further information on the full specification, options and accessories, please refer to the order form.
All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.





Sunrise Medical GmbH

Kahlbachring 2-4

69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-222

E-Mail: international@sunrisemedical.de

www.SunriseMedical.eu

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.eu.

