

Pull Device

Empulse F35



Rev:07/2024
Until:11/2024

ORDER FORM

Errors and omissions excepted.

Billing Address Or Customer-ID			Shipping Address		Identical To Billing Address	
C.-ID:			C.-ID:			
Name:			Name:			
Street:			Street:			
Town/PC:			Town / PC:			
E-Mail:			E-Mail:			
Phone:			Phone:			
Order Information			User Weight			
Order	Calculation	Built-4-Me®	Max. user weight: 110 kg			
Purchase Order No.:			User Weight:		kg	
Art. No.: 4770000X			Reimbursement Code: xxx			
Order Date:	Day	Month	Year			

Base Price € X.XXX,-

Technical Data Overview

Max. user weight	110 kg	Max. speed forward	6 kph; optional 10 kph or 15 kph
Front frame width	260 – 480 mm Frame width measured from center frame tube to center frame tube at the mounting position. Use only the measured dimension and not the seat width to check compatibility	Wheel size	6,5" solid tyre
Total weight with subframe	from 9 kg According to ISO 7176-4.	Battery type	Li-Ion, 36 V, 4 Ah, intelligent battery management system
Total weight with battery	7,5 kg	Brake type	Dual-brake-system: mechanical drum-brake and separate proportional electrical motor brake
Charger	36 V, 2 A	Operating Voltage	36 Volt
Range of action	15 km max	Motor power	6,5" 250 watt (Nominal)
Colour	black	Climbing power	up to 11% (6°)

Base		€
	Pull-device F35 250 W (Nominal) motor, 4 Ah Li-ion battery, folding handles with lock in folded position, forward- and reverse gear, assistive cruise control, dual-brake-system (drum-brake, e-brake), mudguard, front-light	–
Ordering/Bundling		€
EMD090213	Without wheelchair order	–
EMD090217	With wheelchair order, not mounted to the wheelchair Wheelchair order on separate order form. Separate shipment of wheelchair and F35.	–
EMD090218	With wheelchair order, mounted to the wheelchair Wheelchair-order on separate order form. Shipment of readily mounted wheelchair/F35-unit	XXX,-
EMD090182	StVZO Kit incl. Attachment and TÜV approval Germany only	XXX,-
Docking		€
EMD010027	Docking for fixed-front wheelchairs	–
EMD010028	Docking for swing-away wheelchairs	–
Docking-Clamps Always check installation space for the docking clamps before ordering. Depending on configuration and option, a particular configuration of a chair might not offer enough space.		
EMD010030	Docking-clamps for oval tubes (Quickie Xenon ² FF, Helium, Nitrum)	–
EMD010031	Docking-clamps for Quickie Xenon ² SA	–
EMD010029	Docking-clamps for round tube frames (adaptable to 19.5 mm, 23 mm, 25 mm, 28.6 mm, 30mm tubes & QS5 X) Always check combination agreements for non-Sunrise wheelchairs.	–
EMD010062	Docking-clamps for round tube-frames 32 mm	–
Adjustability (Please order the smaller version when a distance is measured between mentioned values)		€
Frame width, from 260mm to 480mm, measured from center frame tube to center frame tube at the mounting position. Use only the measured dimension and not the seat width to check compatibility.		
FW26	Frame width 260 mm	–
FW27	Frame width 270 mm	–
FW28	Frame width 280 mm	–
FW29	Frame width 290 mm	–
FW30	Frame width 300 mm	–
FW31	Frame width 310 mm	–
FW32	Frame width 320 mm	–
FW33	Frame width 330 mm	–
FW34	Frame width 340 mm	–
FW35	Frame width 350 mm	–
FW36	Frame width 360 mm	–
FW37	Frame width 370 mm	–
FW38	Frame width 380 mm	–
FW39	Frame width 390 mm	–
FW40	Frame width 400 mm	–
FW41	Frame width 410 mm	–
FW42	Frame width 420 mm	–
FW43	Frame width 430 mm	–
FW44	Frame width 440 mm	–
FW45	Frame width 450 mm	–
FW46	Frame width 460 mm	–
FW47	Frame width 470 mm	–
FW48	Frame width 480 mm	–

Speed Setting		€
EMD000025	6 kph	–
EMD000027	10 kph Check your local laws and regulations for max. permitted speed in your country. Speed-settings > 6 km/h can only be ordered via Fax with a signed waiver. Max. speed-settings can only be altered by Sunrise Medical.	–
EMD000028	15 kph Check your local laws and regulations for max. permitted speed in your country. Speed-settings > 6 km/h can only be ordered via Fax with a signed waiver. Max. speed-settings can only be altered by Sunrise Medical.	–

Display		€
	Display: Battery charging level, current speed, acceleration mode, speed level, odometer, estimated remaining range	–

Accessories		€
	LED front light	–
	Fender	–
EMD010033	Padding for sub-frame (enclosed)	–
EMD090206	Smart phone holder	XXX,-
EMD090208	Car charging convertor (12V to 230V)	XXX,-
EMD110155	Additional charger	XXX,-
EMD090207	Depth stops including safety locks for docking tubes	XXX,-
EMD090165	Theft-protection, cable-lock	XXX,-

Docking-clamps



// EMD010029 Docking-clamps for round tubes frames (adaptable to 19.5 mm, 23 mm, 25 mm, 28.6 mm & 30 mm tubes & QS5 X)

Accessories



// EMD090206 Smart phone holder



// EMD090208 Car charging converter (12V to 230V)



// EMD090207 Depth stops including safety locks for docking tubes

Letter of Acknowledgement

We herewith advise, that using the Empulse F35 with max. speed-settings higher than 6 kph is only permitted on public roads in certain countries and otherwise may only be used on private property.

In the absence of approval, in accordance with the road traffic laws, the Empulse F35 with speed-settings > 6 kph may not be driven on public roads, bicycle-lanes and pavements. You need to inform yourself about the country specific legislation regarding usage and your legal obligation to insure the Empulse F35 speed settings higher than 6 kph (10/15 kph). When you're travelling with your F35, always check the specific legislation in the respective country.

When placing an order for a max. speed-setting > 6 kph, please confirm with your signature below, that you have understood the above information and that you only intend the usage of your F35 on private property.

I herewith confirm, that I am aware of the information above, that I will follow strictly the country specific legislation.

Date**Customer Signature**

Day

Month

Year

Comments

F35

€

Base Price	00,000.00
Total Of Additional Equipment	
Packaging And Transport	
Exclusive VAT	
VAT	
Total Price	

Date

Signature

Day Month Year

**Sunrise Medical GmbH**

Kahlbachring 2-4
D-69254 Malsch/HD
Phone: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-222
E-Mail: export@sunrisemedical.de
www.SunriseMedical.eu