



QUICKIE® Tango®

January 2021
Rev 01_2021

Invoice to: Name: Street: Town: Postal code: Country: Tel: Fax:		Deliver to: Name: Street: Town: Postal code: Country: Tel: Fax:	
Order date: Purchase order: Number required: <input type="checkbox"/> Order <input type="checkbox"/> Calculation <input type="checkbox"/> Built 4 me Reimbursement code: _____ Art. no.: 034000-000		Marked for: Sunrise order no: <input checked="" type="radio"/> = Standard <input type="checkbox"/> = Cost option <input type="checkbox"/> L = Extended delivery time <input type="checkbox"/> = No cost option Customer service fax number... Max. user weight: 140 kg USER WEIGHT <input type="text" value=""/> kg Mandatory information to ensure correct configuration Basic Price: €	

SEAT DIMENSION AND OPTION PACKS

Seat Frame Width		Depth setting	
TANW033	<input type="radio"/> 33 cm (ajust. 25 a 37cm) (1)	TAND0350	<input type="radio"/> 35 cm
TANW038	<input type="radio"/> 38 cm (ajust. 38 a 42cm) (1)	TAND0375	<input type="radio"/> 37,5 cm
TANW043	<input type="radio"/> 43 cm (ajust. 43 a 48cm)	TAND0400	<input type="radio"/> 40 cm
TANW049	<input type="radio"/> 49 cm (ajust. 49 a 52cm)	TAND0425	<input type="radio"/> 42,5 cm
		TAND0450	<input type="radio"/> 45 cm
		TAND0475	<input type="radio"/> 47,5 cm
		TAND0500	<input type="radio"/> 50 cm

PVP

SEAT

TAN020001	<input checked="" type="radio"/> 5 cm (2") Standard Foam Cushion		
TAN020011	<input type="checkbox"/> Without cushion		
TAN020010	<input type="checkbox"/> Jay Comfort Standard Seat cushion & backrest (2)	(only available in a pack)(Pro Skin cover) (It is not available with 33&38cm seat width)	Jay Comfort Standard Seat cushion & backrest

SEAT

TAN020005	<input type="checkbox"/> Fixed tilt (0° - 9° adjustable, factory setting 6°)	(with tools) (0°, 3°, 6° & 9°)	
TAN020006	<input type="checkbox"/> Manual tilt (3° - 25° adjustable with a locking gas cylinder)		
TAN020007	<input type="checkbox"/> Powered tilt (3° to 25°) (2)		
TAN020008	<input type="checkbox"/> Seat Height Kit (increases STF height by 30mm)	(seat height is going to be increased in 3 cm)	
TAN020009	<input type="checkbox"/> Detachable seat		

BACKRESTS

TAN030000	<input checked="" type="radio"/> Std. Sling Backrest (48cm (19") height)	(3)	(height adjustable, top or bottom in 2,5 cm)
TAN030001	<input checked="" type="radio"/> Mechanical fixed recline Back (6° factory setting)		(with tools) (From factory in 6°) (De -3° a 12°: adjustable in 3°. From 12° to 30°: Adjustable in 6°)
TAN030002	<input type="checkbox"/> Manual adjustable recline Back (-3° to 30°)		
TAN030003	<input type="checkbox"/> Powered recline Back (-3° to 30°) (2)		

HEADREST

TAN030004A	<input type="checkbox"/> Standard Headrest (4)	(Only with std backrest upholstery)	Standard Headrest
TAN030004B	<input type="checkbox"/> Comfort Headrest (4)	(Only with Comfort Pro Seating seat and backrest) (It is not available with 33cm seat width)	Comfort Headrest

ADJUSTABLE LATERALS

TAN030005	<input type="checkbox"/> Adjustable Laterals		
-----------	--	--	--

ARMREST TYPE

TAN040002	<input checked="" type="radio"/> Standard ajustables en anchura 3,5 cm por cada lado, y desmontables		
TAN040004	<input type="checkbox"/> Small armrest necessary for reaching 25 to 32 sizes		
TAN040003	<input type="checkbox"/> Small armrest recommended for 33 to 40 sizes	(Recommended for 33 to 40 cm seat width)	
TAN040001	<input type="checkbox"/> Sideguard reduction pads (5)	(total seat width will be reduced in 5 cm)	

STANDARD LEGRESTS (RubiX2 type)

TAN050000	<input type="checkbox"/> 70° s/a legrest with angle adjustable composite footplate Black Colour		
TAN050004	<input checked="" type="radio"/> 70° s/a legrest, with angle adjustable composite footplate Frame coloured		
TAN050005	<input type="checkbox"/> 70° s/a legrest, with adjustable angle and lenght footplate		
TAN050006	<input type="checkbox"/> 70° s/a short legrest, with adjustable angle and lenght footplate	(distance from seat to footplate from 150mm to 300mm, regulable)	
TAN050001	<input type="checkbox"/> Manual Elevating Legrest with angle adjustable composite footplate		
TAN050002	<input type="checkbox"/> Powered Elevating Legrest with angle adjustable composite footplate (6)	(It needs 2 actuator controls)	
TAN050003	<input type="checkbox"/> 80° Uni style hangers and footplates		
TAN050007	<input type="checkbox"/> Manual Elevating Legrest with angle adjustable footplates		





(1) This option not available with depth setting of 45cm, 47,5cm y 50cm. 33 & 38 cm seat width are not available with Jay Comfort Standard Seat cushion & backrest
 (2) with Jay comfort back & seat we recommended to select bigger seat depth to compensate the thickness.
 Not available with 33 & 38 cm width
 (3) With 33 cm seat width the backrest height is going to be 44.5 cm
 (4) with std backrest upholstery the headrest will be mounted into the pushhandles. With comfort backrest is going to be mounted in the backrest. Not available with 33 cm seat width.
 (5) This option is not available with small armrest
 (6) This option is not available with detachable seat



www.SunriseMedical.es
© Sunrise Medical 2021 1/2

Quickie Tango

● = Standard ○ = Opción sin coste □ = Opción con coste

CONTROLS		PVP
<i>(retraíble en paralelo)</i>		
VR2 controls		 VR2 controls  VR2 remote with actuator control  VR2 remote for lights & indicators
TAN011000	● VR2 remote <i>(without buttons for power actuators)</i>	
TAN011008	□ VR2 remote with actuator control <i>(with 2 buttons for two power actuators)</i>	
TAN011001	□ VR2 remote for lights & indicators <i>(with light and indicators and 2 buttons for power actuators)</i>	
TAN011002	□ VR2 remote Dual <i>(Attendant control always mounted into right side) (you need to select the user control)</i>	
TAN011012	□ VR2 attendant only, front JS mounted to the rear <i>(always mounted in the right side)</i>	
TAN011013	□ VR2 attendant only, front JS with L&I mounted to the rear <i>(light and indicators and 2 buttons for power actuators) (always mounted in the right side)</i>	
TAN011014	□ Power swing away chin control (1) <i>(We recommend to select headrest)</i>	
Mounted side		
TAN011003 ○ right TAN011004 ○ left		
VR2 motor controllers		
TAN011005	● VR2 basic motor controller 50A	
TAN011006	□ VR2 motor controller 60A with 2 act. Outputs <i>(with 2 outputs actuators for 6kph speed motors)</i>	
TAN011011	□ VR2 Heavy-Duty motor controller 90A with 2 act. Outputs <i>(with 2 output actuators) (always for 10kph speed motors)</i>	
Mando acompañante direccional <i>(only attendant is not available with user control)</i>		 Tiller Control with actuator
TAN011009	□ VR2 attendant only, Tiller Control with actuator control	
TAN011010	□ VR2 attendant only, Tiller Control with lights & indicators	
VR2 actuator control box		
TAN011007	□ 4 Actuator Control Box (2) <i>(4 actuator control box for 4 power actuators)</i>	
LIGHTS & INDICATORS		
TAN090001	● Omit Lights and Indicators	 Lights & Indicators
TAN090000	□ Lights & Indicators (3) <i>(always with VR2 remote for lights and indicators)</i>	
KERB CLIMBER		
TAN090003	● Omit Kerb climber	
TAN090002	□ Kerb climber	
FRAME COLOUR		
TAN090004	○ Brilliant Silver	
TAN090005	○ Brilliant Blue	
TAN090006	○ Brilliant Red	
WHEEL SIZE		
TAN070008	● 9" solid castor, 14" solid Knobbly Drive wheels	 Drive wheel mud guards
TAN070007	○ 9" solid castor, 14" pneumatic Knobbly Drive wheels	
TAN070003	○ 9" solid castor, 14" solid V-Type Drive wheels	
TAN070002	○ 9" solid castor, 14" pneumatic V-Type Drive wheels	
TAN070001	○ 8" solid castor, 12" solid Drive wheels	
TAN070000	○ 8" solid castor, 12" pneumatic Drive wheels	
TAN070004	□ Drive wheel mud guards	
EMERGENCY FREE WHEEL MECHANISM		
TAN070005	□ Emergency free wheel mechanism with 8" solid castor and 12" solid rear wheel	
MOTORES		
TAN010011	● 6 Kph	
TAN010010	□ 10 Kph <i>(Always with 90A motor control box.)</i>	
BATTERIES		
TAN016001	● 50Ah Gel Batteries	
TAN016000	○ 40Ah Gel Batteries	
TAN016002	○ Omit Batteries	
CHARGER		
TAN016003	○ 8A Charger UK	
TAN016005	● 8A Charger EU	
TAN016004	○ Omit Charger	
ACCESSORIES		
TAN090007	● Transportation tie down points	 Transportation tie down points
TAN090010	● Lap strap	
TAN090008	□ Crutch holder <i>(It is not available with power recline back and power tilt)</i>	
TAN090009	□ Swing- away tray (4) <i>this option, the control box is going to be mounted in the same side as the remote controller)</i>	
TAN090013	□ Rear View Mirror	
TAN090014	□ Detachable height adjustable anti-pits	
TAN090015	□ control support adjustable in length and height	
TAN090011	□ P&G PP1 hand held programmer	
TAN090012	□ P&G PC mobility programmer	
(1) Not available with dual controler. always with Omni screen. Not available with electrical reclining system (2) 4 Actuator Control Box is going to be mounted in the oposite side of the controler. Is not available with small armrest. (3) This option is not available with detachable seat, it is not available with 35cm depth (4) this option is not available with 33 cm width. Not available with small armrests.		© Sunrise Medical 2021 2/2 Para ampliar información visita nuestra página Web www.SunriseMedical.es