



Edition: 01.21
Valid until: 12.21



WHEELDRIVE

Deliver to:

Name:

Street:

Town:

Postal code:

E-Mail:

Tel:

Fax:

Order date:

Purchase order:

Number required:

Deliver to:

Name:

Street:

Town:

Postal code:

E-Mail:

Tel:

Fax:

Marked for:

Sunrise order no:

☐ Order ☐ Calculation ☐ Built 4 me

• = Standard ☐ = Cost option
○ = No cost option

Customer service fax number: + 31 (0) 492593888

Max. user weight: 130 kg

TECHNICAL DATA

WheelDrive // WHEELCHAIR POWER ASSIST

Max. load capacity	130 kg
Total width compared to manual wheels (set)	43 mm
Weight of batteries (set)	3.6 kg
Total weight including battery (one wheel)	11.5 kg
Total weight excluding battery (one wheel)	9.7 kg
Wheel diameter	24"

Max. speed forward (assist rim)	10 km/h
Max. speed forward (drive rim)	6 km/h
Max. range (assist rim)	20 km
Max. range (drive rim)	12 km
Max. kerb climbing, driving backwards	50 mm
Charging time (full charge)	Approx. 2 hours

Intended use and environment: For people with limited mobility; self-propelling; indoor and outdoor use. Wheelchair type class B

Price in €

WHEELDRIVE SET

1019498	<input type="radio"/>	WheelDrive with 24" Marathon Plus Evolution tyre, pneumatic, puncture proof - right			
1019499	<input type="radio"/>	WheelDrive with 24" Marathon Plus Evolution tyre, pneumatic, puncture proof - left			

WheelDrive mounting (Choose applicable mounting set)

9011415	<input type="radio"/>	Mounting set, quick-release axle diameter 12.7 mm (25-75mm), (only HeliX)			
9011414	<input type="radio"/>	Mounting set, quick-release axle diameter 12.7 mm (for axle plate 25 - 75 mm)			
9012380	<input type="radio"/>	Mounting set, quick-release axle diameter 12.7 mm (25-65mm)			
9011416	<input type="radio"/>	Mounting set, quick-release axle diameter 12 mm, short (for axle plate 37 mm)			
9011417	<input type="radio"/>	Mounting set, quick-release axle diameter 12.7 mm, short (for axle plate 37 mm)			

ACCESSORIES

1019491	<input type="checkbox"/>	Spoke guards 24", set			
1018562	<input type="checkbox"/>	Push rim cover 24", silicone, set (for assist rim)			
1016015	<input type="checkbox"/>	Additional batteries 3.8 Ah, set			

COMMENTS

Wheeldrive

Option price

€.....

Packaging and transport

€.....

Exclusive VAT

€.....

VAT.

€.....

Total Price

€.....

Sunrise Medical Customer service

Tel. +31(0)492 593 888
Fax. +31 (0)492 546 497

E-mail: customerservice@sunrisemedical.nl
Website: sunrisemedical.eu



www.SunriseMedical.eu
© Sunrise Medical 2021