



**QUICKIE**  
Live without limits.

# Easy<sup>max</sup>

Lightweight  
performer...



# Easy Max

## Lightweight performance

"Inside cross brace" ensures maximum rigidity and driving characteristics with a lightweight design



Fix Front version



Swing away version

# Ready for the active users



Active range of frames and adjustments, sideguards, backrests, wheels and footrests, designed to meet all active user needs

## Easy transportation

Easy to fold with retaining strap, removable legrests and rounded rear frame ensures compact transportation



Back angle is adjustable in 3° steps from -9° to +9°.



Triple cross-brace with "inside cross brace" fixation for rigid driving characteristics








New carbon sideguards offers modern style and light weight design



Swing in and out legrests with wide range of divided and platform footrests options



# The lightweight Performer



## Technical data sheet

26 individual frame colours		Maximum user weight	120 kg
6 individual sling colours		Weight of heaviest part	6.40kg
Total Weight	from 10.9kg	Folded width	290 mm

For further information on the full specification, options and accessories, please refer to the order form.  
All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.  
If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com) or alternatively is available on request in large text.



PB\_EASYMAX\_EN\_EU\_REV1

Sunrise Medical GmbH

Kahlbachring 2-4

69254 Malsch/HD

Tel.: +49 (0) 72 53 / 9 80-0

Fax: +49 (0) 72 53 / 9 80-111

[www.SunriseMedical.eu](http://www.SunriseMedical.eu)